

## Laboratory Filtration Products



turning science **into solutions**

Filtration and ultrafiltration are essential process steps in nearly all environmental, chemistry and bioscientific laboratory applications.

Sartorius supplies a wide range of individual filter papers, microporous membranes, filtration devices, ultrafiltration units and protein purification devices to suit these applications. This catalog provides a condensed overview of the Sartorius Lab Filtration product range. Please contact us directly for specialty catalogs – available for in-depth technical information.

## ■ Table of Contents

### □ Ultrafiltration

- 9 Protein Concentration
- 36 DNA Concentration
- 42 Endotoxin Concentration
- 44 Protein Purification
- 50 Virus Purification and Concentration


### □ Filtration Devices

- 60 Minisart® Syringe Filters
- 76 Sartolab® P20 and Sartolab® P20 Plus
- 79 Sartolab® 150v
- 81 Sartolab® RF | BT
- 84 Chemical Compatibility

### □ Basic Filtration

- 91 Filter Papers
- 95 Glass and Quartz Microfiber Filters
- 98 Membrane Filters
- 106 Blotting | Chromatography Papers & Blotting Membranes
- 108 Filtration Equipment
- 118 Combisart® Manifolds













## ■ Ultrafiltration

### Table of Contents

#### ☐ Protein Concentration

- 8 Vivaspin® 500
- 9 Vivaspin® 2
- 11 Centriscart® I
- 12 Vivaspin® Turbo 4
- 13 Vivaspin® 6
- 14 Vivaspin® 15R
- 15 Vivaspin® Turbo 15
- 16 Vivaspin® 20
- 18 Vivaclear Centrifugal Filters
- 19 Vivacell 70
- 22 Vivacell 100
- 24 Vivacell 250
- 26 Vivaflow® 50
- 28 Vivaflow® 50R
- 30 Vivaflow® 200
- 32 Vivapore® Solvent Absorption Concentrators
- 34 Ultrafiltration Membrane Filters

#### ☐ DNA Concentration

- 37 Vivacon® 500
- 40 Vivacon® 2

#### ☐ Endotoxin Concentration

- 43 Vivaspin® Endotest

#### ☐ Protein Purification

- 45 Vivapure® Ion Exchange Protein Purification Products
- 47 Vivapure® Mini and Maxiprep Purification Kits

#### ☐ Virus Purification and Concentration

- 51 Vivapure® Virus Purification and Concentration Kits
- 52 Adenovirus Purification with AdenoPACK Kits
- 54 Lentivirus Purification with LentiSELECT Kit



Ultrafiltration is a convective process using anisotropic semi-permeable membranes to separate macromolecular species and solvents – primarily on the basis of size. By allowing solvents and salts to pass the ultrafiltration membrane while macromolecules are hindered to pass it, they are being concentrated. However, ultrafiltration can also be applied in solvent exchange applications. Multiple concentration and refilling steps will progressively lead to a buffer exchange, replacing lengthy techniques like dialysis. Although ultrafiltration is not a standard method for separating and fractionating macromolecules, it can be used as such if the macromolecules differ at least 10 times in size. Ultrafiltration is a gentle, non-denaturing method that is more efficient and flexible than other processes.

### Ultrafiltration Methods

Sartorius offers you a comprehensive range of ultrafiltration process methods for the concentration of your biological samples. Below is a guide to selecting the most suitable one as this decision depends on your sample volume and the available equipment, as well as your desired filtration speed and process control.

#### – Centrifugal Concentration

##### (100 $\mu$ l to 100 ml Starting Volume)

Centrifugation provides the vector to clear solvent and micro-molecules through the ultrafiltration membrane and into a filtrate container positioned below. This is a gentle process that is quick to set up and offers fast filtration speeds for most solutions. Sartorius offers seven Vivaspin® devices and two Vivacell® and Centrosart® I for protein concentration, as well as Vivacon® devices for DNA and peptide concentration | fractionation.

#### – Pressure Ultrafiltration

##### (5 ml to 250 ml Starting Volume)

Pressurized air or inert gas provide the filtration vector for pressure filtration. For the fastest filtration results, pressurized Vivacell® products can be placed on an orbital laboratory shaker. Agitation prevents macromolecules from polarizing on the membrane surface and thus reducing filtration speed. Vivaspin® 20, Vivacell® 70, Vivacell® 100 and Vivacell® 250 can be run with gas pressure. Pressure-fugation is a unique Sartorius method combining gas pressure with centrifugation. This is the fastest of all methods,

providing process times that are typically 30 to 50 percent faster than centrifugation. You can use Vivaspin® 20 and Vivacell® 70 for pressure-fugation (5 to 50 ml starting volumes).

– **Crossflow (Tangential Flow) (100 ml to 5l Starting Volumes)**

The sample is pumped across an ultrafiltration membrane and then returned to the original reservoir by building up pressure at the outlet of the device. The solution is progressively concentrated as solvent and micro-molecules pass through the membrane into a separate filtrate vessel. Reusable Vivaflow® 50R and Vivaflow® 200, as well as disposable Vivaflow® 50, are offered for your dedicated laboratory crossflow filtration.

– **Static Absorption (3 ml to 20 ml Starting Volume)**

This technique uses an absorbent cellulose pad mounted behind the ultrafiltration membrane to draw solvents and micromolecules through the membrane. The retained macromolecules thus concentrate at the bottom of the sample container. You will not need any additional equipment. These devices are ideal for clinical applications like urine concentration prior to further analysis. Both Vivapore® 5 and Vivapore® 20 offer this procedure.

### Typical Applications for Ultrafiltration

- Concentration | desalting of proteins, enzymes, DNA, monoclonal antibodies, immunoglobulins, viruses and nanoparticles
- Bence Jones Protein concentration from urine samples prior to capillary electrophoresis
- Forensic DNA sample concentration prior to sequencing reaction
- FASP (filter-based sample preparation) – peptide fractionation
- Free drug | hormone assays
- Removal of primers from PCR amplified DNA
- Removal of labeled amino acids and nucleotides
- Deproteinization of samples
- General purpose laboratory concentration and desalting of proteins, enzymes, DNA, biomolecules, viruses, antibodies and immunoglobulins

### Process Optimization

When the highest recoveries are crucial, particularly with solute quantities in the microgram range, Sartorius recommends considering the following tips for optimal ultrafiltration results:

- Select the lowest MWCO membrane that suits your application. For the highest recovery, choose a membrane MWCO which is at least half of the molecular weight of the solute to be retained.
- Avoid over-concentration. The smaller the final concentrate volume, the more difficult it is to achieve complete recovery. If feasible, rinse the device with one or more drops of buffer after the first concentration cycle and then recover it again.
- Pretreat the device overnight in distilled water with a passivation solution such as 5% SDS, Tween 20 or Triton X. Rinse thoroughly before use.





## ■ Vivaspin® 500



### 100 µl to 500 µl Samples


Vivaspin® 500 µl centrifugal filter units offer a simple, one-step procedure for sample preparation. They can effectively be used in fixed-angle rotors accepting 2.2 ml centrifuge tubes.

The patented vertical membrane design and thin channel filtration chamber (US 5,647,990) minimizes membrane fouling and provides high-speed concentrations – even with particle-laden solutions.

### □ Specifications

#### Vivaspin® 500

|                           |                                      |                     |
|---------------------------|--------------------------------------|---------------------|
| Concentrator capacity     | Swing-bucket rotor                   | do not use          |
|                           | Fixed-angle rotor                    | 500 µl              |
| Dimensions                | Total length                         | 50 mm               |
|                           | Width                                | 11 mm               |
|                           | Active membrane area                 | 0.5 cm <sup>2</sup> |
|                           | Hold-up volume, membrane and support | < 5 µl              |
|                           | Dead-stop volume                     | 5 µl                |
| Materials of construction | Body                                 | Polycarbonate       |
|                           | Filtrate vessel                      | Polypropylene       |
|                           | Concentrator cap                     | Polycarbonate       |
|                           | Membrane                             | Polyethersulfone    |




### □ Ordering Information

| Vivaspin® 500 Polyethersulfone                           | Qty./Pkg. | Prod. No. |
|--|-----------|-----------|
| 3,000 MWCO   | 25        | VS0191    |
| 3,000 MWCO   | 100       | VS0192    |
| 5,000 MWCO   | 25        | VS0111    |
| 5,000 MWCO   | 100       | VS0112    |
| 10,000 MWCO  | 25        | VS0101    |
| 10,000 MWCO  | 100       | VS0102    |
| 30,000 MWCO  | 25        | VS0121    |
| 30,000 MWCO  | 100       | VS0122    |
| 50,000 MWCO  | 25        | VS0131    |
| 50,000 MWCO  | 100       | VS0132    |
| 100,000 MWCO   | 25        | VS0141    |
| 100,000 MWCO   | 100       | VS0142    |
| 300,000 MWCO   | 25        | VS0151    |
| 300,000 MWCO   | 100       | VS0152    |
| 1,000,000 MWCO   | 25        | VS0161    |
| 1,000,000 MWCO   | 100       | VS0162    |
| 0.2 µm   | 25        | VS0171    |
| 0.2 µm   | 100       | VS0172    |
| Starter pack<br>(5 of each 5 k, 10 k, 30 k, 50 k, 100 k) | 25        | VS01S1    |

## Vivaspin® 2

### Choice of Membranes



#### 0.4 ml to 2 ml Samples

The Vivaspin® 2 bridges the gap between the 500 µl and 4 ml centrifugal concentrators. This device combines the speed of the classic Vivaspin® products with low internal surface and membrane area for superior recoveries from very dilute solutions.

Available with a choice of PES, Cellulose Triacetate and Hydrosart® membranes, Vivaspin® 2 offers the highest flexibility for process optimization.

Also unique to the Vivaspin® 2 is the choice of directly pipetting the concentrate from the dead-stop pocket built into the bottom of the concentrator or alternatively reverse spinning into the concentrate recovery cap, which can then be sealed for storage. Both methods result in nearly total concentrate recoveries.

#### Specifications

##### Vivaspin® 2


|                           |                            |                     |
|---------------------------|----------------------------|---------------------|
| Concentrator capacity     | Swing-bucket rotor         | 3 ml                |
|                           | Fixed-angle rotor          | 2 ml                |
| Dimensions                | Total length               | 126 mm              |
|                           | Width                      | 17 mm               |
|                           | Active membrane area       | 1.2 cm <sup>2</sup> |
|                           | Hold-up volume of membrane | < 10 µl             |
|                           | Dead-stop volume           | 8 µl                |
| Materials of construction | Body                       | Polycarbonate       |
|                           | Filtrate vessel            | Polycarbonate       |
|                           | Concentrator cap           | Polycarbonate       |
|                           | Membrane                   | PES, CTA, HY        |

#### Ordering Information


| Vivaspin® 2 Polyethersulfone                          | Qty./Pkg. | Prod. No. |
|---|-----------|-----------|
| 3,000 MWCO  | 25        | VS0291    |
| 3,000 MWCO  | 100       | VS0292    |
| 5,000 MWCO  | 25        | VS0211    |
| 5,000 MWCO  | 100       | VS0212    |
| 10,000 MWCO   | 25        | VS0201    |
| 10,000 MWCO   | 100       | VS0202    |
| 30,000 MWCO   | 25        | VS0221    |
| 30,000 MWCO   | 100       | VS0222    |
| 50,000 MWCO   | 25        | VS0231    |
| 50,000 MWCO   | 100       | VS0232    |
| 100,000 MWCO  | 25        | VS0241    |
| 100,000 MWCO  | 100       | VS0242    |
| 300,000 MWCO  | 25        | VS0251    |
| 300,000 MWCO  | 100       | VS0252    |
| 1,000,000 MWCO  | 25        | VS0261    |
| 1,000,000 MWCO  | 100       | VS0262    |
| 0.2 µm  | 25        | VS0271    |
| 0.2 µm  | 100       | VS0272    |
| Starter pack (5 of each 5 k, 10 k, 30 k, 50 k, 100 k) | 25        | VS02S1    |

Integral deadstop  
avoids risk of  
concentrating  
to dryness

PES, CTA, or Hydrosart®  
membranes; Filtrate container  
fits standard 15 ml tube  
carriers



Direct pipette recovery or  
choice of reverse spinning  
concentrate into sample cap



Filtrate and concentrate can  
be sealed for storage

#### Vivaspin® 2 Cellulose Triacetate

|             | Qty./Pkg. | Prod. No. |
|-------------|-----------|-----------|
| 5,000 MWCO  | 25        | VS02U1    |
| 5,000 MWCO  | 100       | VS02U2    |
| 10,000 MWCO | 25        | VS02V1    |
| 10,000 MWCO | 100       | VS02V2    |
| 20,000 MWCO | 25        | VS02X1    |
| 20,000 MWCO | 100       | VS02X2    |

#### Vivaspin® 2 Hydrosart®


|             | Qty./Pkg. | Prod. No. |
|-------------|-----------|-----------|
| 2,000 MWCO  | 25        | VS02H91   |
| 2,000 MWCO  | 100       | VS02H92   |
| 5,000 MWCO  | 25        | VS02H11   |
| 5,000 MWCO  | 100       | VS02H12   |
| 10,000 MWCO | 25        | VS02H01   |
| 10,000 MWCO | 100       | VS02H02   |
| 30,000 MWCO | 25        | VS02H21   |
| 30,000 MWCO | 100       | VS02H22   |

#### Ordering Tips

- Choose a membrane pore size at least 50% smaller than the size of the molecule to be retained.
- It is usually best to select polyethersulfone membranes to achieve the fastest concentrations.
- Usually choose cellulose triacetate for protein removal or ultrafiltrate recovery.
- Usually choose Hydrosart® membranes for the highest recovery with Ig fractions.



## Centrisart® I



Centrisart® I\* is a ready-to-use unit for small-volume, centrifugal ultrafiltration to separate proteins from low molecular weight substances in biological samples.

Centrisart® I features a unique design that enables ultrafiltration in the direction opposite to centrifugal force. This is so effective in preventing premature blockage of the filter that even whole blood samples can be deproteinized.

The ultrafiltrate is collected in the floater insert tube, where it is readily accessible without removing the tube.

Centrisart® I is ideal for the following applications:

- Drug binding studies
- Determination of metabolites in serum
- Protein removal from blood samples
- Cleaning of liposomes
- Virus removal

### Specifications

#### Centrisart® I

|                           |                            |                      |
|---------------------------|----------------------------|----------------------|
| Concentrator capacity     | Swing-bucket rotor         | 2.5 ml               |
|                           | Fixed-angle rotor          | 2.5 ml               |
| Dimensions                | Total length               | 93 mm                |
|                           | Width                      | 14 mm                |
|                           | Active membrane area       | 0.79 cm <sup>2</sup> |
|                           | Hold-up volume of membrane | < 5 µl               |
|                           | Dead-stop volume           | 100 µl               |
| Materials of construction | Centrifuge tube            | Polystyrene          |
|                           | Floater tube               | Cellulose propionate |
|                           | Cap                        | Polyethylene         |
|                           | Membrane                   | CTA, PES             |

### Ordering Information

|   | Qty./Pkg. | Prod. No. |
|---|-----------|-----------|
| 5,000 MWCO CTA  | 12        | 13229-E   |
| 10,000 MWCO CTA                                       | 12        | 13239-E   |
| 20,000 MWCO CTA                                       | 12        | 13249-E   |
| 100,000 MWCO PES                                      | 12        | 13269-E   |
| 300,000 MWCO PES                                      | 12        | 13279-E   |
| Starter pack (3 units each of 5 k, 10 k, 20 k, 100 k) | 12        | 13209-E   |

#### References


P. Nebinger and M. Koel:  
Determination of acyclovir by  
ultrafiltration and high-performance  
liquid chromatography.  
J. Chromatography 619, 342-344 (1993)

F. da Fonseca-Wollheim, K.-G. Heinze,  
K. Lomsky and H. Schreiner:  
Serum ultrafiltration for the elimination  
of endogenous interfering substances  
in creatinine determination.  
J. Clin. Chem. Clin. Biochem. 26, 523-525  
(1988)

R. H. Christenson, S. D. Studenberg,  
S. Beck-Davis and F. A. Sedor:  
Digoxin-like immunoreactivity  
eliminated from serum by centrifugal  
ultrafiltration before fluorescence  
polarization immunoassay of digoxin.  
Clinical Chemistry 33, 606-608 (1987)

\* Centrisart is a registered  
trademark in the U.S. and the  
European Union

## ■ Vivaspin® Turbo 4



### 2 ml to 4 ml Samples

Vivaspin® Turbo 4 is the newest member of the ultrafiltration family and allows the fastest sample concentration with the highest recoveries.

This device can handle up to 4 ml sample volumes in swing-bucket rotors and in fixed-angle rotors that accept 15 ml centrifuge tubes.

The Vivaspin® Turbo 4 optimized design and sleek internal profile ensure maximum process speeds all the way down to the last few microliters, resulting in more than 100-fold concentration.

The UV joining technology ensures smooth joint transition between the membrane and the plating housing – allowing removal of the entire sample concentrated in the unique, pipette-friendly dead-stop pocket.

### □ Specifications

#### Vivaspin® Turbo 4

|                           |  |                              |
|---------------------------|--|------------------------------|
| Concentrator capacity     | Swing-bucket rotor                       | 4 ml                         |
|                           | Fixed-angle rotor                        | 4 ml                         |
| Dimensions                | Total length                             | 122.5 mm                     |
|                           | Width                                    | 17 mm                        |
|                           | Active membrane area                     | 3.2 cm <sup>2</sup>          |
|                           | Hold-up volume of membrane               | < 10 µl                      |
|                           | Dead-stop volume swing-bucket rotor      | 40 µl                        |
|                           | Dead-stop volume fixed-angle rotor (25°) | 60 µl                        |
| Materials of construction | Body                                     | Styrene butadiene copolymere |
|                           | Filtrate vessel                          | Polypropylene                |
|                           | Concentrator cap                         | Polypropylene                |
|                           | Membrane                                 | Polyethersulfone             |

### □ Ordering Information

| Vivaspin® Turbo 4 Polyethersulfone | Qty./Pkg. | Prod. No. |
|------------------------------------|-----------|-----------|
| 3,000 MWCO                         | 25        | VS04T91   |
| 3,000 MWCO                         | 100       | VS04T92   |
| 5,000 MWCO                         | 25        | VS04T11   |
| 5,000 MWCO                         | 100       | VS04T12   |
| 10,000 MWCO                        | 25        | VS04T01   |
| 10,000 MWCO                        | 100       | VS04T02   |
| 30,000 MWCO                        | 25        | VS04T21   |
| 30,000 MWCO                        | 100       | VS04T22   |
| 50,000 MWCO                        | 25        | VS04T31   |
| 50,000 MWCO                        | 100       | VS04T32   |
| 100,000 MWCO                       | 25        | VS04T41   |
| 100,000 MWCO                       | 100       | VS04T42   |

## ■ Vivaspin® 6




### 2 ml to 6 ml Samples

Vivaspin® 6 ml concentrators have been developed to offer increased volume flexibility and performance.

Vivaspin® 6 can process an impressive 6 ml in either swing-bucket or fixed-angle rotors accepting standard 15 ml conical bottom test tubes.

The Vivaspin® 6 features twin vertical membranes for unparalleled filtration speeds and 100× plus concentrations. The remaining volume is easy to read off the printed scale on the side of the concentrator and the modified dead-stop pocket further simplifies direct pipette recovery of the final concentrate.



### □ Specifications

#### Vivaspin® 6

|                           |                            |                     |
|---------------------------|----------------------------|---------------------|
| Concentrator capacity     | Swing-bucket rotor         | 6 ml                |
|                           | Fixed-angle rotor          | 6 ml                |
| Dimensions                | Total length               | 122 mm              |
|                           | Width                      | 17 mm               |
|                           | Active membrane area       | 2.5 cm <sup>2</sup> |
|                           | Hold-up volume of membrane | < 10 µl             |
|                           | Dead-stop volume           | 30 µl               |
| Materials of construction | Body                       | Polycarbonate       |
|                           | Filtrate vessel            | Polycarbonate       |
|                           | Concentrator cap           | Polypropylene       |
|                           | Membrane                   | Polyethersulfone    |

### □ Ordering Information

| Vivaspin® 6 Polyethersulfone                             | Qty./Pkg. | Prod. No. |
|--|-----------|-----------|
| 3,000 MWCO   | 25        | VS0691    |
| 3,000 MWCO   | 100       | VS0692    |
| 5,000 MWCO   | 25        | VS0611    |
| 5,000 MWCO   | 100       | VS0612    |
| 10,000 MWCO  | 25        | VS0601    |
| 10,000 MWCO  | 100       | VS0602    |
| 30,000 MWCO  | 25        | VS0621    |
| 30,000 MWCO  | 100       | VS0622    |
| 50,000 MWCO  | 25        | VS0631    |
| 50,000 MWCO  | 100       | VS0632    |
| 100,000 MWCO   | 25        | VS0641    |
| 100,000 MWCO   | 100       | VS0642    |
| 300,000 MWCO   | 25        | VS0651    |
| 300,000 MWCO   | 100       | VS0652    |
| 1,000,000 MWCO   | 25        | VS0661    |
| 1,000,000 MWCO   | 100       | VS0662    |
| 0.2 µm   | 25        | VS0671    |
| 0.2 µm   | 100       | VS0672    |
| Starter pack<br>(5 of each 5 k, 10 k, 30 k, 50 k, 100 k) | 25        | VS06S1    |



## Vivaspin® 15R




### 4 to 15 ml Samples

Vivaspin® 15R is designed for the volume segment of 2 to 15 ml and features a modified regenerated cellulose membrane; Hydrosart®. This membrane is ideal where extremely high recovery with very low adsorption is needed. Examples of these applications include desalting and concentration of Ig fractions.

### Advantages


- Ultimate recovery with low adsorption (95-98%)
- Exceptionally fast concentration time (30 × in 15 min.)
- Convenient application protocol with easy handling
- Easy scale-up to Vivaflow 200 with Hydrosart® membrane for volumes up to 5 liters
- Very low hold-up volume (< 20 µl)




### Specifications

#### Vivaspin® 15R

|                           |                            |                     |
|---------------------------|----------------------------|---------------------|
| Concentrator capacity     | Swing-bucket rotor         | 15 ml               |
|                           | Fixed-angle rotor          | 12.5 ml             |
| Dimensions                | Total length               | 116 mm              |
|                           | Width                      | 30 mm               |
|                           | Active membrane area       | 3.9 cm <sup>2</sup> |
|                           | Hold-up volume of membrane | < 20 µl             |
|                           | Dead-stop volume           | 30 µl               |
| Materials of construction | Body                       | Polycarbonate       |
|                           | Filtrate vessel            | Polypropylene       |
|                           | Concentrator cap           | Polycarbonate       |
|                           | Membrane                   | Hydrosart®          |



Spin




Recover

### Ordering Information

| Vivaspin® 15R Hydrosart® | Qty./Pkg. | Prod. No. |
|--------------------------|-----------|-----------|
| 2,000 MWCO               | 12        | VS15RH91  |
| 2,000 MWCO               | 48        | VS15RH92  |
| 5,000 MWCO               | 12        | VS15RH11  |
| 5,000 MWCO               | 48        | VS15RH12  |
| 10,000 MWCO              | 12        | VS15RH01  |
| 10,000 MWCO              | 48        | VS15RH02  |
| 30,000 MWCO              | 12        | VS15RH21  |
| 30,000 MWCO              | 48        | VS15RH22  |

## ■ Vivaspin® Turbo 15



### 4 to 15 ml Samples

Vivaspin® Turbo 15 enables the fastest sample concentration with the highest recoveries. This device can handle a sample volume of up to 15 ml in swing-bucket rotors and of up to 11 ml in fixed-angle rotors that accept 50 ml centrifuge tubes.

The UV joining technology ensures smooth joint transition between the membrane and the plating housing – allowing removal of the entire sample concentrated in the unique, pipette-friendly dead-stop pocket.

The optimized design and sleek internal profile of the Vivaspin® Turbo 15 ensure maximum process speeds all the way down to the last few microliters, which results in more than 100-fold concentration.

### □ Specifications

#### Vivaspin® Turbo 15

|  |   |                             |
|--|---|-----------------------------|
| Materials                                | Body                                    | Styrene butadiene copolymer |
|  | Filtrate vessel                         | Polypropylene               |
|  | Concentrator cap                        | Polypropylene               |
|  | Membrane                                | Polyethersulfone (PES)      |
| Dimensions                               | Total length (concentrator insert)      | 77 mm                       |
|  | Total length (in tube with cap)         | 118 mm                      |
|  | Diameter (concentrator insert)          | 27 mm                       |
|  | Active membrane area                    | 7.2 cm <sup>2</sup>         |
|  | Hold-up volume of membrane              | <10 µl                      |
|  | Dead-stop volume for swing-bucket rotor | 100 µl                      |
|  | Dead-stop volume for fixed-angle rotor  | 60 µl                       |
| Concentrator capacity                    | Swing-bucket rotor                      | 15 ml                       |
|  | Fixed-angle rotor (25°)                 | 11 ml                       |
| Maximum speed                            | 4,000 × g                               | 4,000 × g                   |
| Sterilization                            | ETO or 70% EtOH                         |                             |
| Removal of endotoxins [Depyrogenization] | Flushing with 1N NaOH                   |                             |

### □ Ordering Information

| Vivaspin® Turbo 15 Polyethersulfone | Qty./Pkg. | Prod. No. |
|-------------------------------------|-----------|-----------|
| 3,000 MWCO                          | 12        | VS15T91   |
| 3,000 MWCO                          | 48        | VS15T92   |
| 5,000 MWCO                          | 12        | VS15T11   |
| 5,000 MWCO                          | 48        | VS15T12   |
| 10,000 MWCO                         | 12        | VS15T01   |
| 10,000 MWCO                         | 48        | VS15T02   |
| 30,000 MWCO                         | 12        | VS15T21   |
| 30,000 MWCO                         | 48        | VS15T22   |
| 50,000 MWCO                         | 12        | VS15T31   |
| 50,000 MWCO                         | 48        | VS15T32   |
| 100,000 MWCO                        | 12        | VS15T41   |
| 100,000 MWCO                        | 48        | VS15T42   |

## ■ Vivaspin® 20



### 5 to 20 ml Samples

Vivaspin® 20 ml centrifugal concentrators have been developed to offer increased volume flexibility and performance.

Vivaspin® 20 handles up to 20 ml in swing-bucket centrifuges and 14 ml in 25° fixed-angle rotors that accept 50 ml centrifuge tubes.



Featuring twin vertical membranes for unparalleled filtration speeds, the Vivaspin® 20 can achieve 100x plus concentrations.


The remaining volume is easy to read off the printed scale on the side of the concentrator and the modified dead-stop pocket further simplifies direct pipette recovery of the final concentrate.

### More Process Flexibility

Vivaspin® 20 is available with unique accessories and operating methods that are designed to provide more process flexibility and further time savings.

### Gas Pressure Filtration

When an appropriate centrifuge is unavailable or for single sample processing, Vivaspin® 20 can be filled with up to 15 ml and then pressurized for bench-top concentration. For even faster processing, gas pressure can be combined with centrifugal force. "Pressure-fugation" is particularly suitable for difficult or viscous samples, such as serum, or for use of a low process temperature, which reduces filtration speed, and generally when minimum process time is essential.



## □ Specifications

### Vivaspin® 20

|                           |                            |                           |
|---------------------------|----------------------------|---------------------------|
| Concentrator capacity     | Swing-bucket rotor         | 20 ml                     |
|                           | Fixed-angle rotor          | 14 ml                     |
|                           | With pressure head         | 15 ml                     |
| Dimensions                | Total length               | 116 mm                    |
|                           |                            | 125 mm with pressure head |
|                           | Width                      | 30 mm                     |
|                           | Active membrane area       | 6.0 cm <sup>2</sup>       |
|                           | Hold-up volume of membrane | < 20 µl                   |
|                           | Dead-stop volume           | 50 µl                     |
| Materials of construction | Body                       | Polycarbonate             |
|                           | Filtrate vessel            | Polycarbonate             |
|                           | Concentrator cap           | Polypropylene             |
|                           | Pressure head              | Acetal   Aluminum         |
|                           | Membrane                   | Polyethersulfone          |


## □ Ordering Information

| Vivaspin® 20 Polyethersulfone | Qty./Pkg. | Prod. No. |
|-------------------------------|-----------|-----------|
| 3,000 MWCO                    | 12        | VS2091    |
| 3,000 MWCO                    | 48        | VS2092    |
| 5,000 MWCO                    | 12        | VS2011    |
| 5,000 MWCO                    | 48        | VS2012    |
| 10,000 MWCO                   | 12        | VS2001    |
| 10,000 MWCO                   | 48        | VS2002    |



| <b>Vivaspin® 20 Polyethersulfone</b>  | <b>Qty./Pkg.</b> | <b>Prod. No.</b> |
|---|------------------|------------------|
| 30,000 MWCO   | 12               | VS2021           |
| 30,000 MWCO   | 48               | VS2022           |
| 50,000 MWCO   | 12               | VS2031           |
| 50,000 MWCO   | 48               | VS2032           |
| 100,000 MWCO  | 12               | VS2041           |
| 100,000 MWCO  | 48               | VS2042           |
| 300,000 MWCO  | 12               | VS2051           |
| 300,000 MWCO  | 48               | VS2052           |
| 1,000,000 MWCO  | 12               | VS2061           |
| 1,000,000 MWCO  | 48               | VS2062           |
| 0.2 µm  | 12               | VS2071           |
| 0.2 µm  | 48               | VS2072           |
| Starter pack<br>(2 each of 5 kDa, 10 kDa, 30 kDa, 50 kDa,<br>100 kDa, 0.2 µm) | 12               | VS20S1           |

| <b>Vivaspin® 20 Accessories</b> | <b>Qty./Pkg.</b> | <b>Prod. No.</b> |
|---------------------------------|------------------|------------------|
| Air pressure controller (APC)   | 1                | VCA002           |
| Charge valve for pressure head  | 1                | VCA005           |
| Diafiltration cups              | 12               | VSA005           |
| Female connector                | 1                | VCA010           |
| Male connector                  | 1                | VCA011           |
| 4 mm OD pneumatic tube (3 m)    | 1                | VCA012           |
| Vivaspin® 20 pressure head      | 1                | VCA200           |



Using the Vivaspin® 20 pressure cap

## ■ Vivaclear Centrifugal Filters



Vivaclear centrifugal filters are disposable microfiltration devices for the fast and reliable clarification|filtration of biological samples in the range 100 µl to 500 µl. They can be used in fixed-angle rotors accepting 2.2 ml centrifuge tubes.

### Product Features

- High-flux polyethersulfone membrane
- 0.8 µm pore size
- Low hold-up volume (<5 µl)
- Fast and reproducible performance

### Applications

- Clarification of samples before loading in Vivapure® protein purification spin columns
- Removal of particles and precipitates
- Filtration of plasma and serum
- Filtration of cells or cell debris

## □ Specifications

### Vivaclear Centrifugal Filters

|                           |                                       |                      |
|---------------------------|---------------------------------------|----------------------|
| Rotor                     | 40–45° fixed-angle rotor, 500 µl      |                      |
| Pore size                 | 0.8 µm                                |                      |
| Dimensions                | Total length                          | 43 mm                |
|                           | Filtrate collection tube diameter     | 11 mm                |
|                           | Active membrane area                  | 0.34 cm <sup>2</sup> |
|                           | Hold-up volume, membrane plus support | < 5 µl               |
|                           | Maximum RCF                           | 2,000 × g            |
| Materials of construction | Body                                  | Polypropylene        |
|                           | Membrane                              | Polyethersulfone     |
|                           | Filtrate collection tube              | Polypropylene        |

## □ Ordering Information

|                           | Qty./Pkg. | Prod. No. |
|---------------------------|-----------|-----------|
| Vivaclear Mini 0.8 µm PES | 100       | VK01P042  |

## ■ Vivacell 70



### 10 ml to 70 ml Samples

Vivacell 70 combines the ease of use of centrifugal devices with the flexibility and control provided by pressurized ultrafiltration cells. Vivacell 70 is inexpensive, quick and easy to assemble, requires no tubing connections or stirring mechanisms and can be adapted to equipment availability or to specific user preferences.

For more process control or for single samples, combine gas pressure with a gentle orbital shake, or you can even pressurize and then leave standing on a bench top or in a refrigerator for the highest simplicity with minimum equipment requirements.

The longitudinal membrane inhibits fouling, while the built-in dead stop will hinder further concentration when the residual volume drops below 150  $\mu$ l.



For convenience, simply spin in a large capacity centrifuge (rotors accepting 250 ml bottles). For the highest speeds, particularly with difficult samples, pressurize the device with air or inert gas before centrifuging.

### □ Specifications

#### Vivacell 70

|                           |                                       |   |
|---------------------------|---------------------------------------|---|
| Concentrator capacity     | Swing-bucket rotor                    | 70 ml   |
|                           | Fixed-angle rotor                     | 50 ml   |
|                           | With pressure head                    | 70 ml   |
|                           | With pressure-fuge head               | 50 ml   |
| Dimensions                | Total length                          | 119 mm standard centrifugal<br>185 mm with pressure head<br>125 mm with pressure<br>fuge head |
|                           | Width                                 | 62 mm   |
|                           | Active membrane area                  | 20 cm <sup>2</sup>  |
|                           | Hold-up volume of membrane            | < 200 $\mu$ l   |
|                           | Dead-stop volume                      | 150 $\mu$ l   |
|                           |                                       |   |
| Operating requirements    | Rotor type                            | Swing-bucket or fixed angle   |
|                           | Minimum rotor angle                   | 25°   |
|                           | Rotor cavity                          | To fit 250 ml (62 mm)<br>centrifuge bottles   |
|                           | Maximum speed                         | 1,000 g   |
|                           | Maximum pressure                      | 5 bar (75 psi)  |
| Materials of construction | Body                                  | Polycarbonate   |
|                           | Filtrate vessel                       | Polycarbonate   |
|                           | Concentrator cap                      | Santoprene  |
|                           | Pressure head   pressure<br>fuge head | Acetal  |
|                           | Membrane                              | Polyethersulfone  |

### Ordering Information


#### **Vivacell 70 Polyethersulfone – Concentrator Bodies With Polycarbonate Filtrate Bottles**

|              | <b>Qty./Pkg.</b> | <b>Prod. No.</b> |
|--------------|------------------|------------------|
| 5,000 MWCO   | 2                | VS6011           |
| 10,000 MWCO  | 2                | VS6001           |
| 30,000 MWCO  | 2                | VS6021           |
| 50,000 MWCO  | 2                | VS6031           |
| 100,000 MWCO | 2                | VS6041           |
| 0.2 µm       | 2                | VS6071           |


#### **Vivacell 70 Polyethersulfone – Concentrator Body Only**

|              | <b>Qty./Pkg.</b> | <b>Prod. No.</b> |
|--------------|------------------|------------------|
| 5,000 MWCO   | 10               | VS6012           |
| 10,000 MWCO  | 10               | VS6002           |
| 30,000 MWCO  | 10               | VS6022           |
| 50,000 MWCO  | 10               | VS6032           |
| 100,000 MWCO | 10               | VS6042           |
| 0.2 µm       | 10               | VS6072           |


| <b>Vivacell 70 Accessories</b>   | <b>Qty./Pkg.</b> | <b>Prod. No.</b> |
|--|------------------|------------------|
| Air pressure controller (APC), complete with pressure gauge, regulator, over-pressure safety valve, female connector to Sartorius Stedim Biotech pressure products and 1 m extension line (4 mm pneumatic tubing) with male and female connectors and 1 m of 6 mm inlet tubing | 1                | VCA002           |
| 250 ml centrifuge bottle – standard caps   | 4                | VSA003           |
| Modified caps for use in fixed-angle rotors with 250 ml centrifuge bottles   | 2                | VCA004           |
| Charge valve for pressure-fuge head  | 1                | VCA005           |
| Replacement seals for pressure-fuge head (VCA701)  | 10               | VCA007           |
| Female connector   | 1                | VCA010           |
| Male connector   | 1                | VCA011           |
| 4 mm pneumatic tubing (3 m)  | 1                | VCA012           |
| Vivacell 70 pressure head with reservoir and filtrate bottle (bench-top use)   | 1                | VCA700           |
| Vivacell 70 pressure-fuge head (for use in centrifuge)   | 2                | VCA701           |

**Centrifuge**


- Process convenience
- Low shear, no foaming
- Less visual control

**Pressure**

- Simplicity and highest process control
- Ideal for refrige-rated use
- Slower concentrations

**Pressure-Shake**

- Speed and process control
- Ideal for single samples
- If left unattended can concentrate to dryness

**Pressure-Fuge**

- Fastest processing
- Ideal with low MWCO or with difficult solutions
- Less visual control



# Vivacell 100



## 20 ml to 100 ml Samples



Vivacell 100 is the latest member of the Vivacell family and bridges the volume range between the Vivacell 70 and the Vivacell 250.

The patented vertical membrane design ensures the highest performance and unmatched flexibility.



Vivacell 100 is a unique and innovative concentrator for volumes from 20 ml to 100 ml, and utilizes pressure, centrifuge or pressure-shake to rapidly concentrate even samples with very high-particle loading.

Vivacell 100 is designed for centrifugal concentration of samples up to 100 ml, which makes it the largest centrifugal unit available. At the same time, its new design allows for maximum centrifugal force of  $4,000 \times g$  to be used for even faster concentration.



Filtrate container fits standard 250 ml rotors

## Vivacell 100 Utilizes:

- Pressure
- Centrifuge
- Pressure-shake

Like the smaller Vivacell 70 unit, Vivacell 100, when used as a centrifugal device, fits only into rotors that accept 250 ml bottles.

Vivacell 100 units can also be used for single or extremely sensitive samples in the pressurized mode only and left on a bench or placed on a swing-bucket laboratory shaker for faster concentration. It can also be kept in a pressurized mode in the refrigerator. Handling is facilitated by the use of quick connectors. In whichever mode Vivacell 100 is used, the vertical membrane design inhibits membrane fouling, while the built-in dead stop impedes concentration to dryness and loss of sample.


## Specifications

### Vivacell 100

|                           |                            |  |
|---------------------------|----------------------------|--|
| Concentrator capacity     | Swing-bucket rotor         | 90 ml  |
|                           | With pressure head         | 98 ml  |
| Dimensions                | Total length               | 123 mm centrifugal<br>197 mm with pressure head                        |
|                           | Width                      | 62 mm  |
|                           | Active membrane area       | 23.5 cm <sup>2</sup>   |
|                           | Hold-up volume of membrane | < 250 µl   |
|                           | Dead-stop volume           | 350 µl   |
| Operating requirements    | Rotor type                 | Swing-bucket   |
|                           | Rotor cavity               | To fit 250 ml (62 mm) centrifuge bottles (maximum cavity depth 105 mm) |
|                           | Maximum speed              | 2,000 g  |
|                           | Maximum pressure           | 5 bar (75 psi)   |
| Materials of construction | Body                       | Polycarbonate  |
|                           | Filtrate vessel            | Polycarbonate  |
|                           | Concentrator cap           | Santoprene   |
|                           | Pressure head              | Acetal   |
|                           | Membrane                   | Polyethersulfone   |


### □ Ordering Information

| <b>Vivacell 100 Polyethersulfone<br/>with Polypropylene Concentrator Cap</b>   | <b>Qty./Pkg.</b> | <b>Prod. No.</b> |
|--|------------------|------------------|
| 5,000 MWCO   | 2                | VC1011           |
| 5,000 MWCO   | 10               | VC1012           |
| 10,000 MWCO  | 2                | VC1001           |
| 10,000 MWCO  | 10               | VC1002           |
| 30,000 MWCO  | 2                | VC1021           |
| 30,000 MWCO  | 10               | VC1022           |
| 50,000 MWCO  | 2                | VC1031           |
| 50,000 MWCO  | 10               | VC1032           |
| 100,000 MWCO   | 2                | VC1041           |
| 100,000 MWCO   | 10               | VC1042           |
| 300,000 MWCO   | 2                | VC1051           |
| 300,000 MWCO   | 10               | VC1052           |
| 1,000,000 MWCO   | 2                | VC1061           |
| 1,000,000 MWCO   | 10               | VC1062           |
| 0.2 µm   | 2                | VC1071           |
| 0.2 µm   | 10               | VC1072           |
| <b>Vivacell 100 Accessories</b>  | <b>Qty./Pkg.</b> | <b>Prod. No.</b> |
| Air pressure controller (APC), complete with pressure gauge, regulator, over-pressure safety valve, female connector, 1 m extension line (4 mm pressure tubing) with male and female connectors and 1 m of 6 mm inlet tubing | 1                | VCA002           |
| Female connector   | 1                | VCA010           |
| Male connector   | 1                | VCA011           |
| 4 mm pressure tubing (3 m)   | 1                | VCA012           |
| Santoprene replacement seals   | 10               | VCA014           |
| Vivacell 100 pressure head with replacement seals (5)  | 1                | VCA800           |




#### **Centrifuge**

- Process convenience
- Low shear, non-foaming
- Less visual control



#### **Pressure**

- Simplicity and the highest process control
- Ideal for refrigerated use
- Slower concentrations



#### **Pressure-Shake**

- Speed and process control
- Ideal for single samples

# Vivacell 250



### 50 ml to 250 ml Samples

The Vivacell 250 is a totally new concept for the concentration of relatively high-volume biological samples. This product offers numerous advantages when compared to stirred cells.

- One size handles a volume range from under 50 ml to 250 ml.
- Use free standing on a bench top or in a refrigerator for maximum simplicity, or use a laboratory shaker for fastest concentrations.
- The unique, conical dead stop built into the bottom of the membrane insert allows concentrations to under 1 ml.


- The gentle vortex action controls membrane polarization, while greatly reducing the shear effects that are typical of stirring mechanisms.
- Setup or membrane replacement takes just a few seconds. Quick-connect fittings and a simple screw closure further enhance ease of use.

This unequaled membrane module takes only seconds to replace. The concentrate can be easily monitored through the graduated inspection window.

## Specifications

### Vivacell 250

|                           |                                   |                    |
|---------------------------|-----------------------------------|--------------------|
| Concentrator capacity     | 250 ml                            |                    |
| Max pressure              | 4 bar (60 psi)                    |                    |
| Dimensions                | Width                             | 116 mm             |
|                           | Height (incl. pressure indicator) | 235 mm             |
|                           | Active membrane area              | 40 cm <sup>2</sup> |
|                           | Hold-up vol. memb. & support      | < 200 µl           |
|                           | Dead-stop volume                  | 600 µl             |
| Materials of construction | Screw closure                     | Acetal             |
|                           | Pressure head                     | Acetal             |
|                           | Quick release connector           | Acetal             |
|                           | Concentrator body   sleeve        | Polycarbonate      |
|                           | Filtrate container                | Polycarbonate      |



## □ Ordering Information

| <b>Vivacell 250</b>   | <b>Qty./Pkg.</b> | <b>Prod. No.</b> |
|---|------------------|------------------|
| Vivacell 250 complete with pressure head, pressure indicator   over-pressure safety valve, quick-release connection to APC, 2 sample reservoirs, filtrate container and insert tool   | 1                | VCA250           |
| <b>Vivacell 250 Polyethersulfone Inserts</b>  | <b>Qty./Pkg.</b> | <b>Prod. No.</b> |
| 5,000 MWCO  | 5                | VC2511           |
| 10,000 MWCO   | 5                | VC2501           |
| 30,000 MWCO   | 5                | VC2521           |
| 50,000 MWCO   | 5                | VC2531           |
| 100,000 MWCO  | 5                | VC2541           |
| 0.2 µm  | 5                | VC2571           |
| Starter kit (1 each of 5 kDa, 10 kDa, 30 kDa, 50 kDa, 100 kDa)  | 5                | VC25S1           |
| <b>Vivacell 250 Accessories</b>   | <b>Qty./Pkg.</b> | <b>Prod. No.</b> |
| Air pressure controller (APC), complete with pressure gauge, regulator, over-pressure safety valve, female connector to Sartorius Stedim Biotech pressure products and 1 m extension line (4 mm pneumatic tubing) with male and female connector and 1 m of 6 mm inlet tubing | 1                | VCA002           |
| Replacement pressure indicator   over-pressure safety valve   | 1                | VCA008           |
| Vivacell 250 maintenance kit (includes one sample reservoir and filtrate container and O-ring seals for pressure head)  | 1                | VCA009           |
| Female connector  | 1                | VCA010           |
| Male connector  | 1                | VCA011           |
| 4 mm OD pressure tubing (3 m)   | 1                | VCA012           |
| Replacement pressure head and screw closure   | 1                | VCA015           |

# Vivaflow<sup>®</sup> 50



## 100 ml to 3 l Samples

The unique, patented Vivaflow<sup>®</sup> 50 system provides ease of use, performance, flexibility and economy that are unrivaled by any laboratory or pilot-scale filtration system on the market.


## Unique Features

- The thin-channel, flip-flow recirculation path provides high crossflow velocities with minimum pump requirements
- No need for pressure holders
- Crystal clear for simple checking of remaining hold-up volume and membrane status
- Unique interlocking modules with series connectors for easy scale up
- Disposable

## Unique Performance

- A single 50 cm<sup>2</sup> module will typically reduce 500 ml to less than 15 ml in under 50 minutes
- Less than 10 ml minimum system recirculation for the highest concentrations
- Less than 500 µl non-recoverable hold-up volume
- Nearly total recoveries achievable with a single 10 ml rinse


Unique "flip-flow" thin channel flow path results in high turbulence and linear velocity for exceptional flux even at high concentrations




## Specifications

### Vivaflow<sup>®</sup> 50

|                           |                              |                          |
|---------------------------|------------------------------|--------------------------|
| Dimensions                | Overall L×H×W                | 107 mm×84 mm×25 mm       |
|                           | Channel W×H                  | 15 mm×0.3 mm             |
|                           | Active membrane area         | 50 cm <sup>2</sup>       |
|                           | Hold-up volume (module)      | 1.5 ml                   |
|                           | Minimum recirculation volume | < 10 ml                  |
|                           | Non-recoverable hold-up      | < 0.5 ml                 |
| Operating conditions      | Pump flow                    | 200 ml/min to 400 ml/min |
|                           | Maximum pressure             | 3 bar (45 psi)           |
|                           | Maximum temperature          | 60°C                     |
| Materials of construction | Main housing                 | Polycarbonate            |
|                           | Flow channel                 | TPX (PMP)                |
|                           | Membrane support             | TPX (PMP)                |
|                           | Seals and O-rings            | Silicone                 |
|                           | Pressure indicator           | Polypropylene, SS spring |
|                           | Flow restrictor              | Polypropylene            |
|                           | Fittings                     | Nylon                    |
|                           | Tubing                       | PVC (medical grade)      |



Multiple modules



Single module



## □ Ordering Information

| <b>Vivaflow® 50*</b> | <b>Qty./Pkg.</b> | <b>Prod. No.</b> |
|----------------------|------------------|------------------|
| 3,000 MWCO PES       | 2                | VF05P9           |
| 5,000 MWCO PES       | 2                | VF05P1           |
| 10,000 MWCO PES      | 2                | VF05P0           |
| 30,000 MWCO PES      | 2                | VF05P2           |
| 50,000 MWCO PES      | 2                | VF05P3           |
| 100,000 MWCO PES     | 2                | VF05P4           |
| 1,000,000 MWCO PES   | 2                | VF05P6           |
| 0.2 µm PES           | 2                | VF05P7           |
| 100,000 MWCO RC      | 2                | VF05C4           |

\* Vivaflow® 50 modules include filtrate tube, size 16 peristaltic tubing, flow restrictors and fittings.

### **Vivaflow® 50 Complete System Comprises**

|  |   |        |
|--|---|--------|
| Pump (240 V), Easy Load pump head (size 16), tubing, 500 ml sample   diafiltration reservoir, module stand, pressure indicator, T-connectors, series interconnectors | 1 | VFS502 |
| Pump (115 V), Easy Load pump head (size 16), tubing, 500 ml sample   diafiltration reservoir, module stand, pressure indicator, T-connectors, series interconnectors | 1 | VFS504 |

### **PVC Tubing and Fittings**

|  |        |
|--|--------|
| Size 16 PVC pump tubing (3 m, 3.2 × 1.6 mm)  | VFA004 |
| Flow restrictor set (2 × 0.4 mm, 0.6 mm, 0.8 mm)   | VFA009 |
| T-connectors for running 2 stacks (2 units)  | VFA030 |
| Series interconnectors (6 units)   | VFA031 |
| Female luer fittings (10 units)  | VFA032 |
| VF50 tubing kit (2 × 1 m size 16 PVC tubing with inlet fittings, 2 × 50 cm size 16 PVC tubing with 0.6 mm flow restrictors, 1 × series interconnector) | VFA034 |
| Flow restrictor 0.6 mm (6 units)   | VFA035 |

### **Accessories**

|  |        |
|--|--------|
| Masterflex economy drive variable speed peristaltic pump (240 V) | VFP001 |
| Masterflex economy drive variable speed peristaltic pump (115 V) | VFP002 |
| 500 ml sample and   or diafiltration reservoir                   | VFA006 |
| Masterflex Easy Load pump head – size 16                         | VFA012 |
| Vivaflow® 50 stand   | VFA016 |
| Pressure indicator (1 to 3 bar [15 psi to 45 psi])               | VFA020 |


# Vivaflow® 50R




**0.1 ml to 1 l Samples**  
Concentrate 100 ml to under 20 ml in just a few minutes or concentrate one liter 50 times in less than 60 minutes. Alternatively, speed up your process by using two Vivaflow® 50R units in parallel and concentrate 1 liters in under 30 min.

- Fast and easy protein sample concentration
- Reusable
- Concentrates volumes from 0.1 l to 1 l
- Optimal for concentration of culture supernatants and viruses
- The most compact crossflow cassette with a premium Hydrosart® membrane

Vivaflow® 50R is a plug-and-play laboratory crossflow cassette for concentrating up to 1 L aqueous samples. The active membrane area per device is 50 cm². One unit comes with all the necessary accessories for running the device with a laboratory pump and a size 16 pump head. For speeding up concentration, two cassettes can be run simultaneously.



Vivaflow® 50R – Single module



Vivaflow® 50R – Two modules

## Specifications

### Vivaflow® 50R

|                           |                           |                          |
|---------------------------|---------------------------|--------------------------|
| Dimensions                | Overall L×H×W             | 100 mm×100 mm×24 mm      |
|                           | Channel W×H               | 7.5 mm×0.4 mm            |
|                           | Active membrane area      | 50 cm²                   |
|                           | Hold-up volume (module)   | 1.7 ml                   |
|                           | Min. recirculation volume | 10 ml                    |
|                           | Non-recoverable hold-up   | < 0.5 ml                 |
| Operating conditions      | Pump flow                 | 200 ml/min to 400 ml/min |
|                           | Maximum pressure          | 4 bar (60 psi)           |
|                           | Maximum temperature       | 60°C                     |
| Materials of construction | Main housing              | Acrylic                  |
|                           | Flow channel              | Acrylic                  |
|                           | Membrane support          | Polypropylene            |
|                           | Seals and O-rings         | Silicone                 |
|                           | Pressure indicator        | Polypropylene, SS spring |
|                           | Flow restrictor           | Polypropylene            |
|                           | Fittings                  | Nylon                    |
|                           | Tubing                    | PVC (medical grade)      |

## □ Ordering Information

| <b>Vivaflow® 50R*</b>   | <b>Qty./Pkg.</b> | <b>Prod. No.</b> |
|-------------------------|------------------|------------------|
| 5,000 MWCO Hydrosart®   | 1                | VF05H1           |
| 10,000 MWCO Hydrosart®  | 1                | VF05H0           |
| 30,000 MWCO Hydrosart®  | 1                | VF05H2           |
| 100,000 MWCO Hydrosart® | 1                | VF05H4           |

\* Vivaflow® 50R modules include pressure indicator, flow restrictor and size 16 pvc peristaltic tubing and fittings.

### **Vivaflow® 50R/200 Complete System**

|  |   |        |
|--|---|--------|
| Pump (240 V), Easy Load pump head (size 16), tubing, 500 ml sample   diafiltration reservoir | 1 | VFS202 |
| Pump (115 V), Easy Load pump head (size 16), tubing, 500 ml sample   diafiltration reservoir | 1 | VFS204 |


### **Tubing and Fittings**

|   |    |        |
|---|----|--------|
| Size 16 pvc pump tubing and Luer fittings (3 m, 3.2 × 1.6 mm) | 1  | VFA004 |
| T-connectors for running 2 units                              | 2  | VFA030 |
| Flow restrictor set (2 × 0.4 mm, 0.6 mm, 0.8 mm)              | 6  | VFA009 |
| Female luer fittings – size 16                                | 10 | VFA032 |
| Flow restrictors 0.6 mm                                       | 6  | VFA035 |
| Female luer fittings – size 15                                | 10 | VFA036 |

### **Accessories**

|  |   |        |
|--|---|--------|
| Masterflex economy drive variable speed peristaltic pump (240 V) | 1 | VFP001 |
| Masterflex economy drive variable speed peristaltic pump (115 V) | 1 | VFP002 |
| 500 ml sample and   or diafiltration reservoir                   | 1 | VFA006 |
| Masterflex Easy Load pump head – size 16                         | 1 | VFA012 |

# Vivaflow<sup>®</sup> 200




### 0.5 ml to 5 l Samples

Concentrate 250 ml to under 20 ml in just a few minutes or concentrate one litre 50 times in less than 30 minutes. Alternatively, use two Vivaflow<sup>®</sup> 200 units in parallel and concentrate 5 liters in under 75 minutes.

Nearly total sample recoveries can be expected with most solutions.

The standard economy package comes complete with tubing, pressure indicator, flow restrictor and high-pressure pump tubing. All you need is a peristaltic pump capable of handling 6.4 mm OD (size 16) tubing. Should your pump head require larger tubing, just use the interconnector provided to attach your own peristaltic tubing to this standard product.

Two modules in parallel will concentrate 5 liters in under 75 minutes



Vivaflow 200 setup for diafiltration

## Specifications

### Vivaflow 200

|                           |                           |                          |
|---------------------------|---------------------------|--------------------------|
| Dimensions                | Overall L×H×W             | 126 mm×138 mm×38 mm      |
|                           | Channel W×H               | 10 mm×0.4 mm             |
|                           | Active membrane area      | 200 cm <sup>2</sup>      |
|                           | Hold-up volume (module)   | 5.3 ml                   |
|                           | Min. recirculation volume | < 20 ml                  |
|                           | Non-recoverable hold-up   | < 2 ml                   |
| Materials of construction | Main housing              | Acrylic                  |
|                           | Flow channel              | Acrylic                  |
|                           | Membrane support          | Polypropylene            |
|                           | Seals and O-rings         | Silicone                 |
|                           | Pressure indicator        | Polypropylene, SS spring |
|                           | Flow restrictor           | Polypropylene            |
|                           | Fittings                  | Nylon                    |
|                           | Tubing                    | PVC (medical grade)      |
| Operating conditions      | Pump flow                 | 200 ml/min to 400 ml/min |
|                           | Maximum pressure          | 4 bar (60 psi)           |
|                           | Maximum temperature       | 60°C                     |

## □ Ordering Information

| <b>Vivaflow® 200*</b>   | <b>Qty./Pkg.</b> | <b>Prod. No.</b> |
|-------------------------|------------------|------------------|
| 3,000 MWCO PES          | 1                | VF20P9           |
| 5,000 MWCO PES          | 1                | VF20P1           |
| 10,000 MWCO PES         | 1                | VF20P0           |
| 30,000 MWCO PES         | 1                | VF20P2           |
| 50,000 MWCO PES         | 1                | VF20P3           |
| 100,000 MWCO PES        | 1                | VF20P4           |
| 0.2 µm PES              | 1                | VF20P7           |
| 2,000 MWCO Hydrosart®   | 1                | VF20H9           |
| 5,000 MWCO Hydrosart®   | 1                | VF20H1           |
| 10,000 MWCO Hydrosart®  | 1                | VF20H0           |
| 30,000 MWCO Hydrosart®  | 1                | VF20H2           |
| 100,000 MWCO Hydrosart® | 1                | VF20H4           |

\* Vivaflow® 200 modules include pressure indicator, flow restrictor and size 16 PVC peristaltic tubing and fittings.

### **Vivaflow® 200 Complete System Comprises**

|  |   |        |
|--|---|--------|
| Pump (240 V), Easy Load pump head (size 16), tubing, 500 ml sample   diafiltration reservoir | 1 | VFS202 |
| Pump (115 V), Easy Load pump head (size 16), tubing, 500 ml sample   diafiltration reservoir | 1 | VFS204 |

### **Accessories**

|  |        |
|--|--------|
| Masterflex economy drive variable speed peristaltic pump (240 V) | VFP001 |
| Masterflex economy drive variable speed peristaltic pump (115 V) | VFP002 |
| 500 ml sample and   or diafiltration reservoir                   | VFA006 |
| Masterflex Easy Load pump head – size 16                         | VFA012 |
| Masterflex Easy Load pump head – size 15                         | VFA013 |

### **Tubing and Fittings**

|  |        |
|--|--------|
| Size 15 PVC pump tubing and Luer fittings (3 m, 4.8 mm × 2.6 mm) | VFA003 |
| Size 16 PVC pump tubing and Luer fittings (3 m, 3.2 mm × 1.6 mm) | VFA004 |
| Y-connetor (size 15 to 2 × size 16)                              | VFA005 |
| Flow restrictor set (2 × 0.4 mm, 0.6 mm, 0.8 mm)                 | VFA009 |
| Female luer fittings, size 16 (10 units)                         | VFA032 |
| Flow restrictors 0.6 mm (6 units)                                | VFA035 |
| Female luer fittings, size 15 (10 units)                         | VFA036 |



# Vivapore® Solvent Absorption Concentrators




**3 ml to 20 ml Samples**  
With no need for additional equipment, pressure or vacuum, solvent absorption is the most economic and user-friendly concentration technique available to the clinician and research scientist.


Just fill the unit with the solution to be concentrated, wait for the desired concentration level to be achieved and then pipette the concentrated sample from the bottom of the reservoir.

Vivapore® is ideal for general-purpose laboratory concentration and purification prior to further analysis. It is particularly suited for labile solutions that can denature with alternative shear- or pressure-inducing methods or that require processing in a cold room environment.

Vivapore® concentrators extend the solvent absorption technique to a totally new level of performance, application potential and ease of use.



Vivapore® 5



Vivapore® 10 | 20

## Specifications

|                              | Vivapore® 5        | Vivapore® 10   20      |
|------------------------------|--------------------|------------------------|
| Membrane material            | PES                | PES                    |
| Membrane MWCO                | 7,500              | 7,500                  |
| Membrane surface area        | 20 cm <sup>2</sup> | 28 cm <sup>2</sup>     |
| Reservoir material           | SAN                | SAN                    |
| Volume range                 | 1 ml to 5 ml       | 2 ml to 10 ml   20 ml* |
| Minimum concentrate volume   | 50 µl              | 50 µl                  |
| Vivapore® overall dimensions |                    |                        |
| Width (mm)                   | 42                 | 46                     |
| Height (mm)                  | 82                 | 100                    |

## □ Ordering Information

| <b>Vivapore® 5*</b> | <b>Qty./Pkg.</b> | <b>Prod. No.</b> |
|---------------------|------------------|------------------|
| 7,500 MWCO PES      | 4                | VP0503           |
| 7,500 MWCO PES      | 30               | VP0501           |

\* includes stand and recovery pipettes

### **Requires Stand**

|                |     |        |
|----------------|-----|--------|
| 7,500 MWCO PES | 100 | VP0502 |
|----------------|-----|--------|

### **Vivapore® 10 | 20\***

|                |    |        |
|----------------|----|--------|
| 7,500 MWCO PES | 4  | VP2003 |
| 7,500 MWCO PES | 30 | VP2001 |

\* includes stand and recovery pipettes

### **Requires Stand**


|                |     |        |
|----------------|-----|--------|
| 7,500 MWCO PES | 100 | VP2002 |
|----------------|-----|--------|

### **Accessories**

|   |     |        |
|---|-----|--------|
| Disposable stands for 4 units                 | 6   | VPA002 |
| Plastic recovery pipettes (Vivapore® 10   20) | 100 | VPA005 |
| 10 ml expansion reservoir (Vivapore® 10   20) | 10  | VPA006 |
| Plastic recovery pipettes (Vivapore® 5)       | 100 | VPA007 |

## ■ Ultrafiltration Membrane Filters

PES 146, CTA 145 and Hydrosart® 144



### Polyethersulfone (PES)

This is a general purpose membrane that provides excellent performance with most solutions when retentate recovery is of primary importance. Polyethersulfone membranes exhibit no hydrophobic or hydrophilic interactions and are usually preferred for their low fouling characteristics, exceptional flux and broad pH range.

### Cellulose Triacetate (CTA)

High hydrophilicity and exceptionally low non-specific binding characterize this membrane. Cast without any membrane

support that could trap or bind passing microsolute, these membranes are to be ideal for sample cleaning and protein removal and when high recovery of the filtrate solution is of primary importance.

### Hydrosart®

These membranes are also highly hydrophilic and are often preferred for their high protein recovery when processing very dilute solutions. Resistance to autoclaving, ease of cleaning and extended chemical resistance also characterize this type of membrane.

## □ Specifications

### Typical Performance for Polyethersulfone, Type 146

|                   |              |                            |
|-------------------|--------------|----------------------------|
| Thickness         | 120 µm       |                            |
| pH range          | 1–14         |                            |
| Water flux        | MWCO 10,000  | 0.2 ml/min/cm <sup>2</sup> |
| Protein retention | Cytochrome C | 95%                        |

### Specifications for Cellulose Triacetate, Type 145

|                   |              |                             |
|-------------------|--------------|-----------------------------|
| Thickness         | 120 µm       |                             |
| pH range          | 4–8          |                             |
| Water flux        | MWCO 10,000  | 0.11 ml/min/cm <sup>2</sup> |
| Protein retention | Cytochrome C | 90%                         |

### Specifications for Hydrosart®, Type 144

|                   |              |                             |
|-------------------|--------------|-----------------------------|
| Thickness         | 180 µm       |                             |
| pH range          | 1–13         |                             |
| Water flux        | MWCO 10,000  | 0.08 ml/min/cm <sup>2</sup> |
| Protein retention | Cytochrome C | 99%                         |

## □ Ordering Information

### Polyethersulfone Membrane Filters, Type 146

| Diameter in mm | MWCO           | Qty./Pkg. | Prod. No.       |
|----------------|----------------|-----------|-----------------|
| 47             | 1,000 daltons  | 10        | 14609--47-----D |
| 63             | 1,000 daltons  | 10        | 14609--63-----D |
| 76             | 1,000 daltons  | 10        | 14609--76-----D |
| 25             | 5,000 daltons  | 10        | 14629--25-----D |
| 47             | 5,000 daltons  | 10        | 14429--47-----D |
| 63             | 5,000 daltons  | 10        | 14629--63-----D |
| 76             | 5,000 daltons  | 10        | 14629--76-----D |
| 25             | 10,000 daltons | 10        | 14639--25-----D |

**Polyethersulfone Membrane Filters, Type 146**

| Diameter in mm | MWCO            | Qty./Pkg. | Prod. No.       |
|----------------|-----------------|-----------|-----------------|
| 63             | 10,000 daltons  | 10        | 14639--63-----D |
| 76             | 10,000 daltons  | 10        | 14639--76-----D |
| 150            | 10,000 daltons  | 10        | 14639-150-----D |
| 25             | 30,000 daltons  | 10        | 14659--25-----D |
| 63             | 30,000 daltons  | 10        | 14659--63-----D |
| 76             | 30,000 daltons  | 10        | 14659--76-----D |
| 25             | 50,000 daltons  | 10        | 14650--25-----D |
| 47             | 50,000 daltons  | 10        | 14650--47-----D |
| 76             | 50,000 daltons  | 10        | 14650--76-----D |
| 25             | 300,000 daltons | 10        | 14679--25-----D |
| 47             | 300,000 daltons | 10        | 14679--47-----D |
| 76             | 300,000 daltons | 10        | 14679--76-----D |

**Cellulose Triacetate Membrane Filters, Type 145**

| Diameter in mm | MWCO           | Qty./Pkg. | Prod. No.       |
|----------------|----------------|-----------|-----------------|
| 25             | 5,000 daltons  | 10        | 14529--25-----D |
| 47             | 5,000 daltons  | 10        | 14529--47-----D |
| 25             | 10,000 daltons | 10        | 14539--25-----D |
| 47             | 10,000 daltons | 10        | 14539--47-----D |
| 50             | 10,000 daltons | 10        | 14539--50-----D |
| 25             | 20,000 daltons | 10        | 14549--25-----D |
| 43             | 20,000 daltons | 10        | 14549--43-----D |
| 47             | 20,000 daltons | 10        | 14549--47-----D |
| 47             | 20,000 daltons | 100       | 14549--47-----N |
| 63             | 20,000 daltons | 10        | 14549--63-----D |
| 25             | 30,000 daltons | 10        | 14459--25-----D |
| 76             | 30,000 daltons | 10        | 14459--76-----D |

**Hydrosart® Membrane Filters, Type 144**

| Diameter in mm | MWCO           | Qty./Pkg. | Prod. No.       |
|----------------|----------------|-----------|-----------------|
| 25             | 5,000 daltons  | 10        | 14429--25-----D |
| 44             | 5,000 daltons  | 10        | 14429--44-----D |
| 63             | 5,000 daltons  | 10        | 14429--63-----D |
| 76             | 5,000 daltons  | 10        | 14429--76-----D |
| 25             | 10,000 daltons | 10        | 14439--25-----D |
| 47             | 10,000 daltons | 10        | 14439--47-----D |
| 63             | 10,000 daltons | 10        | 14439--63-----D |
| 76             | 10,000 daltons | 10        | 14439--76-----D |
| 25             | 30,000 daltons | 100       | 14459--25-----D |
| 47             | 30,000 daltons | 10        | 14459--47-----D |
| 63             | 30,000 daltons | 10        | 14459--63-----D |
| 76             | 30,000 daltons | 10        | 14459--76-----D |





## ■ Vivacon® 500

For DNA Sample Desalting and Concentration



### Reproducible DNA and Protein Sample Desalting and Concentration

Vivacon® 500 centrifugal concentrators offer the optimal solution for DNA and protein concentration and buffer exchange applications. For optimal performance with highly dilute samples, e.g. forensic samples, Vivacon® 500 incorporates the patented regenerated cellulose membrane Hydrosart®.

High recoveries and excellent reproducibilities are combined with convenience offered by the molecular weight cutoff printed on the individual Vivacon® 500 units.

As Vivacon® 500 can be respun after sample processing, this ensures complete concentrate recovery, which is especially important when working with low-sample concentrations.

### Vivacon® 500-PCR Grade

To use DNA amplification technologies, any traces of DNA originating from the equipment need to be eliminated.

Vivacon® 500-PCR Grade units are treated with ethylene oxide (ETO) in a validated process to deactivate all traces of DNA that might interfere with subsequent amplification procedures.

References: K. Shaw et al., Int. J. Legal Med. (2008) 122: 29–33

## □ Specifications

### Vivacon® 500

|                           |  |                      |
|---------------------------|--|----------------------|
| Concentrator capacity     | Fixed-angle rotor                      | 0.5 ml               |
| Dimensions                | Total length (concentration)           | 45 mm                |
|                           | Total length (backspin)                | 47.5 mm              |
|                           | Width                                  | 12.4 mm              |
|                           | Active membrane area                   | 0.32 cm <sup>2</sup> |
|                           | Hold-up volume of membrane and support | < 5 µl               |
|                           | Dead-stop volume (40° rotor)           | 5 µl                 |
| Materials of construction | Body                                   | Polycarbonate        |
|                           | Filtrate vessel                        | Polypropylene        |
|                           | Membrane                               | Hydrosart®           |

### Conversion Table for Hydrosart® MWCO to Nucleotide Cutoff

| Membrane          | MWCO    | Double-Stranded Nucleotide Cutoff (bp) |
|-------------------|---------|--|
| Hydrosart®        | 2 kDa   | > 10                                   |
| Hydrosart®        | 10 kDa  | > 30                                   |
| Hydrosart®        | 30 kDa  | > 50                                   |
| Hydrosart®        | 50 kDa  | > 300                                  |
| Hydrosart®        | 100 kDa | > 600                                  |
| Cellulose Acetate | 125 kDa | > 650                                  |

### ☐ Ordering Information

| <b>Vivacon® 500</b>  | <b>Qty./Pkg.</b> | <b>Prod. No.</b> |
|--|------------------|------------------|
| 2,000 MWCO   | 25               | VN01H91          |
| 2,000 MWCO   | 100              | VN01H92          |
| 10,000 MWCO  | 25               | VN01H01          |
| 10,000 MWCO  | 100              | VN01H02          |
| 30,000 MWCO  | 25               | VN01H21          |
| 30,000 MWCO  | 100              | VN01H22          |
| 50,000 MWCO  | 25               | VN01H31          |
| 50,000 MWCO  | 100              | VN01H32          |
| 100,000 MWCO   | 25               | VN01H41          |
| 100,000 MWCO   | 100              | VN01H42          |
| 125,000 MWCO   | 25               | VN01H81          |
| 125,000 MWCO   | 100              | VN01H82          |
| 125,000 MWCO   | 500              | VN01H83          |
| <b>Vivacon® 500</b>  | <b>Qty./Pkg.</b> | <b>Prod. No.</b> |
| Sample kit L<br>(4 units each of 2 kDa, 10 kDa,<br>30 kDa)   | 12               | VN01HL12         |
| Sample kit H<br>(4 units each of 30 kDa, 50 kDa,<br>100 kDa) | 12               | VN01HH12         |

| <b>Vivacon® 500-PCR Grade</b> | <b>Qty./Pkg.</b> | <b>Prod. No.</b> |
|-------------------------------|------------------|------------------|
| 30,000 MWCO                   | 25               | VN01H21ETO       |
| 30,000 MWCO                   | 100              | VN01H22ETO       |
| 50,000 MWCO                   | 25               | VN01H31ETO       |
| 50,000 MWCO                   | 100              | VN01H32ETO       |
| 100,000 MWCO                  | 25               | VN01H41ETO       |
| 100,000 MWCO                  | 100              | VN01H42ETO       |
| 125,000 MWCO                  | 25               | VN01H81ETO       |
| 125,000 MWCO                  | 100              | VN01H82ETO       |
| <b>Accessories</b>            | <b>Qty./Pkg.</b> | <b>Prod. No.</b> |
| Tubes                         | 100              | VNCT01           |

## Vivacon® 2

For DNA Sample Desalting and Concentration



### Reproducible DNA Sample Desalting and Concentration

Vivacon® 2 centrifugal concentrators offer the optimal solution for DNA and protein concentration and buffer exchange applications. For optimal performance with highly dilute samples, e.g. forensic samples, Vivacon® 2 incorporates the patented regenerated cellulose membrane Hydrosart®.

High recoveries and excellent reproducibilities are combined with the convenience provided by the volume graduation and molecular weight cutoff printed on the individual Vivacon® 2 units.

As Vivacon® 2 can be respun after sample processing, this ensures complete concentrate recovery, which is especially important when working with low-sample concentrations.

### Vivacon® 2-PCR Grade

Vivacon® 2-PCR Grade units are treated with ethylene oxide (ETO) in a validated process to deactivate all traces of DNA that might interfere with subsequent amplification procedures.

## Specifications

|                           |                                     |                      |
|---------------------------|-------------------------------------|----------------------|
| Concentrator capacity     | Fixed-angle rotor                   | 2 ml                 |
| Dimensions                | Total length (concentration)        | 125 mm               |
|                           | Total length (backspin)             | 115 mm               |
|                           | Width                               | 16 mm                |
|                           | Active membrane area                | 0.95 cm <sup>2</sup> |
|                           | Hold-up volume membrane and support | 10 µl                |
|                           | Dead-stop volume (25° rotor)        | 55 µl                |
| Materials of construction | Body                                | Polycarbonate        |
|                           | Filtrate vessel                     | Polypropylene        |
|                           | Backspin vial                       | Polypropylene        |
|                           | Concentrator cap                    | Polypropylene        |
|                           | Membrane                            | Hydrosart®           |

### Conversion Table for Hydrosart® MWCO to Nucleotide Cutoff

| Membrane          | MWCO    | Double-Stranded Nucleotide Cutoff (bp) |
|-------------------|---------|--|
| Hydrosart®        | 2 kDa   | > 10                                   |
| Hydrosart®        | 10 kDa  | > 30                                   |
| Hydrosart®        | 30 kDa  | > 50                                   |
| Hydrosart®        | 50 kDa  | > 300                                  |
| Hydrosart®        | 100 kDa | > 600                                  |
| Cellulose Acetate | 125 kDa | > 650                                  |

### ☐ Ordering Information

| <b>Vivacon® 2</b>           | <b>Qty./Pkg.</b> | <b>Prod. No.</b> |
|-----------------------------|------------------|------------------|
| 2,000 MWCO                  | 25               | VN02H91          |
| 2,000 MWCO                  | 100              | VN02H92          |
| 10,000 MWCO                 | 25               | VN02H01          |
| 10,000 MWCO                 | 100              | VN02H02          |
| 30,000 MWCO                 | 25               | VN02H21          |
| 30,000 MWCO                 | 100              | VN02H22          |
| 50,000 MWCO                 | 25               | VN02H31          |
| 50,000 MWCO                 | 100              | VN02H32          |
| 100,000 MWCO                | 25               | VN02H41          |
| 100,000 MWCO                | 100              | VN02H42          |
| 125,000 MWCO                | 25               | VN02H81          |
| 125,000 MWCO                | 100              | VN02H82          |
| <b>Vivacon® 2-PCR Grade</b> | <b>Qty./Pkg.</b> | <b>Prod. No.</b> |
| 30,000 MWCO                 | 25               | VN02H21ETO       |
| 30,000 MWCO                 | 100              | VN02H22ETO       |
| 50,000 MWCO                 | 25               | VN02H31ETO       |
| 50,000 MWCO                 | 100              | VN02H32ETO       |
| 100,000 MWCO                | 25               | VN02H41ETO       |
| 100,000 MWCO                | 100              | VN02H42ETO       |
| 125,000 MWCO                | 25               | VN02H81ETO       |
| 125,000 MWCO                | 100              | VN02H82ETO       |



## ■ Vivaspin® Endotest



Vivaspin® Endotest is a single-use ultrafiltration device for endotoxin concentration and removal of interfering substances from liquid samples for LAL testing. The Vivaspin® Endotest is certified as endotoxin-free and is available with a cellulose triacetate membrane that has a molecular weight cutoff (MWCO) of 20,000.

Due to its centrifugal design, Vivaspin® Endotest enables parallel testing to minimize hands-on time.

Vivaspin® Endotest can effectively be used in either swing-bucket or fixed-angle rotors that accept 50 ml conical bottom centrifuge tubes.

### □ Specifications

|                           |                                  |                      |
|---------------------------|----------------------------------|----------------------|
| Concentrator capacity     | Swing-bucket rotor               | 15.0 ml              |
|                           | Fixed-angle rotor                | 12.5 ml              |
| Dimensions                | Total length                     | 116 mm               |
|                           | Diameter                         | 30 mm                |
|                           | Active membrane area             | 3.9 cm <sup>2</sup>  |
|                           | Dead-stop volume                 | 30 µl                |
| Materials of construction | Non-recoverable (hold-up) volume | ≤ 10 µl              |
|                           | Body                             | Polycarbonate        |
|                           | Filtrate vessel                  | Polycarbonate        |
|                           | Concentrator cap                 | Polypropylene        |
|                           | Membrane                         | Cellulose triacetate |

### □ Ordering Information


| Vivaspin® Endotest | Qty./Pkg. | Prod. No. |
|--------------------|-----------|-----------|
| 20,000             | 12        | VS15RXETO |





## Vivapure®

### Ion Exchange Protein Purification Products



Chromatography gel beads (right) are shown on top of a membrane adsorber in this SEM. The membrane adsorber pores are more than 50 times larger than bead pores.

#### Fast and Easy-to-Use Spin Columns

Vivapure® Ion Exchange (IEX) spin columns are centrifugal devices that incorporate Sartobind® membrane adsorber technology as their chromatography matrix. Vivapure® IEX spin columns make protein purification as easy as filtration. The devices are ready to use and eliminate the risk of running dry. For many protein purification applications, they replace time-consuming and tedious column chromatography.

The rapid Vivapure® IEX 1-2-3 bind-wash-elute protocol especially lends itself to screening applications, where many different samples are processed in parallel.

#### The Sartobind® Membrane Adsorber Matrix

Sartobind® IEX membrane adsorbers are based on stabilized regenerated cellulose and exhibit a microporous structure with a pore size of  $>3\ \mu\text{m}$ , which is orders of magnitude larger than conventional chromatographic gel materials. This allows molecules to be transported to the ligands immobilized on the membrane adsorber by convective flow, resulting in exceptionally high flow rates.

By contrast, gel chromatography is slowed down due to diffusion limitations, as the molecules need to enter the small bead pores in order to be bound by the ligands. The porous membrane adsorber enables fast, reproducible and scalable protein purification.

#### Fast and Simple-to-Use Spin Columns

- Devices are ready to use
- They make protein purification as simple as filtration

#### Reproducible Results

- No column packing necessary as the devices are ready to use
- Membrane adsorber spin columns cannot crack or run dry

#### Centrifugal Devices


- Offer the possibility of working in parallel

#### Low Bed Volume

- Small membrane adsorber bed volumes allow working with lower buffer amounts, resulting in concentrated elution fractions

#### Scalable Product Range

- Process-scale modules are available with the same Sartobind® IEX membrane adsorber matrix



Fast and easy protein purification with Vivapure® spin columns



Vivapure® Mini-400 | 500 µl  
Binding capacities:  
1 mg to 4 mg



Vivapure® Maxi-19 | 20 ml  
Binding capacities:  
15 mg to 80 mg

## Specifications

### Available Formats

| Vivapure® IEX Products      | Application  |
|-----------------------------|--|
| Vivapure® Mini spin columns | <ul style="list-style-type: none"> <li>– Sample fractionation</li> <li>– Purification condition scouting</li> <li>– Small-scale purification</li> </ul>                                  |
| Vivapure® Maxi spin columns | <ul style="list-style-type: none"> <li>– Large-scale sample fractionation</li> <li>– One-step protein purification   concentration</li> <li>– Polishing of His-tagged protein</li> </ul> |

### Membrane Availability

| Functional Groups       | Ion Exchanger Type              |   |
|-------------------------|---------------------------------|---|
| Sulfonic acid (S)       | Strong acidic cation exchanger: | $\text{R-CH}_2\text{-SO}_3^-\text{Na}^+$              |
| Quaternary ammonium (Q) | Strong basic anion exchanger:   | $\text{R-CH}_2\text{-N}^+\text{-(CH}_3)_3\text{Cl}^-$ |
| Diethylamine (D)        | Weak basic anion exchanger:     | $\text{R-CH}_2\text{-NH}^+\text{-(CH}_2\text{H}_5)_2$ |

### Performance Characteristics

| Vivapure® Spin Columns | Protein Binding Capacity* (mg) | Max. Volume per Centrifuge Run Using a Swing-Out Rotor (ml) | Max. Volume per Centrifuge Using a Fixed-Angle Rotor Run (ml) |
|------------------------|--------------------------------|---|---|
| Vivapure® Mini H       | 4                              | 0.4   |   |
| Vivapure® Maxi H       | 60 to 80                       | 19  | 10.5  |

\* Actual yields depend on specific protein sample and selected pH and salt conditions.  
Yields established using 1 mg/ml BSA in 25 mM Tris/HCL pH 8.0 with Vivapure® Q & D spin columns and 1 mg/ml cytochrome c in 25 mM sodium acetate buffer pH 5.5 with Vivapure® S spin columns.

## Ordering Information

| Description  | Spin Columns | Centrifuge Tubes | Prod. No.   |
|--|--------------|------------------|-------------|
| <b>Vivapure® Mini Ion Exchange Spin Columns (up to 0.5 ml)</b> |              |                  |             |
| Vivapure® Mini S&Q H starter kit                               | 16           | 32               | VS-IX01SQ16 |
| Vivapure® D Mini H   | 24           | 48               | VS-IX01DH24 |
| Vivapure® Q Mini H   | 24           | 48               | VS-IX01QH24 |
| Vivapure® S Mini H   | 24           | 48               | VS-IX01SH24 |
| <b>Vivapure® Maxi Ion Exchange Spin Columns (up to 20 ml)</b>  |              |                  |             |
| Vivapure® D Maxi H   | 8            | 16               | VS-IX20DH08 |
| Vivapure® Q Maxi H   | 8            | 16               | VS-IX20QH08 |
| Vivapure® S Maxi H   | 8            | 16               | VS-IX20SH08 |

# ■ Vivapure® Miniprep and Maxiprep Purification Kits

For Fast Antibody and His-Tagged Protein Purification



## Rapid Purification with High Yields

Vivapure® Miniprep and Maxiprep kits are spin-column-based kits for fast and effective purification of His-tagged proteins and antibodies.

Spin columns have the advantage of speed over gravity drip columns and batch protocols.


With the patented FlowGo regulator, the sample residence time is extended to assure adequate sample binding to the resin. As a result, Vivapure® Miniprep and Maxiprep spin column kits combine the merits of spin columns and gravity drip columns. This results in rapid purification with up to 95% protein recovery and purity.




All spin columns can be conveniently used in a centrifuge. For processing large sample volumes, e.g. from diluted cell culture supernatants, Vivapure® Maxiprep spin columns can additionally be run with a peristaltic pump collar (VS-PPCSC).

Vivapure® Miniprep and Maxiprep columns come in a variety of different configurations for your convenience. They can be purchased as completely ready-to-use kits with buffers and ultrafiltration devices but just as well as stand-alone spin columns in small or large packages for frequent users.

## Vivapure® Miniprep



## Vivapure® Maxiprep



## □ Specifications

### Protein A and G for Antibody Purification

#### Protein A and G

| Miniprep                 | Centrifuge      |
|--------------------------|-----------------|
| Sample size              | 0.65 ml         |
| Typical binding capacity | 1 mg IgG/column |
| Number of reuses         | 3               |

#### Protein A&G

| Maxiprep                 | Centrifuge <sup>1</sup> |
|--------------------------|-------------------------|
| Sample size              | 20 ml                   |
| Typical binding capacity | 20 mg IgG/column        |
| Number of reuses         | 5                       |

<sup>1</sup> Use the peristaltic pump accessory (VS-PPCSC) for larger volumes

## □ Ordering Information

| Product Name                            | Qty./Pkg.           | Prod. No.   |
|---|---------------------|-------------|
| Vivapure® A Starter Pack*               | 2 Miniprep A units  | VS-ARSTPKA2 |
| Vivapure® Miniprep A Kit*               | 16 Miniprep A units | VS-ARAMINIK |
| Vivapure® Miniprep A Bulk Pack          | 48 Miniprep A units | VS-ARAMINIB |
| Vivapure® Miniprep A Kit*               | 4 Maxiprep A units  | VS-ARAMAXIK |
| Vivapure® Maxiprep A Bulk Pack          | 12 Maxiprep A units | VS-ARAMAXIB |
| Vivapure® A Buffer Pack                 |                     | VS-ARABUFPK |
| Vivapure® G Starter Pack*               | 2 Miniprep G units  | VS-ARSTPKG2 |
| Vivapure® Miniprep G Kit*               | 16 Miniprep G units | VS-ARGMINIK |
| Vivapure® Miniprep G Bulk Pack          | 48 Miniprep G units | VS-ARGMINIB |
| Vivapure® Maxiprep G Kit*               | 4 Maxiprep G units  | VS-ARGMAXIK |
| Vivapure® Maxiprep G Bulk Pack          | 12 Maxiprep G units | VS-ARGMAXIB |
| Vivapure® G Buffer Pack                 |                     | VS-ARGBUFPK |
| Sealing cap and peristaltic pump collar | 1                   | VS-PPCSC    |

\* including ultrafiltration concentrators and buffers

## □ Specifications

| <b>Protein MC<br/>Miniprep Kits</b> |         | <b>Centrifuge</b>  |
|-------------------------------------|---------|--------------------|
| Sample size                         | 0.65 ml |                    |
| Typical binding capacity            | 1 mg    | His-tagged protein |
| Number of reuses                    | 2       |                    |

| <b>Protein MC<br/>Maxiprep Kits</b> |       | <b>Centrifuge<sup>1</sup></b> |
|-------------------------------------|-------|-------------------------------|
| Sample size                         | 20 ml |                               |
| Typical binding capacity            | 10 mg | His-tagged protein            |
| Number of reuses                    | 2     |                               |

<sup>1</sup> Use the peristaltic pump accessory (VS-PPCSC) for larger volumes

## □ Ordering Information

| <b>Product Name</b>                                    | <b>Qty./Pkg.</b> | <b>Prod. No.</b> |
|--|------------------|------------------|
| Vivapure <sup>®</sup> Metal Chelate (MC) Starter Pack* | 4                | VS-MCST04        |
| Vivapure <sup>®</sup> Miniprep MC Kit*                 | 24               | VS-MCMINI24      |
| Vivapure <sup>®</sup> Miniprep MC Bulk Pack            | 72               | VS-MCMINIB       |
| Vivapure <sup>®</sup> Maxiprep MC Kit*                 | 8                | VS-MCMAXIK       |
| Vivapure <sup>®</sup> Maxiprep MC Bulk Pack            | 24               | VS-MCMAXIB       |
| Vivapure <sup>®</sup> Metal Chelate Buffer Pack        |                  | VS-MCBUFPK       |

\* including ultrafiltration concentrators and buffers





## ■ Vivapure<sup>®</sup> Virus Purification and Concentration Kits



Recombinant virus vectors are the preferred method for a wide range of gene delivery applications. Especially adenovirus type 5 and VSV-G pseudotyped lentivirus are two frequently utilized viral vectors for in vitro and in vivo applications.

### **Recombinant Adenovirus Vectors**

Recombinant adenovirus vectors are versatile tools in research and therapeutic applications for gene transfer and protein expression in cell lines that have low transfection efficiency with liposomes. After entering cells, the virus remains epichromosomal – i.e., does not integrate into the host chromosome, leaving the host genome unaffected. The delivery of RNAi into cells is becoming a major application for adenovirus vectors.

### **Lentivirus Vectors**

Lentivirus vectors are frequently used in gene transfer studies, due to their ability of gene transfer and integration into dividing and non-dividing cells. The pseudotyped envelope with vesicular stomatitis virus envelope G (VSV-G) protein broadens their target cell range. Lentiviral vectors have been shown to deliver genes into cell types (e.g. neurons, lymphocytes and macrophages) which other retrovirus vectors could not be used for. The lentivirus vector is increasingly used to integrate siRNA efficiently in a wide variety of cell lines and primary cells, both in vitro and in vivo.

### **Rapid Virus Purification by Membrane Chromatography**

The Sartobind<sup>®</sup> ion exchange membrane adsorber technology used in AdenoPACK and LentiSELECT is unique in its capability to efficiently and rapidly capture and recover large virus particles. Compared with chromatography media, membrane adsorbers provide large 3,000 nm pores, allowing unrestricted access and recovery of virus from the charged adsorber surface. Convective flow through the syringe filter devices provides high-speed separations not possible with traditional chromatography, cesium chloride density gradients and ultracentrifugation methods. Sartorius membrane adsorbers with porous matrices, high capacities, low differential pressures, high flow rates and low unspecific adsorption show excellent performance in small-scale virus purification. In addition, these syringe filter devices are scalable and comply with cGMP requirements for large-volume, high-performance separation, reducing final process time ten-fold.

## ■ Adenovirus Purification with Vivapure® AdenoPACK Kits

### AdenoPACK 20|100|500

The AdenoPACK adenovirus purification and concentration kits offer researchers who need to recover up to  $3 \times 10^{13}$  purified recombinant adenovirus particles for in vitro transfection a fast, safe and easy-to-use solution. The kits include all reagents and devices necessary for clarification, purification and concentration of adenovirus type 5 from HEK293 cell cultures – all within just two hours. These straightforward kits replace time-consuming and labor-intensive 48-hour CsCl density gradients.

AdenoPACK kits are offered as AdenoPACK 20, AdenoPACK 100 and AdenoPACK 500 for the purification and concentration of adenovirus type 5 from 20 ml to 500 ml cell cultures, resulting in  $1 \times 10^{11}$  to  $3 \times 10^{13}$  purified viral particles. For each sample volume, the most convenient handling method is provided for ultimate convenience.

To this end, preparations using AdenoPACK 20 are supplied in a spin column format for centrifuges. AdenoPACK 100 is a manually operated kit in a syringe filter format\* and AdenoPACK 500 is a pump-driven kit.

\*Vivapure® AdenoPACK 100 can be alternatively be operated with a laboratory pump or an infusion pump, for which protocols are provided on our web page at [www.sartorius-stedim.com](http://www.sartorius-stedim.com). Additionally, the tubes and adaptors needed for these operating modes can be ordered.

### AdenoPACK Advantages

#### Fast and Easy Virus Purification

- Purification completed in just 2 hours
- Convenient, over 10× faster alternative to CsCl density gradient

#### Quantitative Yields

- In contrast to CsCl density gradient, the complete cell culture is used for virus purification and not only the viral pellet

#### Flexible Product Range

- Applicable from initial construct screening to large-scale virus production

#### Complete Kit

- Including filtration devices, AdenoPACK units for virus purification, Vivaspin® and all buffers

#### Low Endotoxin Levels

- High cell viability and infection rates due to endotoxin levels of < 0.025 EU/ml

## □ Specifications

### Adenovirus Purification Kit Specifications

| Product                     | AdenoPACK 20                                   | AdenoPACK 100                      | AdenoPACK 500                      |
|-----------------------------|--|------------------------------------|------------------------------------|
| Sample size                 | 20 ml cell culture                             | 20 ml to 200 ml of cell culture    | 500 ml of cell culture             |
| Number of purifications     | 6 × 20 ml                                      | 2 × 20 ml to 60 ml<br>1 × 200 ml   | 1 × 500 ml                         |
| Virus particles (VP) per ml | Typically up to $1 \times 10^{11}$ – $10^{12}$ | Typically up to $1 \times 10^{13}$ | Typically up to $3 \times 10^{13}$ |
| VP/IU                       | 50 to 100                                      | 20 to 50                           | 20 to 50                           |
| Processing time             | Typically one hour                             | Typically two hours                |                                    |
| Endotoxin level             | <0.025 EU/ml                                   | <0.025 EU/ml                       | <0.025 EU/ml                       |



## □ Ordering Information

### **Vivapure® AdenoPACK 20**

|                        |            |
|------------------------|------------|
| Vivapure® AdenoPACK 20 | VS-AVPQ020 |
|------------------------|------------|

|                            |            |
|----------------------------|------------|
| Vivapure® AdenoPACK 20 RT* | VS-AVPQ022 |
|----------------------------|------------|

\* AdenoPACK 20 RT does not contain Benzonase®\*\*



### **Vivapure® AdenoPACK 100**

|                         |            |
|-------------------------|------------|
| Vivapure® AdenoPACK 100 | VS-AVPQ101 |
|-------------------------|------------|

|                             |            |
|-----------------------------|------------|
| Vivapure® AdenoPACK 100 RT* | VS-AVPQ102 |
|-----------------------------|------------|

### **AdenoPACK 100 Accessories**

|            |  |
|------------|--|
| VS-AVPA001 | Pump tubing set for<br>Vivapure® AdenoPACK 100 |
|------------|--|

\* AdenoPACK 100 RT does not contain Benzonase®\*\*



### **Vivapure® AdenoPACK 500**

|                         |            |
|-------------------------|------------|
| Vivapure® AdenoPACK 500 | VS-AVPQ501 |
|-------------------------|------------|

|                             |            |
|-----------------------------|------------|
| Vivapure® AdenoPACK 500 RT* | VS-AVPQ502 |
|-----------------------------|------------|

\* AdenoPACK 500 RT does not contain Benzonase®\*\*

\*\* Benzonase® is a registered trademark of Merck

## ■ Lentivirus Purification with the Vivapure® LentiSELECT Kit

### LentiSELECT 40|500|1000

The LentiSELECT kits for lentivirus purification and concentration offer researchers who need to recover up to  $5 \times 10^9$  infective lentivirus particles per ml for in vitro transfection or animal studies a fast and easy-to-use solution.

These straightforward kits replace time-consuming ultracentrifugation protocols, which typically take approximately one day for large sample volumes. Vivapure® LentiSELECT thus reduces purification time to just a few hours.

LentiSELECT kits are offered as LentiSELECT 40, LentiSELECT 500 and LentiSELECT 1000 for the purification and concentration of VSV-G pseudotyped lentivirus from 40 ml to 1,000 ml cell cultures, resulting in  $8 \times 10^8$  to  $1 \times 10^{10}$  purified infective particles. The most convenient handling method is provided for each sample volume. To this end, 40 ml sample volumes are processed manually with LentiSELECT 40, while LentiSELECT 500 and 1000 are pump-driven kits.

### LentiSELECT Advantages

#### Fast and Easy Virus Purification

- Purification completed in less than one to six hours, depending on sample volume
- Kit is as easy to use as filtration

#### No Need for Expensive Instruments

- Lentivirus purification with LentiSELECT is independent of equipment, such as ultracentrifuges

#### High Virus Purity

- Achieve pure virus based on a chromatographic method for your experiments instead of a crude and variable cell culture supernatant pellet

#### Optimal for Multiple Virus Construct Screening

With LentiSELECT 40, four purification runs can be conducted in parallel with one kit

#### Complete Kits

- Including LentiSELECT units for virus purification, Vivaspin® units for concentration| buffer exchange and all buffers and syringes necessary

#### Low Endotoxin Levels

- High cell viability and infection rates due to endotoxin levels of  $< 0.025$  EU/ml

## □ Specifications

| Product                     | LentiSELECT 40                  | LentiSELECT 500                     | LentiSELECT 1000                     |
|-----------------------------|---------------------------------|-------------------------------------|--------------------------------------|
| Sample size                 | 40 ml cell culture              | 500 ml of cell culture              | 1,000 ml of cell culture             |
| Number of purifications     | 4 × 40 ml                       | 1 × 500 ml                          | 1 × 1,000 ml                         |
| Virus particles (VP) per ml | Typically up to $3 \times 10^9$ | Typically up to $2-5 \times 10^9$ * | Typically up to $4-6 \times 10^{13}$ |
| VP/IU                       | 5 to 15                         | 5 to 15                             | 20 to 50                             |
| Processing time             | Typically up to 45 min          | Typically up to 3 hours             | Typically up to 6 hours              |
| Endotoxin level             | $< 0.025$ EU/ml                 | $< 0.025$ EU/ml                     | $< 0.025$ EU/ml                      |



### ☐ Ordering Information

**Vivapure® LentiSELECT 40**

Vivapure® LentiSELECT 40

VS-LVPQ040

**Vivapure® LentiSELECT 500**

Vivapure® LentiSELECT 500

VS-LVPQ500

**Vivapure® LentiSELECT 1000**

Vivapure® LentiSELECT 1000

VS-LVPQ1000



A gloved hand is shown pouring a yellow liquid from a beaker into a funnel. The background is a light, neutral color.

## ■ Filtration Devices

### Table of Contents

60 Minisart® Syringe Filters

76 Sartolab® P20 and Sartolab® P20 Plus

79 Sartolab® 150v

81 Sartolab® RF | BT

84 Chemical Compatibility





Syringe filters are used for many routine preparation steps in laboratories all over the world. They are convenient, ready-to-use disposables for sterile filtration of liquids and removal of particles from solutions and gases. Depending on the reagents filtered, syringe filters have to fulfill certain requirements to best serve customer's application. Sartorius offers Minisart® syringe filters and filters optimized for a wide range of relatively large volumes. The filters are clean and safe as they are virtually free of leachables and extractables and reliably remove particles and microorganisms with no leakage. If you need to rely on the quality of your filtrate – whether it needs to be sterile prior to injection or particulate-free before analysis – field-proven, high-quality Sartorius filter syringes are the No. 1 choice for reliable, convenient preparation steps.

#### **Our Product Range**

Minisart® syringe filters feature a housing made of medical-grade acrylic (MBS), which makes them the perfect choice for sterile filtration and clarification of additives, buffers, reagents, drugs and gases. Their effective filtration area of 6.2 cm<sup>2</sup> for liquids is the largest among premium syringe filters. The combination of the MBS housing material and surfactant-free cellulose acetate, PES or PTFE membranes add comma provides the lowest non-specific binding. As a result, Minisart® delivers excellent performance in filtering medical drugs or gases for sterilization or particle removal prior to use.

Minisart® syringe filters with a polypropylene housing are optimized for filtration prior to analysis and withstand even harsh solvents and chemicals. Due to the typical range of volumes from less than 1 ml to 100 ml, these filters are available in three different diameters with an effective filtration area of 0.07 cm<sup>2</sup>, 1.7 cm<sup>2</sup> and 4.8 cm<sup>2</sup> for the optimum filter capacity and the lowest hold-up volume.

Sartolab® vacuum filtration devices with 0.1 µm and 0.22 µm PES membranes for convenient filtration of 150 ml up to 1 L are ready to use and sterile. Sartolab® RF is a complete system that includes a receiver flask. Sartolab® BT Sterile is a bottle top filter without a receiver flasks. This enables customers to use a receiver bottle of their choice and to even expand filtration capacity, depending on the particle load of the filtered liquid by filling more than one receiver flask. Sartolab® 150V is a disposable vacuum filter with a pleated 0.22 µm PES membrane, which is suitable for up to 15 L of liquid.

Sartolab® pressure filtration devices with a 0.2 µm SFCA or PES membrane are available with or without a glass fiber prefilter, depending on your needs. Sartolab P20 is designed for up to 3 L volumes and can also be used in-line. The polycarbonate housing and membrane components are ideal for filtering liquid medical drugs. The glass fiber prefilter types are ideal for filtering environmental samples that have a high particle load prior to analyzing such samples.

### Typical Applications for Filtration Devices

- Sterile filtration of liquids and gases with virtually no effect on the ingredients
- Particle removal from liquids and gases prior to downstream processes
- Venting of vials, bottles, containers, bags and bioreactors and fermenters
- Removal of precipitates and coagulates from medical solutions prior to use



## ■ Minisart® Syringe Filters

### Removal of Particles and Microorganisms from Liquids and Gases

Sartorius offers Minisart® syringe filters for a wide range of applications. The filters are clean and safe as they are virtually free of leachables and extractables. In addition, they reliably remove particles and microorganisms – without any leakage. Minisart® with a polypropylene housing is optimized for filtration prior to analysis and withstands even harsh solvents and chemicals. Minisart® with a housing made of medical-grade acrylic (MBS) is the perfect choice for sterile filtration and clarification of additives, buffers, reagents, drugs and gases.

### ■ A Full Range of Filters Dedicated for Various Filtration Applications



#### Sample Preparation HPLC | UHPLC | Analytics

Elimination of particles from your samples prior to HPLC or other chromatographic analysis is essential in order to maintain the integrity of your chromatography column and to maximize its operating lifetime. Minisart® syringe filters optimized for sample preparation consist of a

polypropylene housing and membrane components featuring maximum chemical compatibility and minimum extractables to ensure excellent results. Due to the typical range of volumes from less than 1 ml to 100 ml, these filters are available in three different diameters with an effective filtration area of 0.07 cm<sup>2</sup>, 1.7 cm<sup>2</sup> and 4.8 cm<sup>2</sup>. For a selection guide, please see page 61.



#### Filtration of Aqueous Liquids Clarification | Sterile Filtration

For clarification and sterilization of liquids, filtration is the optimal method. It removes all microorganisms and particles reliably, without any effects on the ingredients due to adsorption or decomposition. For optimal results, Minisart® syringe filters with an MBS housing provide a choice of

membranes with pore sizes ranging from 0.1 µm to 5 µm for high flow rates and the lowest adsorption characteristics. The effective filtration area of 6.2 cm<sup>2</sup> for the fastest filtration is the largest among premium syringe filters available, and the MBS housing is color-coded for easy pore size identification. For a list of the types offered, please see page 65.



#### Medical Use and Special Sterile Venting Applications

Minisart® syringe filters are ideal for clarification of liquids laden with particles, e.g. for preparation of pharmaceuticals or infusion solutions. For sterilization and removal of particles from air and other gases, syringe filters are optimal for

sterile venting of containers, bioreactors, fermenters and tubing systems in medical devices. Many Minisart® syringe filters have the CE mark of conformity (European directive) and are available with a wide choice of membranes, connectors and housing materials. For a list of the types offered, please see page 68.

## Minisart® Selection Guide





Please refer to Minisart® RC, NY or SRP for the highest chemical compatibility on page 63.

Please refer to Minisart® NML or Minisart® High Flow on page 66.

Several additional Minisart® units for venting and special purposes are shown on page 70.

| Sample Composition | Aqueous  |   | Aqueous   Solvents                            |                                     |   |
|--------------------|--|---|---|-------------------------------------|---|
|                    | ✓<br><b>All Aqueous Solutions</b><br>Buffers, Protein Analysis | ✓<br><b>All Aqueous Solutions</b><br>Tissue Culture Media | ✓<br>Aqueous   Solvent<br>Mixtures   Solvents | ✓<br>Solvent Mixtures  <br>Solvents | ✓<br>Solvents   Gases  <br>Acids   Bases    |
|                    | ✓<br><b>SFCA</b><br>Surfactant-Free<br>Cellulose Acetate       | ✓<br><b>PES</b><br>Polyethersulfone                       | ✓<br><b>RC</b><br>Regenerated<br>Cellulose    | ✓<br><b>NY</b><br>Polyamide, Nylon  | ✓<br><b>PTFE</b><br>Polytetrafluoroethylene |
|                    | Hydrophilic  |   |   |                                     | Hydrophobic                                 |

| Pore Sizes | Sterilization   |  | Sample Preparation   Clarification   Particle Removal            |                               |                               |  |  | Prefiltration   |
|------------|---|--|--|-------------------------------|-------------------------------|--|--|---|
|            | ✓<br>Small<br>Bacteria<br>Mycoplasma<br>Colloids<br>> 0.1 | ✓<br>UHPLC, etc.<br>(Columns<br>< 3 µm<br>Particles)<br>Bacteria | ✓<br>HPLC, etc.<br>(Columns<br>> 3 µm<br>Particles)<br>Particles | ✓<br>Particles<br>Yeast Cells | ✓<br>Particles<br>Yeast Cells | ✓<br>Particles<br>Yeast Cells<br>Platelets | ✓<br>Large<br>Particles<br>Grit, Cells | ✓<br>Glass Prefilter<br>Glass+Membrane<br>Highly<br>Particle-laden<br>Samples |
|            | ✓<br><b>0.1 µm</b>  | ✓<br><b>0.2 µm</b>   | ✓<br><b>0.45 µm</b>  | ✓<br><b>0.65 µm</b>           | ✓<br><b>0.8 µm</b>            | ✓<br><b>1.2 µm</b>                         | ✓<br><b>5 µm</b>                       | ✓<br><b>GF (Glass Fiber)</b>  |

| Sample Volumes |  |  |  |  |
|----------------|---|---|--|---|
|                | ✓<br>1 ml to 200 ml   | ✓<br>1 ml to 100 ml   | ✓<br>0.5 ml to 15 ml   | ✓<br>0.05 ml to 1 ml  |
|                | ✓<br><b>28 mm for up to 200 ml</b>  | ✓<br><b>25 mm for up to 100 ml</b>  | ✓<br><b>15 mm for up to 15 ml</b>  | ✓<br><b>4 mm for up to 1 ml</b>   |

4 mm packages are color-coded

Male Spike Outlet

Male Luer Slip Outlet

Minisart® RC 15 mm

Minisart® NY 15 mm

Minisart® SRP 15 mm

Minisart® RC 25 mm

Minisart® NY 25 mm

Minisart® SRP 25 mm

## Sample Preparation for Analytical Procedures

### Reliable Removal of Particles and Microorganisms from Liquids and Gases




Particle removal by filtration before analysis substantially increases the lifetime of your columns. Minisart® RC is optimized for aqueous liquids and solvents and is compatible with DMSO, other amides, ketones, esters and ethers. Minisart® NY is exceptionally pure compared with other common polyamide (nylon) filters and competitor products. For this product raw materials are used which do not interfere with standard analytical methods.

Our coating-free hydrophobic PTFE membrane used in Minisart® SRP is suitable for venting applications and leachables-free clarification of exceptionally harsh chemicals.

#### Minisart® Features

- Low adsorption of analytes
- Maximum chemical compatibility
- Minimum extractables or leachables
- Bidirectional use possible
- 100% optical integrity test

### HPLC Certification



### HPLC Procedure

Column C18: 250 × 4.6 mm, Flow Rate: 1 ml/min, Wavelength: 220 nm

HPLC Injection Volume: 20 µl, Analysis Time: 65 min, Temperature: 40°C, Mobile Phases:

A) Acetonitrile | B) Water, Gradient: Hold 60% A for 10 min, 60% to 95% A in 20 min, 95% to 100% A in 35 min

## Sample Preparation for Chromatography

☐ Ordering Information

**Dia. in mm | EFA<sup>1</sup> Membrane Housing Pore Size Connector Outlet Color | Printing Sterile\* Qty./Pkg. Order No.**

**Minisart® RC (Regenerated Cellulose)**

|       |    |    |         |                |                |     |     |               |
|-------|----|----|---------|----------------|----------------|-----|-----|---------------|
| 25 mm | RC | PP | 0.2 µm  | Male Luer Slip | White, Printed | Yes | 50  | 17764-----ACK |
| 25 mm | RC | PP | 0.2 µm  | Male Luer Slip | White, Printed | No  | 50  | 17764-----K   |
| 25 mm | RC | PP | 0.2 µm  | Male Luer Slip | White, Printed | No  | 200 | 17764-----S   |
| 25 mm | RC | PP | 0.2 µm  | Male Luer Slip | White, Printed | No  | 500 | 17764-----Q   |
| 25 mm | RC | PP | 0.45 µm | Male Luer Slip | White, Printed | No  | 50  | 17765-----K   |
| 25 mm | RC | PP | 0.45 µm | Male Luer Slip | White, Printed | No  | 200 | 17765-----S   |
| 25 mm | RC | PP | 0.45 µm | Male Luer Slip | White, Printed | No  | 500 | 17765-----Q   |
| 15 mm | RC | PP | 0.2 µm  | Male Luer Slip | White, Printed | Yes | 50  | 17761-----ACK |
| 15 mm | RC | PP | 0.2 µm  | Male Luer Slip | White, Printed | No  | 50  | 17761-----K   |
| 15 mm | RC | PP | 0.2 µm  | Male Luer Slip | White, Printed | No  | 500 | 17761-----Q   |
| 15 mm | RC | PP | 0.45 µm | Male Luer Slip | White, Printed | No  | 50  | 17762-----K   |
| 15 mm | RC | PP | 0.45 µm | Male Luer Slip | White, Printed | No  | 500 | 17762-----Q   |
| 4 mm  | RC | PP | 0.2 µm  | Male Luer Slip | Blue Tray      | No  | 50  | 17821-----K   |
| 4 mm  | RC | PP | 0.2 µm  | Male Luer Slip | Blue Tray      | No  | 500 | 17821-----Q   |
| 4 mm  | RC | PP | 0.45 µm | Male Luer Slip | Yellow Tray    | No  | 50  | 17822-----K   |
| 4 mm  | RC | PP | 0.45 µm | Male Luer Slip | Yellow Tray    | No  | 500 | 17822-----Q   |

**Minisart® SRP (Hydrophobic PTFE)**

|       |      |    |         |                |                |     |     |               |
|-------|------|----|---------|----------------|----------------|-----|-----|---------------|
| 25 mm | PTFE | PP | 0.2 µm  | Male Luer Slip | White, Printed | Yes | 50  | 17575-----ACK |
| 25 mm | PTFE | PP | 0.2 µm  | Male Luer Slip | White, Printed | No  | 50  | 17575-----K   |
| 25 mm | PTFE | PP | 0.2 µm  | Male Luer Slip | White, Printed | No  | 200 | 17575-----S   |
| 25 mm | PTFE | PP | 0.2 µm  | Male Luer Slip | White, Printed | No  | 500 | 17575-----Q   |
| 25 mm | PTFE | PP | 0.45 µm | Male Luer Slip | White, Printed | No  | 50  | 17576-----K   |
| 25 mm | PTFE | PP | 0.45 µm | Male Luer Slip | White, Printed | No  | 200 | 17576-----S   |
| 25 mm | PTFE | PP | 0.45 µm | Male Luer Slip | White, Printed | No  | 500 | 17576-----Q   |
| 15 mm | PTFE | PP | 0.2 µm  | Male Spike     | White, Printed | No  | 50  | 17558-----K   |
| 15 mm | PTFE | PP | 0.2 µm  | Male Spike     | White, Printed | No  | 500 | 17558-----Q   |
| 15 mm | PTFE | PP | 0.2 µm  | Male Luer Slip | White, Printed | Yes | 50  | 17573-----ACK |
| 15 mm | PTFE | PP | 0.2 µm  | Male Luer Slip | White, Printed | No  | 50  | 17573-----K   |
| 15 mm | PTFE | PP | 0.2 µm  | Male Luer Slip | White, Printed | No  | 500 | 17573-----Q   |
| 15 mm | PTFE | PP | 0.45 µm | Male Spike     | White, Printed | No  | 50  | 17559-----K   |
| 15 mm | PTFE | PP | 0.45 µm | Male Spike     | White, Printed | No  | 500 | 17559-----Q   |
| 15 mm | PTFE | PP | 0.45 µm | Male Luer Slip | White, Printed | No  | 50  | 17574-----K   |
| 15 mm | PTFE | PP | 0.45 µm | Male Luer Slip | White, Printed | No  | 500 | 17574-----Q   |
| 4 mm  | PTFE | PP | 0.2 µm  | Male Luer Slip | Blue Tray      | No  | 500 | 17844-----Q   |
| 4 mm  | PTFE | PP | 0.45 µm | Male Luer Slip | Yellow Tray    | No  | 50  | 17820-----K   |
| 4 mm  | PTFE | PP | 0.45 µm | Male Luer Slip | Yellow Tray    | No  | 500 | 17820-----Q   |



**Dia. in mm | EFA<sup>1</sup> Membrane Housing Pore Size Connector Outlet Color | Printing Sterile\* Qty./Pkg. Order No.**

**Minisart® NY (Nylon) and NY25 Plus (Glass Fiber 0.7 µm<sup>2</sup> + Nylon)**

|       |          |    |         |                |                |     |     |               |
|-------|----------|----|---------|----------------|----------------|-----|-----|---------------|
| 25 mm | Nylon    | PP | 0.2 µm  | Male Luer Slip | White, Printed | Yes | 50  | 17845-----ACK |
| 25 mm | Nylon    | PP | 0.2 µm  | Male Luer Slip | White, Printed | No  | 500 | 17845-----Q   |
| 25 mm | Nylon    | PP | 0.45 µm | Male Luer Slip | White, Printed | Yes | 50  | 17846-----ACK |
| 25 mm | Nylon    | PP | 0.45 µm | Male Luer Slip | White, Printed | No  | 500 | 17846-----Q   |
| 15 mm | Nylon    | PP | 0.2 µm  | Male Luer Slip | White, Printed | No  | 50  | 1776B-----K   |
| 15 mm | Nylon    | PP | 0.2 µm  | Male Luer Slip | White, Printed | No  | 500 | 1776B-----Q   |
| 15 mm | Nylon    | PP | 0.45 µm | Male Luer Slip | White, Printed | No  | 50  | 1776C-----K   |
| 15 mm | Nylon    | PP | 0.45 µm | Male Luer Slip | White, Printed | No  | 500 | 1776C-----Q   |
| 25 mm | GF+Nylon | PP | 0.2 µm  | Male Luer Slip | White, Printed | No  | 50  | 1784B-----K   |
| 25 mm | GF+Nylon | PP | 0.2 µm  | Male Luer Slip | White, Printed | No  | 500 | 1784B-----Q   |
| 25 mm | GF+Nylon | PP | 0.45 µm | Male Luer Slip | White, Printed | No  | 50  | 1784C-----K   |
| 25 mm | GF+Nylon | PP | 0.45 µm | Male Luer Slip | White, Printed | No  | 500 | 1784C-----Q   |

\* Sterile Minisart® syringe filters are individually packaged. If not stated otherwise, Minisart® units have been sterilized by ethylene oxide.

Non-sterilized Minisart® units: RC, PTFE and nylon can be sterilized by autoclaving at 121°C for 30 min/or by using ethylene oxide (EO).

<sup>1</sup> Diameter of EFA – Effective Filtration Area

<sup>2</sup> 0.7 µm = GF particle retention ≠ pore size!

For technical product specifications, please see page 72.





Minisart® High Flow with PES



Minisart® NML with SFCA

## ■ Filtration of Aqueous Liquids – Clarification | Sterilization

### Filtration Is the Optimal Method for Clarification and Sterilization of Liquids




Sterilization by filtration is the fastest method for removal of bacterial cells from liquids, while minimizing the effects on ingredients. Minisart® NML with surfactant-free cellulose acetate (SFCA) is the best choice for all aqueous solutions with a pH of 4 to 8. It combines fast flow rates with an exceptionally pure and leachables-free membrane and is available in many different pore sizes – also for the removal of larger particles. Minisart® High Flow with polyethersulfone (PES) is optimal for delivering the highest flow rates and for a broad pH compatibility range from 1 to 13. Due to the asymmetric membrane structure, the PES surface virtually behaves like a prefilter. Both Minisart® types – NML and High Flow – can be sterilized by

ethylene oxide (EO) or gamma irradiation and are suitable for gluing and assembling on tubing. In addition, Minisart® NML has the CE mark of conformity.

#### Minisart® Features

- The largest effective filtration area (EFA) of 6.2 cm<sup>2</sup>
- The lowest adsorption
- Superior flow rate
- High total throughput
- Low hold-up volume
- Minimum extractables
- PVC-free
- Gamma-irradiated or EO-sterilized
- Bidirectional use possible
- 100% optical integrity test

### Water Flow Rates at 1 bar (ml/min); 15.4 psi; 0.45 µm Hydrophilic Membranes



28 mm EFA  
33 mm housing diameter

## Preparation of Aqueous Liquids

## □ Ordering Information

| Dia. in<br>mm   EFA <sup>1</sup>                                | Membrane | Housing | Pore<br>Size | Connector<br>Outlet | Color  <br>Printing | Sterile* | Qty./<br>Pkg. | Order No.     | CE Marking |
|---|----------|---------|--------------|---------------------|---------------------|----------|---------------|---------------|------------|
| <b>Minisart® High Flow</b> (PES – Polyethersulfone)             |          |         |              |                     |                     |          |               |               |            |
| 28 mm   | PES      | MBS     | 0.1 µm       | Male Luer Lock      | Dark Red            | Yes      | 50            | 16553-----K   |            |
| 28 mm   | PES      | MBS     | 0.22 µm      | Male Luer Lock      | Royal Blue          | Yes#     | 50            | 16532-----GUK |            |
| 28 mm   | PES      | MBS     | 0.22 µm      | Male Luer Lock      | Royal Blue          | Yes      | 50            | 16532-----K   |            |
| 28 mm   | PES      | MBS     | 0.22 µm      | Male Luer Slip      | Royal Blue          | Yes      | 50            | 16541-----K   |            |
| 28 mm   | PES      | MBS     | 0.22 µm      | Male Luer Lock      | Royal Blue          | No       | 500           | 16532-----Q   |            |
| 28 mm   | PES      | MBS     | 0.22 µm      | Male Luer Slip      | Royal Blue          | No       | 500           | 16541-----Q   |            |
| 28 mm   | PES      | MBS     | 0.45 µm      | Male Luer Lock      | Amber               | Yes      | 50            | 16537-----K   |            |
| 28 mm   | PES      | MBS     | 0.45 µm      | Male Luer Lock      | Amber               | No       | 500           | 16537-----Q   |            |
| 28 mm   | PES      | MBS     | 0.45 µm      | Male Luer Slip      | Amber               | Yes#     | 50            | 16533-----GUK |            |
| 28 mm   | PES      | MBS     | 0.45 µm      | Male Luer Slip      | Amber               | Yes      | 50            | 16533-----K   |            |
| 28 mm   | PES      | MBS     | 0.45 µm      | Male Luer Slip      | Amber               | No       | 500           | 16533-----Q   |            |
| <b>Minisart® NML</b> (SFCA – Surfactant-free Cellulose Acetate) |          |         |              |                     |                     |          |               |               |            |
| 28 mm   | SFCA     | MBS     | 0.2 µm       | Male Luer Lock      | Blue                | Yes      | 50            | 16534-----K   | CE-marked  |
| 28 mm   | SFCA     | MBS     | 0.2 µm       | Male Luer Lock      | Blue                | Yes#     | 50            | 16534-----GUK | CE-marked  |
| 28 mm   | SFCA     | MBS     | 0.2 µm       | Male Luer Lock      | Blue                | No       | 500           | 16534-----Q   | CE-marked  |
| 28 mm   | SFCA     | MBS     | 0.2 µm       | Male Luer Slip      | Blue                | Yes      | 50            | 17597-----K   | CE-marked  |
| 28 mm   | SFCA     | MBS     | 0.2 µm       | Male Luer Slip      | Blue                | No       | 500           | 17597-----Q   | CE-marked  |
| 28 mm   | SFCA     | MBS     | 0.45 µm      | Male Luer Lock      | Yellow              | Yes      | 50            | 16555-----K   | CE-marked  |
| 28 mm   | SFCA     | MBS     | 0.45 µm      | Male Luer Lock      | Yellow              | Yes#     | 50            | 16555-----GUK | CE-marked  |
| 28 mm   | SFCA     | MBS     | 0.45 µm      | Male Luer Lock      | Yellow              | No       | 500           | 16555-----Q   | CE-marked  |
| 28 mm   | SFCA     | MBS     | 0.45 µm      | Male Luer Slip      | Yellow              | Yes      | 50            | 17598-----K   | CE-marked  |
| 28 mm   | SFCA     | MBS     | 0.45 µm      | Male Luer Slip      | Yellow              | No       | 500           | 17598-----Q   | CE-marked  |
| 28 mm   | SFCA     | MBS     | 0.65 µm      | Male Luer Slip      | Pink                | Yes      | 50            | 16569-----K   |            |
| 28 mm   | SFCA     | MBS     | 0.8 µm       | Male Luer Lock      | Green               | Yes      | 50            | 16592-----K   |            |
| 28 mm   | SFCA     | MBS     | 0.8 µm       | Male Luer Lock      | Green               | Yes#     | 50            | 16592-----GUK |            |
| 28 mm   | SFCA     | MBS     | 0.8 µm       | Male Luer Lock      | Green               | No       | 500           | 16592-----Q   |            |
| 28 mm   | SFCA     | MBS     | 1.2 µm       | Male Luer Lock      | Red                 | Yes      | 50            | 17593-----K   |            |
| 28 mm   | SFCA     | MBS     | 1.2 µm       | Male Luer Lock      | Red                 | No       | 500           | 17593-----Q   |            |
| 28 mm   | SFCA     | MBS     | 5 µm         | Male Luer Lock      | Brown               | Yes      | 50            | 17594-----K   | CE-marked  |
| 28 mm   | SFCA     | MBS     | 5 µm         | Male Luer Lock      | Brown               | No       | 500           | 17594-----Q   |            |

| Dia. in mm   EFA <sup>1</sup>                                      | Membrane | Housing | Pore Size           | Connector Outlet | Color   Printing | Sterile* | Qty./ Pkg. | Order No.   | CE Marking |
|--|----------|---------|---------------------|------------------|------------------|----------|------------|-------------|------------|
| <b>Minisart® NML Plus</b> (Glass Fiber 0.7 µm <sup>2</sup> + SFCA) |          |         |                     |                  |                  |          |            |             |            |
| 28 mm  | GF+SFCA  | MBS     | 0.2 µm              | Male Luer Lock   | Blue             | Yes      | 50         | 17823-----K |            |
| 28 mm  | GF+SFCA  | MBS     | 0.2 µm              | Male Luer Lock   | Blue             | No       | 500        | 17823-----Q |            |
| 28 mm  | GF+SFCA  | MBS     | 0.45 µm             | Male Luer Lock   | Yellow           | Yes      | 50         | 17829-----K |            |
| 28 mm  | GF+SFCA  | MBS     | 0.45 µm             | Male Luer Lock   | Yellow           | No       | 500        | 17829-----Q |            |
| 28 mm  | GF+SFCA  | MBS     | 1.2 µm              | Male Luer Lock   | Red              | No       | 500        | 17825-----Q |            |
| 28 mm  | GF       | MBS     | 0.7 µm <sup>2</sup> | Male Luer Lock   | White            | No       | 50         | 17824-----K |            |
| 28 mm  | GF       | MBS     | 0.7 µm <sup>2</sup> | Male Luer Lock   | White            | No       | 500        | 17824-----Q |            |

\* Sterilized Minisart® units are individually packaged. If not stated otherwise, Minisart® are sterilized by ethylene oxide.

#-Mark indicates sterilization by gamma irradiation.

Non-sterilized Minisart® units: PES, SFCA, GF+SFCA and GF can be sterilized by ethylene oxide or gamma irradiation.

<sup>1</sup> Diameter of EFA – Effective Filtration Area

<sup>2</sup> 0.7 µm = GF particle retention ≠ pore size!

For technical product specifications, please see page 74.





Minisart® available in many configurations

## Special Applications – Medical Uses and Sterile Venting

### Choose from a Broad Range of Pore Sizes, Materials and Formats

Bacterial cell removal or particulate removal from liquids, including medical drugs, can easily be performed with Minisart® syringe filters. Minisart® has a minimal effect on the ingredients of a filtered solution.


Minisart® NML and Ophthalsart with surfactant-free cellulose acetate (SFCA) and Minisart® HY and SRP with hydrophobic PTFE have the CE mark of conformity and are frequently used for sterile filtration of aqueous and oily ear or eye solutions and other drugs. Minisart® NML with a 5 µm pore size removes particulates or coagulates in dissolved medical drugs prior to injection, offering the highest total throughput under sterile conditions – without any early

clogging. Hydrophobic PTFE filters are suitable for venting purposes and are additionally available in special formats with activated carbon.

#### Minisart® Features

- 100% optical integrity test
- Low adsorption
- Minimum extractables
- Particulate-free
- PVC-free
- Gamma-irradiated or EO-sterilized
- Many CE-marked types
- Bidirectional use possible
- Many configurations available

#### Pressure-Hold Test



Method: Pressure-hold tests were performed by connecting pre-wetted syringe filter units to a 10-fold pressure device with a pressure gauge. The pressure applied equaled 4/5 of the bubble point. Units failed the test if they released air bubbles before reaching 4/5 of the minimum bubble point. These units were remeasured to determine at which pressure air bubbles appear.

Result: Testing the pressure-hold capability of syringe filters revealed that many competitor 3's filters are dysfunctional and not intact. Filtration performed with such filters will result in unsterile filtrates or filtrates without an appropriately reduced level of particles.



Biocompatibility  
Certificate Minisart® HY



Biocompatibility  
Certificate Minisart® NML



Declaration of Conformity  
Minisart®

## Minisart® Syringe Filters – Special

## □ Ordering Information

| Dia. in<br>mm   EFA <sup>1</sup>  | Membrane         | Housing | Pore<br>Size | Connector<br>Outlet         | Color  <br>Printing | Sterile* | Qty./<br>Pkg. | Order No.     | CE Marking |
|---|------------------|---------|--------------|-----------------------------|---------------------|----------|---------------|---------------|------------|
| <b>Minisart® NML</b> (SFCA – Cellulose Acetate) Aqueous Filtration  |                  |         |              |                             |                     |          |               |               |            |
| 28 mm   | SFCA             | MBS     | 0.2 µm       | Male Luer Lock              | Blue                | Yes      | 50            | 16534-----K   | CE-marked  |
| 28 mm   | SFCA             | MBS     | 0.2 µm       | Male Luer Lock              | Blue                | Yes#     | 50            | 16534-----GUK | CE-marked  |
| 28 mm   | SFCA             | MBS     | 0.2 µm       | Male Luer Lock              | Blue                | No       | 500           | 16534-----Q   | CE-marked  |
| 28 mm   | SFCA             | MBS     | 0.2 µm       | Male Luer Slip              | Blue                | Yes      | 50            | 17597-----K   | CE-marked  |
| 28 mm   | SFCA             | MBS     | 0.2 µm       | Male Luer Slip              | Blue                | No       | 500           | 17597-----Q   | CE-marked  |
| 28 mm   | SFCA             | MBS     | 0.45 µm      | Male Luer Lock              | Yellow              | Yes      | 50            | 16555-----K   | CE-marked  |
| 28 mm   | SFCA             | MBS     | 0.45 µm      | Male Luer Lock              | Yellow              | Yes#     | 50            | 16555-----GUK | CE-marked  |
| 28 mm   | SFCA             | MBS     | 0.45 µm      | Male Luer Lock              | Yellow              | No       | 500           | 16555-----Q   | CE-marked  |
| 28 mm   | SFCA             | MBS     | 0.45 µm      | Male Luer Slip              | Yellow              | Yes      | 50            | 17598-----K   | CE-marked  |
| 28 mm   | SFCA             | MBS     | 0.45 µm      | Male Luer Slip              | Yellow              | No       | 500           | 17598-----Q   | CE-marked  |
| 28 mm   | SFCA             | MBS     | 5 µm         | Male Luer Lock              | Brown               | Yes      | 50            | 17594-----K   | CE-marked  |
| <b>Minisart® Ophthalmart</b> (SFCA – Cellulose Acetate) Aqueous Filtration                                  |                  |         |              |                             |                     |          |               |               |            |
| 28 mm   | Ophthalmart      | MBS     | 0.2 µm       | Male Luer Slip              | Pink                | Yes      | 50            | 17528-----K   | CE-marked  |
| <b>Minisart® High Flow</b> (PES – Polyethersulfone) Aqueous Filtration                                      |                  |         |              |                             |                     |          |               |               |            |
| 28 mm   | PES              | MBS     | 0.1 µm       | Male Luer Lock              | Dark Red            | Yes      | 50            | 16553-----K   |            |
| <b>Minisart® PES</b> (Polyethersulfone) Aqueous Filtration  |                  |         |              |                             |                     |          |               |               |            |
| 15 mm   | PES              | PP      | 0.22 µm      | Male Luer Slip              | White               | Yes      | 50            | 1776D-----ACK |            |
| <b>Minisart® Air</b> (Hydrophobic PTFE) Venting   |                  |         |              |                             |                     |          |               |               |            |
| 15 mm   | PTFE             | MBS     | 0.2 µm       | Male Luer Slip              | Yellow              | No       | 500           | 1751A-----Q   |            |
| 15 mm   | PTFE             | MBS     | 0.2 µm       | Male Luer Slip +<br>Needle  | Yellow              | Yes#     | 50            | 16596-----HNK |            |
| <b>Minisart® HY</b> (hydrophobic PTFE), CE-Marked, for Venting and Gas Filtration                           |                  |         |              |                             |                     |          |               |               |            |
| 26 mm   | PTFE             | MBS     | 0.2 µm       | Male Luer Lock              | Clear               | Yes      | 50            | 16596-----HYK | CE-marked  |
| 26 mm   | PTFE             | MBS     | 0.2 µm       | Male Luer Lock              | Clear               | No       | 500           | 16596-----HYQ | CE-marked  |
| 26 mm   | PTFE             | MBS     | 0.2 µm       | Male Luer Lock <sup>a</sup> | Clear               | No       | 500           | 16599-----HYQ | CE-marked  |
| 26 mm   | PTFE             | MBS     | 0.2 µm       | Hose Barbs <sup>b</sup>     | Clear               | No       | 500           | 40078-----Q   | CE-marked  |
| 26 mm   | PTFE             | MBS     | 1 µm         | Male Luer Lock              | Clear               | No       | 500           | 1659A-----HYQ |            |
| 26 mm   | PTFE             | MBS     | 1 µm         | Hose Barbs <sup>b</sup>     | Clear               | No       | 500           | 1659B-----HYQ |            |
| <b>Minisart® Acticosart with Dome Reservoir and Hydrophobic PTFE for Venting and Ultracleaning of Gases</b> |                  |         |              |                             |                     |          |               |               |            |
| 26 mm   | Active<br>Carbon | MBS     | 0.45 µm      | Male Luer Slip <sup>a</sup> | Blue                | No       | 500           | 17840-----Q   |            |

| Dia. in mm   EFA <sup>1</sup>  | Membrane | Housing | Pore Size | Connector Outlet | Color   Printing | Sterile* | Qty./ Pkg. | Order No.     | CE Marking |
|--|----------|---------|-----------|------------------|------------------|----------|------------|---------------|------------|
| <b>Minisart® SRP</b> (Hydrophobic PTFE) CE-marked Venting & Gas Filtration |          |         |           |                  |                  |          |            |               |            |
| 25 mm  | PTFE     | PP      | 0.2 µm    | Male Luer Slip   | White, Printed   | Yes      | 50         | 17575-----ACK | CE-marked  |
| 25 mm  | PTFE     | PP      | 0.2 µm    | Hose Barb        | White            | No       | 500        | 1757A-----Q   |            |

\* Sterilized Minisart® units are individually packaged. If not stated otherwise, Minisarts are sterilized by ethylene oxide.

#-mark indicates sterilization by gamma irradiation

Non-sterilized Minisart® units: SFCA can be sterilized by ethylene oxide or gamma irradiation. PTFE can be sterilized by ethylene oxide.

<sup>a</sup> Connector inlet: Male Luer slip (all other Minisart® types have female luer lock inlets)

<sup>b</sup> Hose barbs, inlet and outlet, 5 mm diameter

<sup>1</sup> Diameter of EFA – Effective Filtration Area

For technical product specifications, please see pages 72 and 74.

Do you need other quantities per package?

Are you looking for special versions or do you have requirements?

Do you need other inlet and | or outlet connectors?

Please contact us to learn more about additional Minisart® configurations that are available.





## Minisart® with Polypropylene Housing

### □ Specifications

#### Minisart® RC | SRP | NY | PES with 4 | 15 | 25 mm Diameter Membrane Filtration Area

|                        |  |
|------------------------|--|
| Housing material       | Polypropylene (PP)   |
| Membranes              | RC = Regenerated Cellulose   NY = Polyamide  <br>SRP = Hydrophobic PTFE = Polytetrafluoroethylene   PES = Polyethersulfone |
| Glass fiber prefilter  | NY Plus: Ultrapure quartz, 0.7 µm particle retention   |
| Application limits     | Max. recommended operating pressure 4.5 bar   65 psi   |
| Housing burst pressure | > 7 bar   102 psi  |
| Max. temperature       | 121 °C, 30 min (autoclavable)  |
| Sterilization          | Non-sterile Minisart® can be autoclaved or sterilized by ethylene oxide (EO)   |

| Minisart® Membrane Types   | RC 0.2 µm                      | RC 0.2 µm                      | RC 0.45 µm                     | SRP 0.2 µm                       | SRP 0.45 µm                      |
|--|--------------------------------|--------------------------------|--------------------------------|----------------------------------|----------------------------------|
| Non-sterile packs: 50 (K), 200 (S),<br>500 (Q), 1000 (R)   sterile packs:<br>individual packaged, 50 (ACK) | K   S   Q   R                  | ACK                            | K   S   Q   R                  | K   S   Q   ACK                  | K   S   Q                        |
| <b>Bubble point (≥)</b>  | With water<br>3.0 bar   44 psi | With water<br>4.6 bar   67 psi | With water<br>2.0 bar   29 psi | With ethanol<br>1.4 bar   20 psi | With ethanol<br>0.9 bar   13 psi |

#### Flow rate, 4 mm dia. = 0.07 cm<sup>2</sup> filter area | Hold-up volume<sup>1</sup>: 5 µl to 10 µl (► ml/min)

|                         |                |   |                |                |                |
|-------------------------|----------------|---|----------------|----------------|----------------|
| – For water at 1 bar    | 0.5            | – | 1.5            | – <sup>3</sup> | – <sup>3</sup> |
| – For methanol at 1 bar | 1.5            | – | 3.0            | 2.0            | 4.5            |
| – For air at 0.1 bar    | – <sup>2</sup> | – | – <sup>2</sup> | 30             | 60             |

#### Flow rate, 15 mm dia. = 1.7 cm<sup>2</sup> filter area | Hold-up volume<sup>1</sup>: 30 µl to 100 µl (► ml/min)

|                         |                |                |                |                |                |
|-------------------------|----------------|----------------|----------------|----------------|----------------|
| – For water at 1 bar    | 20             | 10             | 40             | – <sup>3</sup> | – <sup>3</sup> |
| – For methanol at 1 bar | 55             | 25             | 105            | 55             | 150            |
| – For air at 0.1 bar    | – <sup>2</sup> | – <sup>2</sup> | – <sup>2</sup> | 800            | 1600           |

#### Flow rate, 25 mm dia. = 4.8 cm<sup>2</sup> filter area | Hold-up volume<sup>1</sup>: 100 µl to 200 µl (► ml/min)

|                         |                |                |                |                |                |
|-------------------------|----------------|----------------|----------------|----------------|----------------|
| – For water at 1 bar    | 80             | 50             | 160            | – <sup>3</sup> | – <sup>3</sup> |
| – For methanol at 1 bar | 160            | 90             | 325            | 160            | 260            |
| – For air at 0.1 bar    | – <sup>2</sup> | – <sup>2</sup> | – <sup>2</sup> | 1800           | 3000           |

|  |  |     |    |                                 |                                 |
|--|--|-----|----|---------------------------------|---------------------------------|
| <b>Water penetration point</b>   | –  | –   | –  | > 4.0 bar   58 psi <sup>3</sup> | > 3.0 bar   44 psi <sup>3</sup> |
| <b>Sterile filtration capability acc. to the bacteria challenge text</b> | No <sup>5</sup>                                      | Yes | No | Yes                             | No                              |
| <b>Pyrogen-free according to the USP</b>                                 | Yes  |     |    |                                 |                                 |
| <b>Cytotoxicity (17575-ACK)</b>  | No inhibition with MRC-5 (human lung cells) and L929 |     |    |                                 |                                 |

| Minisart® Membrane Types  | NY 0.2 µm  | NY 0.45 µm                     | NY Plus 0.2 µm                 | NY Plus 0.45 µm                | PES 0.2 µm                     |
|---|--|--------------------------------|--------------------------------|--------------------------------|--------------------------------|
| Non-sterile packs: 50 (K), 200 (S), 500 (Q), 1000 (R)   sterile packs: individual packaged, 50 (ACK)                    | K   Q   R   ACK                                      | K   Q   R   ACK                | K   Q                          | K   Q                          | K   Q   ACK                    |
| <b>Bubble point (≥)</b>   | With water<br>3.0 bar   44 psi                       | With water<br>2.0 bar   29 psi | With water<br>3.0 bar   44 psi | With water<br>2.0 bar   29 psi | With water<br>3.2 bar   46 psi |
| <b>Flow rate, 4 mm dia. = 0.07 cm<sup>2</sup> filter area   Hold-up volume<sup>1</sup>: 5 µl to 10 µl (► ml/min)</b>    |  |                                |                                |                                |                                |
| – For water at 1 bar  | –  | –                              | –                              | –                              | 1.5                            |
| – For methanol at 1 bar   | –  | –                              | –                              | –                              | – <sup>4</sup>                 |
| – For air at 0.1 bar  | –  | –                              | –                              | –                              | – <sup>2</sup>                 |
| <b>Flow rate, 15 mm dia. = 1.7 cm<sup>2</sup> filter area   Hold-up volume<sup>1</sup>: 30 µl to 100 µl (► ml/min)</b>  |  |                                |                                |                                |                                |
| – For water at 1 bar  | 20   | 40                             | –                              | –                              | 40                             |
| – For methanol at 1 bar   | 40   | 110                            | –                              | –                              | – <sup>4</sup>                 |
| – For air at 0.1 bar  | – <sup>2</sup>                                       | – <sup>2</sup>                 | –                              | –                              | – <sup>2</sup>                 |
| <b>Flow rate, 25 mm dia. = 4.8 cm<sup>2</sup> filter area   Hold-up volume<sup>1</sup>: 100 µl to 200 µl (► ml/min)</b> |  |                                |                                |                                |                                |
| – For water at 1 bar  | 50   | 100                            | 50                             | 100                            | 100                            |
| – For methanol at 1 bar   | 70   | 200                            | 70                             | 200                            | – <sup>4</sup>                 |
| – For air at 0.1 bar  | – <sup>2</sup>                                       | – <sup>2</sup>                 | – <sup>2</sup>                 | – <sup>2</sup>                 | – <sup>2</sup>                 |
| <b>Water penetration point</b>  | –  | –                              | –                              | –                              | –                              |
| <b>Sterile filtration capability acc. to the bacteria challenge test</b>  | Yes  | No                             | Yes                            | No                             | Yes                            |
| <b>Pyrogen-free according to the USP</b>  |  |                                |                                |                                |                                |
| <b>Cytotoxicity (17575-ACK)</b>   | No inhibition with MRC-5 (human lung cells) and L929 |                                |                                |                                |                                |

<sup>1</sup> Hold-up volume after air purge

<sup>2</sup> Hydrophilic membranes can filter dry air or gas but become impermeable to air or gas when wetted!

<sup>3</sup> Hydrophobic membranes cannot be wetted with aqueous solutions unless you overcome their water penetration point or pre-wet them using an organic solvent (e.g. ethanol)

<sup>4</sup> PES is suitable for solutions only containing up to 30% MeOH

<sup>5</sup> According to the bacterial challenge test (BCT) with 10<sup>7</sup> Brevundimonas diminuta. Non-sterile RC Minisart® types are optimized for sample preparation and are not suitable for sterile filtration according to the bacteria challenge test. All other non-sterile Minisart® types with 0.2 mm pore size can be sterilized by autoclaving or EO before use for sterile filtration.

## Specifications

**Minisart® Air** with 15 mm accessible membrane filtration area diameter, 100 µm hold-up volume<sup>1</sup>

|  |  |     |    |     |     |    |    |     |
|--|--|-----|----|-----|-----|----|----|-----|
| <b>Water penetration point</b>   | –  | –   | –  | –   | –   | –  | –  | –   |
| <b>Sterile filtration capability<sup>4</sup> acc. to the bacteria challenge test</b> | Yes  | Yes | No | Yes | No  | No | No | No  |
| <b>Pyrogen-free according to the USP</b>   |  |     |    | Yes | Yes |    |    | Yes |
| <b>Cytotoxicity</b>  | No inhibition with MRC-5 (human lung cells) and L929 |     |    |     |     |    |    |     |

| Minisart® Membrane Type   | GF + CA<br>0.2 µm                                    | GF + CA<br>0.45 µm                   | GF + CA<br>1.2 µm                    | GF<br>0.7 µm                        | PTFE<br>0.2 µm                         | PTFE<br>1.0 µm                        | Acticosart                             | PTFE (Air)<br>0.2 µm                   |
|---|--|--------------------------------------|--------------------------------------|-------------------------------------|--|---------------------------------------|--|--|
| Non-sterile packages:<br>500 (Q, HYQ), 1000 (R),<br>sterile packs:<br>individually packaged,<br>50 (K, GUK, HYK, HNK) | K   Q  | K   Q                                | Q                                    | K   Q                               | HYK  <br>HYQ   Q                       | HYQ                                   | Q                                      | Q   HNK                                |
| <b>Bubble point (≥)</b>   | With<br>water<br>3.2 bar  <br>46 psi                 | With<br>water<br>2.0 bar  <br>29 psi | With<br>water<br>0.7 bar  <br>10 psi | With<br>water<br>0.5 bar  <br>7 psi | With<br>ethanol<br>1.4 bar  <br>20 psi | With<br>ethanol<br>0.5 bar  <br>7 psi | With<br>ethanol<br>0.9 bar  <br>13 psi | With<br>ethanol<br>0.9 bar  <br>13 psi |
| <b>Flow rate for<sup>2 3</sup> (► ml/min)</b>   |  |                                      |                                      |                                     |  |                                       |  |  |
| 28 mm dia. for water at 1 bar   | 60   | 160                                  | 350                                  | 450                                 | –                                      | –                                     | –                                      | –                                      |
| 15 mm dia. for air at 0.1 bar   | –  | –                                    | –                                    | –                                   | –                                      | –                                     | –                                      | 2000                                   |
| 26 mm dia. for air at 0.1 bar   | –  | –                                    | –                                    | –                                   | 2000                                   | 4000                                  | 2300                                   | –                                      |
| <b>Water penetration point</b>  | –  | –                                    | –                                    | –                                   | > 4.0 bar  <br>58 psi <sup>3</sup>     | > 1.5 bar  <br>22 psi <sup>3</sup>    | N.a.                                   | > 3.0 bar  <br>44 psi <sup>3</sup>     |
| <b>Sterile filtration capability<sup>4</sup><br/>according to the bacteria<br/>challenge test</b>                     | Yes  | No                                   | No                                   | No                                  | Yes                                    | No                                    | N.a.                                   | Yes                                    |
| <b>Pyrogen-free according<br/>to the USP</b>  | Yes  |                                      |                                      |                                     |  |                                       |  |  |
| <b>Cytotoxicity</b>   | No inhibition with MRC-5 (human lung cells) and L929 |                                      |                                      |                                     |  |                                       |  |  |

<sup>1</sup> Hold-up volume after air purge

<sup>2</sup> Hydrophilic membranes can filter dry air or gas but become impermeable to air or gas when wetted!

<sup>3</sup> Hydrophobic membranes cannot be wetted with aqueous solutions unless you overcome their water penetration point.

<sup>4</sup> According to bacterial challenge test (BCT) with 10<sup>7</sup> Brevundimonas diminuta. All non-sterile Minisart® types listed above can be sterilized according to the method recommended in this table.

\* Minisart® Air can be sterilized by gamma rays according to the following parameters: Range 25 – 40 kGy (validated with 50 kGy).

## ■ Sartolab® P20 and Sartolab® P20 Plus

Sterile Filtration of Sample Volumes of Up to 5 Liters



### Using Sartolab® P20 or Sartolab® P20 Plus in Available Systems

#### Systems with Luer Lock Connectors

Sartolab® units with a luer lock inlet fit directly onto the corresponding connectors of tubing from peristaltic pumps. This configuration can also be used with luer lock syringes.



#### Systems with Tubing

If the liquid to be filtered will be pumped through open tubing from a pressure tank or a peristaltic pump, then a Sartolab® unit with a hose nipple inlet is required. The stepwise increase in diameter of the nipple, from 6 mm to 12 mm, makes it suitable for a wide range of tubing.

#### Application

The Sartolab® P20 is a ready-to-use pressure filtration unit for sterile filtration of media and aqueous solutions in batches ranging from 100 ml to 5 l. For media that contain sera and difficult-to-filter solutions, a Sartolab® P20 Plus unit with an incorporated prefilter is also available.


#### Security

The effectiveness of every batch of 0.2 µm cellulose acetate membranes for sterile filtration is confirmed by bacteria challenge tests (HIMA) using *Brevundimonas diminuta*. Only biosafe

material is used in the filtration units; they have been proven non-toxic by passing the USP plastics test for toxicity. Tests with MRC-5 human lung cells on cellulose acetate membranes and glass fiber prefilters showed no cytotoxic effects. Finished units are also tested for their sterile filtration capability and for housing and membrane integrity.

#### Fast Filtration


The combination of a large filtration area (20 cm<sup>2</sup>) and the optimal design of the filter support guarantees high flow rates with high total throughputs. Automatic venting of any trapped air through the PTFE membrane-protected vent ports ensures that the entire filter surface is used for effective filtration.



Choice of inlet connector,  
hose nipple or luer lock

With or without  
binder-free glass fiber  
prefilter

Polycarbonate housing  
and surfactant-free  
cellulose acetate  
membrane filter



20 cm<sup>2</sup> filter area,  
optimal flow path and  
automatic venting

Outlet side with protective  
filling bell

## Specifications

### Technical Specifications for Sartolab® P20 and Sartolab® P20 Plus Units with SFCA Membrane

| Properties                         | Description   |   |
|------------------------------------|---|---|
|                                    | Sartolab® P20 18052 18053   | Sartolab® P20 Plus 18056 18058  |
| Filter material                    | SFCA, type 12587, 0.2 µm pore size and PTFE                                       | SFCA, type 12587 0.2 µm pore size plus GF, 100% free of binding agents and PTFE |
| Housing material                   | Polycarbonate   | Polycarbonate   |
| Color code                         | Transparent   | Transparent   |
| Filter diameter                    | 64 mm   | 64 mm   |
| Connector inlet                    | Female luer lock or stepped hose nipple with 6 mm to 12 mm outer diameter         | Female luer lock or stepped hose nipple with 6 mm to 12 mm outer diameter       |
| Connector outlet                   | Hose nipple   | Hose nipple   |
| Filling bell                       | Yes   | Yes   |
| Filtration area                    | 20 cm <sup>2</sup>  | 20 cm <sup>2</sup>  |
| Hold-up volume before bubble point | 1 ml  | Approx 1.5 ml   |
| Housing burst pressure             | > 5 bar   72.5 psi  | > 5 bar   72.5 psi  |
| Bubble point                       | ≥ 3.2 bar   46.4 psi  | ≥ 3.2 bar   46.4 psi  |
| Max. recommended inlet pressure    | 3 bar   43.5 psi  | 3 bar   43.5 psi  |
| Flow rate for water                | ≥ 250 ml/min at Δp = 1 bar   14.5 psi   | ≥ 250 ml/min at Δp = 1 bar   14.5 psi   |
| Filtration range                   | 100 ml to max. 5 l  | 100 ml to max. 10 l   |
| pH-range                           | 4–8   | 4–8   |
| Non-specific protein adsorption    | No loss of protein detectable (filtration of γ globulin, method acc. to Bradford) | < 80 µg/cm <sup>2</sup> (filtration of γ globulin, method acc. to Bradford)     |
| Sterilization                      | EO sterilization  | EO sterilization  |
| Biosafety                          | Class VI Plastics Test  | Class VI Plastics Test  |
| Operating instructions             | Directions for use included in each box   | Directions for use included in each box   |

### Technical Specifications for Sartolab® P20 and Sartolab® P20 Plus Units with PES Membrane

|                                    | Sartolab® P20 Plus 18068  | Sartolab® P20 18075   |
|------------------------------------|---|---|
| Filter material                    | PES, type 15407 MI 0.2 µm pore size plus GF, 100% free of binding agents and PTFE | PES, type 15407 MI 0.2 µm pore size                                       |
| Housing material                   | Polycarbonate   | Polycarbonate   |
| Color code                         | Transparent   | Transparent   |
| Filter diameter                    | 64 mm   | 64 mm   |
| Connector inlet                    | Female luer lock or stepped hose nipple with 6 mm to 12 mm outer diameter         | Female luer lock or stepped hose nipple with 6 mm to 12 mm outer diameter |
| Connector outlet                   | Hose nipple   | Hose nipple   |
| Filling bell                       | Yes   | No  |
| Filtration area                    | 20 cm <sup>2</sup>  | 20 cm <sup>2</sup>  |
| Hold-up volume before bubble point | Approx 1.5 ml   | 1 ml  |
| Housing burst pressure             | > 5 bar   72.5 psi  | > 5 bar   72.5 psi  |

| Properties                      | Description  |  |
|---------------------------------|--|--|
|                                 | Sartolab® P20 Plus 18068   | Sartolab® P20 18075  |
| Bubble point                    | ≥ 3.2 bar   46.4 psi   | ≥ 3.2 bar   46.4 psi   |
| Max. recommended inlet pressure | 3 bar   43.5 psi   | 3 bar   43.5 psi   |
| Flow rate for water             | 400 ml/min at $\Delta p = 1$ bar   14.5 psi  | 400 ml/min at $\Delta p = 1$ bar   14.5 psi  |
| Filtration range                | 100 ml–max. 10 l   | 100 ml–max. 5 l  |
| pH range                        | 1 to 8   | 1 to 8   |
| Non-specific protein adsorption | < 80 µg/cm <sup>2</sup> (filtration of $\gamma$ globulin, method acc. to Bradford) | No loss of protein detectable (filtration of $\gamma$ globulin, method acc. to Bradford) |
| Sterilization                   | EO sterilization   | EO sterilization   |
| Biosafety                       | Class VI Plastics Test   | Class VI Plastics Test   |
| Operating instructions          | Directions for use included in each box  | Directions for use included in each box  |

## □ Ordering Information

| Type                      | Membrane  | Housing | Pore Size | Inlet       | Outlet      | Sterile | Qty./Pkg. | Order No.   |
|---------------------------|-----------|---------|-----------|-------------|-------------|---------|-----------|-------------|
| <b>Sartolab® P20</b>      |           |         |           |             |             |         |           |             |
|                           | SFCA      | PC      | 0.2 µm    | Hose Nipple | Hose Nipple | Yes     | 10        | 18052-----D |
|                           | SFCA      | PC      | 0.2 µm    | Luer Lock   | Hose Nipple | Yes     | 10        | 18053-----D |
|                           | PES       | PC      | 0.2 µm    | Luer Lock   | Hose Nipple | Yes     | 10        | 18075-----D |
| <b>Sartolab® P20 Plus</b> |           |         |           |             |             |         |           |             |
|                           | SFCA + GF | PC      | 0.2 µm    | Hose Nipple | Hose Nipple | Yes     | 10        | 18056-----D |
|                           | SFCA + GF | PC      | 0.2 µm    | Luer Lock   | Hose Nipple | Yes     | 10        | 18058-----D |
|                           | PES + GF  | PC      | 0.2 µm    | Luer Lock   | Hose Nipple | Yes     | 10        | 18068-----D |

SFCA – Cellulose Acetate, PES – Polyethersulfone, GF – Glass Fiber Prefilter, PC – Polycarbonate



## ■ Sartolab® 150v

Disposable PES Vacuum Filtration Unit for Volumes of Up to 15 Liters



### Description

Sartolab® 150v is a disposable, sterile, ready-to-use membrane filter capsule for highest convenience. Sartolab® 150v capsules are made of a unique hydrophilic polyethersulfone (PES) membrane providing outstanding total throughput, flow rate performance, low extractables and broadest chemical compatibility.

### Applications

Typical applications include sterilizing-grade filtration of:

- Biological fluids
- Purified water
- Media
- Buffers

### Compatibility

The PES membrane is compatible with a pH range from 1 to 14, making Sartolab® 150v ideal for filtration of solutions with either high or low pH values.

### Microbiological Retention

Sartolab® 150v capsules rated to 0.2 µm are fully validated as sterilizing-grade filters according to HIMA and ASTM F-838-05 guidelines.

### Quality Control

Each individual element is integrity tested by the bubble point and diffusion test methods before final assembly. Sartolab® 150v capsules are designed, developed and manufactured in accordance with a ISO 9001:2000 certified quality management system.

### Performance

The unique pleated filter construction, combined with the highly asymmetric pore structure of the PES membrane, provides excellent flow rates and superior total throughput performance, especially in comparison to conventional disc filter systems.


## □ Specifications

### Materials and Data


|                  |                        |
|------------------|------------------------|
| Filter           | Polyethersulfone (PES) |
| Support fleece   | Polypropylene fleece   |
| Housing          | Polypropylene          |
| Vent             | PTFE                   |
| Pore size        | 0.2 µm                 |
| Filtration area  | 150 cm <sup>2</sup>    |
| Flow rate        | 1.5 L/min              |
| Total throughput | 15 L                   |

## □ Ordering Information


| Qty./Pkg. | Order No. |
|-----------|-----------|
| 3         | 18080-M   |

**Easy to Use**


Place Sartolab® 150v on the bottle neck. Vacuum is automatically created when used with a water jet pump or a different vacuum source.



Create vacuum to start filtration.



To fill multiple receiving bottles, break the vacuum by opening the valve on top of the Sartolab® 150v housing.



To continue filtration with the next receiving bottle, close the valve.

## ■ Sartolab® RF | BT

### Disposable Filter Systems and Bottle-Top Filters



Sartolab® single-use sterile filter systems and bottle-top filters are designed for the vacuum filtration of tissue culture media and components, biological fluids, and other aqueous solutions.

The Sartolab® single-use 150 mL, 250 mL, 500 mL and 1,000 mL bottles are designed as storage containers for sterile media, buffers, or other aqueous solutions.



These products are for laboratory use only and not for human parenteral applications.

#### Materials

The filter funnels, dust covers and receiver bottles are manufactured from virgin, heavy metal-free polystyrene. The tubing adapters, filter adapters, and the plug seal caps are made of heavy metal-free polyethylene. Sartolab® filter systems are available with polyethersulfone membranes. All units are sterilized by gamma irradiation.

#### Performance

The filter units contain membranes integrally sealed to a support grid designed to maximize flow and reduce foaming and protein denaturation.

The membrane is compatible with most aqueous solutions and tested for use in cell culture applications.

#### Filter Systems

The filter adapter utilizes a gasket design to ensure a vacuum-tight seal on the receiver | storage bottle. Each filter unit also contains a convenient tubing adapter that will fit most vacuum hoses.

The bottles are single-use containers.

They cannot withstand autoclaving or use at temperatures greater than 70°C. The suitability of the bottles for storage of solutions below 0°C depends both on the solution and the storage conditions. Many aqueous solutions, including culture media, have been successfully frozen and stored at temperatures down to -20°C. However, a trial run under actual conditions is strongly recommended to test the suitability of the bottles for frozen storage.

### □ Specifications

| Pore Size | Membrane Material | Characteristics   |
|-----------|-------------------|---|
| 0.22 µm   | Polyethersulfone  | Very low protein binding and low extractables, fast flow rate |
| 0.1 µm    | Polyethersulfone  | Very low protein binding and low extractables, fast flow rate |

**Bottle-Top Filters**

The filter adapter is available with a 45 mm thread finish, and is designed to ensure a vacuum-tight seal on customer-supplied bottles with the appropriate thread finish. Each filter unit also contains a tubing adapter that will fit most vacuum hoses.

**Chemical Compatibility**

The mechanical strength, color, appearance, and dimensional stability of filter systems, bottle-top filters and plastic bottles are affected to varying degrees by the chemicals with which they come in contact. Specific operating conditions, especially temperature, will also affect their chemical resistance. A table is provided to serve as a general guideline for the chemical resistance of Sartolab® single-use sterile filters and bottles.

**Chemical Resistance of Sartolab® Filters**

| Chemical Class   | Membrane (PES) | Housing (PS) |
|------------------|----------------|--------------|
| Weak Acids       | 3              | 1            |
| Strong Acids     | 3              | 2            |
| Alcohols         | 1              | 2            |
| Aldehydes        | 3              | 3            |
| Aliphatic Amines | 1              | 3            |
| Aromatic Amines  | 3              | 3            |
| Bases            | 3              | 1            |
| Esters           | 3              | 3            |
| Hydrocarbons     | 3              | 3            |
| Ketones          | 3              | 3            |

Key: 1. recommended


2. may be suitable for some applications; a trial run is recommended

3. not recommended. PS, polystyrene; PES, polyethersulfone.

## Ordering Information

| Volume  | Membrane    | Filter Area        | Qty./Pkg. | Order No.   |
|---|-------------|--------------------|-----------|-------------|
| <b>Sartolab® RF Filtration System Including Collection Bottle</b>           |             |                    |           |             |
| 150 ml  | 0.22 µm PES | 18 cm <sup>2</sup> | 12        | 180C1-----E |
| 250 ml  | 0.22 µm PES | 24 cm <sup>2</sup> | 12        | 180C7-----E |
| 500 ml  | 0.22 µm PES | 39 cm <sup>2</sup> | 12        | 180C2-----E |
| 1,000 ml  | 0.22 µm PES | 62 cm <sup>2</sup> | 12        | 180C3-----E |
| 1,000 ml  | 0.1 µm PES  | 62 cm <sup>2</sup> | 12        | 180C8-----E |
| <b>Sartolab® BT Bottle-Top Filters Without Integrated Collection Bottle</b> |             |                    |           |             |
| 150 ml  | 0.22 µm PES | 18 cm <sup>2</sup> | 48        | 180C4-----K |
| 500 ml  | 0.22 µm PES | 39 cm <sup>2</sup> | 12        | 180C5-----E |
| 1,000 ml  | 0.22 µm PES | 62 cm <sup>2</sup> | 12        | 180C6-----E |

## Water Throughput



## Chemical Compatibility

|                           | Material     |               |                |             |                |                 |             |            | Minisart® Types    |               |             |                    |                  |              |               |              |              |                   |               |
|---------------------------|--------------|---------------|----------------|-------------|----------------|-----------------|-------------|------------|--------------------|---------------|-------------|--------------------|------------------|--------------|---------------|--------------|--------------|-------------------|---------------|
|                           | PES membrane | SFCA membrane | PTFE membrane  | RC membrane | Nylon membrane | GF depth filter | Housing MBS | Housing PP | Minisart® HighFlow | Minisart® NML | Ophthalsart | Minisart® NML Plus | Minisart® NML GF | Minisart® HY | Minisart® Air | Minisart® RC | Minisart® NY | Minisart® NY Plus | Minisart® SRP |
| Filter Membrane           | PES          | SFCA          | PTFE           | RC          | PA             |                 |             |            | PES                | SFCA          | SFCA        |                    |                  | PTFE         | RC            | PA           | PA           | PTFE              | PES           |
| Pre-Filter                | GF           |               |                |             |                |                 |             |            | -                  | -             | GF          | GF                 | -                | -            | -             | GF           | -            | -                 |               |
| Housing Material          | MBS PP       |               |                |             |                |                 |             |            | MBS                | MBS           | MBS         | MBS                | MBS              | PP           | PP            | PP           | PP           | PP                |               |
| Sterilization             |              |               |                |             |                |                 |             |            |                    |               |             |                    |                  |              |               |              |              |                   |               |
| Ethylene oxide            | ++           | ++            | ++             | ++          | ++             | ++              | ++          | ++         | ++                 | ++            | ++          | ++                 | ++               | ++           | ++            | ++           | ++           | ++                | ++            |
| Gamma irradiation         | ++           | ++            | - <sup>1</sup> | ++          | -              | ++              | ++          | -          | ++                 | ++            | ++          | ++                 | - <sup>1</sup>   | -            | -             | -            | -            | -                 |               |
| Autoclaving 121°C, 30 min | ++           | ++            | ++             | ++          | ++             | ++              | -           | ++         | -                  | -             | -           | -                  | -                | ++           | ++            | ++           | ++           | ++                |               |
| Solvents                  |              |               |                |             |                |                 |             |            |                    |               |             |                    |                  |              |               |              |              |                   |               |
| Acetone                   | -            | -             | ++             | ++          | ++             | ++              | -           | ++         | -                  | -             | -           | -                  | -                | ++           | ++            | ++           | ++           | -                 |               |
| Acetonitrile              | -            | -             | ++             | ++          | ++             | ++              | -           | ++         | -                  | -             | -           | -                  | -                | ++           | ++            | ++           | ++           | -                 |               |
| Gasoline                  | +            | ++            | ++             | ++          | ++             | ++              | +           | ++         | +                  | +             | +           | +                  | +                | ++           | ++            | ++           | ++           | +                 |               |
| Benzene                   | +            | +             | ++             | ++          | ++             | ++              | -           | ++         | -                  | -             | -           | -                  | -                | ++           | ++            | ++           | ++           | +                 |               |
| Benzyl alcohol            | +            | +             | ++             | ++          | ++             | ++              | -           | +          | -                  | -             | -           | -                  | -                | ++           | ++            | ++           | ++           | +                 |               |
| n-Butyl acetate           | -            | -             | ++             | ++          | ++             | ++              | -           | ++         | -                  | -             | -           | -                  | -                | ++           | ++            | ++           | ++           | -                 |               |
| n-Butanol                 | ++           | ++            | ++             | ++          | ++             | ++              | +           | ++         | +                  | +             | +           | +                  | +                | ++           | ++            | ++           | ++           | ++                |               |
| Cellosolve                | +            | -             | ++             | ++          | ++             | ++              | -           | ++         | -                  | -             | -           | -                  | -                | ++           | ++            | ++           | ++           | +                 |               |
| Chloroform                | -            | -             | ++             | ++          | ++             | ++              | -           | ++         | -                  | -             | -           | -                  | -                | ++           | ++            | ++           | ++           | -                 |               |
| Cyclohexane               | -            | -             | ++             | ++          | ++             | ++              | +           | +          | -                  | -             | -           | +                  | +                | +            | +             | +            | +            | -                 |               |
| Cyclohexanone             | -            | -             | ++             | ++          | ++             | ++              | -           | +          | -                  | -             | -           | -                  | -                | +            | +             | +            | +            | -                 |               |
| Diethylacetamide          | -            | -             | ++             | ++          | ++             | ++              | -           | ++         | -                  | -             | -           | -                  | -                | ++           | ++            | ++           | ++           | -                 |               |
| Diethyl ether             | -            | +             | ++             | ++          | ++             | ++              | -           | ++         | -                  | -             | -           | -                  | -                | ++           | ++            | ++           | ++           | -                 |               |
| Dimethyl formamide        | -            | -             | ++             | +           | +              | ++              | -           | ++         | -                  | -             | -           | -                  | -                | +            | +             | +            | ++           | -                 |               |
| Dimethylsulfoxide         | -            | -             | ++             | ++          | ++             | ++              | -           | ++         | -                  | -             | -           | -                  | -                | ++           | ++            | ++           | ++           | -                 |               |
| Dioxane                   | -            | -             | ++             | ++          | ++             | ++              | -           | ++         | -                  | -             | -           | -                  | -                | ++           | ++            | ++           | ++           | -                 |               |
| Ethanol, 98%              | ++           | ++            | ++             | ++          | ++             | ++              | -           | ++         | -                  | -             | -           | -                  | -                | ++           | ++            | ++           | ++           | ++                |               |
| Ethyl acetate             | -            | -             | ++             | ++          | ++             | ++              | -           | +          | -                  | -             | -           | -                  | -                | +            | +             | +            | +            | -                 |               |
| Ethylene glycol           | ++           | +             | ++             | ++          | ++             | ++              | ++          | ++         | ++                 | +             | +           | ++                 | ++               | ++           | ++            | ++           | ++           | ++                |               |
| Formamide                 | ++           | -             | ++             | +           | ++             | ++              | ++          | ++         | ++                 | -             | -           | ++                 | ++               | +            | ++            | ++           | ++           | ++                |               |
| Glycerin                  | ++           | ++            | ++             | ++          | ++             | ++              | ++          | ++         | ++                 | ++            | ++          | ++                 | ++               | ++           | ++            | ++           | ++           | ++                |               |
| n-Heptane                 | +            | +             | ++             | ++          | ++             | ++              | ++          | +          | +                  | +             | +           | +                  | ++               | +            | +             | +            | +            | +                 |               |
| n-Hexane                  | +            | +             | ++             | ++          | ++             | ++              | ++          | +          | +                  | +             | +           | +                  | ++               | +            | +             | +            | +            | +                 |               |
| Isobutanol                | ++           | +             | ++             | ++          | ++             | ++              | -           | ++         | -                  | -             | -           | -                  | -                | ++           | ++            | ++           | ++           | ++                |               |
| Isopropanol               | ++           | ++            | ++             | ++          | ++             | ++              | -           | ++         | -                  | -             | -           | -                  | -                | ++           | ++            | ++           | ++           | ++                |               |
| Isopropyl acetate         | -            | -             | ++             | ++          | ++             | ++              | -           | ++         | -                  | -             | -           | -                  | -                | ++           | ++            | ++           | ++           | -                 |               |
| Methanol, 98%             | +            | -             | ++             | ++          | ++             | ++              | -           | ++         | -                  | -             | -           | -                  | -                | ++           | ++            | ++           | ++           | +                 |               |
| Methyl acetate            | -            | -             | ++             | ++          | ++             | ++              | -           | +          | -                  | -             | -           | -                  | -                | +            | +             | +            | +            | -                 |               |
| Methylene chloride        | -            | -             | ++             | ++          | ++             | ++              | -           | ++         | -                  | -             | -           | -                  | -                | ++           | ++            | ++           | ++           | -                 |               |
| Methyl ethyl ketone       | -            | +             | ++             | ++          | ++             | ++              | -           | +          | -                  | -             | -           | -                  | -                | +            | +             | +            | +            | -                 |               |
| Methyl isobutyl ketone    | -            | -             | ++             | ++          | ++             | ++              | -           | +          | -                  | -             | -           | -                  | -                | +            | +             | +            | +            | -                 |               |
| Monochlorobenzene         | +            | +             | ++             | ++          | ++             | ++              | -           | +          | -                  | -             | -           | -                  | -                | +            | +             | +            | +            | +                 |               |
| Nitrobenzene              | -            | -             | ++             | ++          | +              | ++              | -           | +          | -                  | -             | -           | -                  | -                | +            | +             | +            | +            | -                 |               |
| n-Pentane                 | ++           | ++            | ++             | ++          | ++             | ++              | +           | +          | +                  | +             | +           | +                  | +                | +            | +             | +            | +            | +                 |               |
| Perchloroethylene         | -            | -             | ++             | ++          | ++             | ++              | -           | +          | -                  | -             | -           | -                  | -                | +            | +             | +            | +            | -                 |               |

|                            | Material     |               |               |             |                |                 |             |            | Minisart® Types    |                           |                    |                  |              |               |              |              |                   |               |               |
|----------------------------|--------------|---------------|---------------|-------------|----------------|-----------------|-------------|------------|--------------------|---------------------------|--------------------|------------------|--------------|---------------|--------------|--------------|-------------------|---------------|---------------|
|                            | PES membrane | SFCA membrane | PTFE membrane | RC membrane | Nylon membrane | GF depth filter | Housing MBS | Housing PP | Minisart® HighFlow | Minisart® NML Ophthalsart | Minisart® NML Plus | Minisart® NML GF | Minisart® HY | Minisart® Air | Minisart® RC | Minisart® NY | Minisart® NY Plus | Minisart® SRP | Minisart® PES |
| Filter Membrane            | PES          | SFCA          | PTFE          | RC          | PA             |                 |             |            | PES                | SFCA                      | SFCA               |                  | PTFE         | RC            | PA           | PA           | PTFE              | PES           |               |
| Prefilter                  |              |               |               |             |                | GF              |             |            | -                  | -                         | GF                 | GF               | -            | -             | -            | GF           | -                 | -             |               |
| Housing Material           |              |               |               |             |                |                 | MBS         | PP         | MBS                | MBS                       | MBS                | MBS              | MBS          | PP            | PP           | PP           | PP                | PP            |               |
| Solvents (continued)       |              |               |               |             |                |                 |             |            |                    |                           |                    |                  |              |               |              |              |                   |               |               |
| Pyridine                   | -            | -             | ++            | ++          | ++             | ++              | -           | ++         | -                  | -                         | -                  | -                | -            | ++            | ++           | ++           | ++                | -             |               |
| Carbon tetrachloride       | -            | -             | ++            | ++          | ++             | ++              | -           | +          | -                  | -                         | -                  | -                | -            | +             | +            | +            | +                 | -             |               |
| Tetrahydrofuran            | -            | -             | ++            | ++          | ++             | ++              | -           | ++         | -                  | -                         | -                  | -                | -            | ++            | ++           | ++           | ++                | -             |               |
| Toluene                    | -            | +             | ++            | ++          | ++             | ++              | -           | +          | -                  | -                         | -                  | -                | -            | +             | +            | +            | +                 | -             |               |
| Trichloroethane            | -            | -             | ++            | ++          | +              | ++              | -           | +          | -                  | -                         | -                  | -                | -            | +             | +            | +            | +                 | -             |               |
| Trichloroethylene          | -            | +             | ++            | ++          | ++             | ++              | -           | +          | -                  | -                         | -                  | -                | -            | +             | +            | +            | +                 | -             |               |
| Xylene                     | -            | +             | ++            | ++          | ++             | ++              | -           | +          | -                  | -                         | -                  | -                | -            | +             | +            | +            | +                 | -             |               |
| Acids                      |              |               |               |             |                |                 |             |            |                    |                           |                    |                  |              |               |              |              |                   |               |               |
| Acetic acid, 25%           | +            | +             | ++            | ++          | -              | ++              | +           | ++         | +                  | +                         | +                  | +                | +            | ++            | -            | -            | ++                | +             |               |
| Acetic acid, 80%           | -            | -             | ++            | +           | -              | ++              | -           | +          | -                  | -                         | -                  | -                | -            | +             | -            | -            | +                 | -             |               |
| Hydrofluoric acid, 25%     | +            | -             | ++            | +           | -              | ++              | +           | +          | +                  | -                         | -                  | +                | +            | +             | -            | -            | +                 | +             |               |
| Hydrofluoric acid, 50%     | +            | -             | ++            | +           | -              | ++              | -           | +          | -                  | -                         | -                  | -                | -            | +             | -            | -            | +                 | +             |               |
| Perchloric acid, 25%       | -            | -             | ++            | -           | -              | ++              | -           | +          | -                  | -                         | -                  | -                | -            | -             | -            | -            | +                 | -             |               |
| Phosphoric acid, up to 10% | +            | +             | ++            | -           | -              | ++              | +           | +          | +                  | +                         | +                  | +                | +            | -             | -            | -            | +                 | +             |               |
| Phosphoric acid, 86%       | +            | +             | ++            | -           | -              | ++              | -           | +          | -                  | -                         | -                  | -                | -            | -             | -            | -            | +                 | +             |               |
| Nitric acid, 30%           | +            | -             | ++            | -           | -              | ++              | +           | +          | +                  | -                         | -                  | +                | +            | -             | -            | -            | +                 | +             |               |
| Nitric acid, conc.         | -            | -             | ++            | -           | -              | ++              | -           | -          | -                  | -                         | -                  | -                | -            | -             | -            | -            | -                 | -             |               |
| Hydrochloric acid, 15%     | ++           | +             | ++            | -           | -              | ++              | +           | +          | +                  | +                         | +                  | +                | +            | -             | -            | -            | +                 | +             |               |
| Hydrochloric acid, 20%     | ++           | -             | ++            | -           | -              | ++              | +           | +          | +                  | -                         | -                  | +                | +            | -             | -            | -            | +                 | +             |               |
| Sulfuric acid, 25%         | +            | -             | ++            | +           | -              | ++              | ++          | ++         | +                  | -                         | -                  | ++               | ++           | +             | -            | -            | ++                | +             |               |
| Sulfuric acid, 98%         | -            | -             | ++            | -           | -              | ++              | -           | -          | -                  | -                         | -                  | -                | -            | -             | -            | -            | -                 | -             |               |
| Trichloroacetic acid, 25%  | -            | -             | ++            | ++          | -              | ++              | -           | +          | -                  | -                         | -                  | -                | -            | +             | -            | -            | +                 | -             |               |
| Bases                      |              |               |               |             |                |                 |             |            |                    |                           |                    |                  |              |               |              |              |                   |               |               |
| Ammonia, 1N                | ++           | +             | ++            | +           | ++             | ++              | +           | ++         | +                  | +                         | +                  | +                | +            | +             | ++           | ++           | ++                | ++            |               |
| Ammonium hydroxide, 25%    | +            | +             | ++            | +           | ++             | +               | -           | +          | -                  | -                         | -                  | -                | -            | +             | +            | +            | +                 | +             |               |
| Potassium hydroxide, 32%   | ++           | -             | ++            | -           | +              | +               | -           | ++         | -                  | -                         | -                  | -                | -            | -             | +            | +            | ++                | ++            |               |
| Sodium hydroxide, 1N       | ++           | -             | ++            | +           | ++             | +               | -           | ++         | -                  | -                         | -                  | -                | -            | +             | ++           | +            | ++                | ++            |               |
| Sodium hydroxide, 32%      | ++           | -             | ++            | -           | +              | -               | -           | +          | -                  | -                         | -                  | -                | -            | -             | +            | -            | +                 | +             |               |
| Aqueous solutions          |              |               |               |             |                |                 |             |            |                    |                           |                    |                  |              |               |              |              |                   |               |               |
| Formaldehyde, 30%          | +            | ++            | ++            | +           | ++             | ++              | +           | +          | +                  | +                         | +                  | +                | +            | +             | +            | +            | +                 | +             |               |
| Sodium hypochlorite, 5%    | ++           | -             | ++            | -           | -              | ++              | +           | +          | +                  | -                         | -                  | +                | +            | -             | -            | -            | +                 | +             |               |
| Hydrogen peroxide, 35%     | ++           | -             | ++            | -           | -              | ++              | +           | ++         | +                  | -                         | -                  | +                | +            | -             | -            | -            | ++                | ++            |               |
| pH range                   |              |               |               |             |                |                 |             |            |                    |                           |                    |                  |              |               |              |              |                   |               |               |
| pH 1 to 14                 | -            | -             | ++            | -           | -              | ++              | -           | ++         |                    |                           |                    |                  |              |               |              |              |                   |               |               |
| pH 1 to 13                 | ++           | -             | ++            | -           | -              | ++              | -           | ++         |                    |                           |                    |                  |              |               |              |              |                   |               |               |
| pH 3 to 14                 | +            | -             | ++            | +           | ++             | ++              | -           | ++         |                    |                           |                    |                  |              |               |              |              |                   |               |               |
| pH 3 to 12                 | ++           | -             | ++            | ++          | ++             | ++              | +           | ++         |                    |                           |                    |                  |              |               |              |              |                   |               |               |
| pH 4 to 8                  | ++           | ++            | ++            | ++          | ++             | ++              | ++          | ++         |                    |                           |                    |                  |              |               |              |              |                   |               |               |

Legend

High compatibility++


Limited compatibility+

Not compatible-

<sup>1</sup> Gamma irradiation feasible for Minisart® Air

Contact time: 24 hours at 20°C. Chemical compatibilities can be influenced by various factors. Therefore, we recommend that you confirm compatibility with the liquid you want to filter by performing a trial filtration run before you start your actual filtration.





## ■ Basic Filtration

### Table of Contents

91 Filter Papers

95 Glass and Quartz Microfiber Filters

98 Membrane Filters

106 Blotting | Chromatography Papers & Blotting Membranes

108 Filtration Equipment

118 Combisart® Manifolds





Filters are indispensable for your routine work in laboratory and industrial applications. Sartorius supplies you with a broad range of filters for a myriad of filtration tasks and supports you with all your filtration challenges.

Our Product Range Covers:

- Filter papers
- Glass and quartz microfiber filters
- Extraction thimbles
- Membrane filters
- Blotting & chromatography papers & membranes
- Filtration equipment

#### **Quality Assurance and Quality Control**

Sartorius pays particular attention to continuous in-process quality control. Regular checks and exact analyses of the raw materials and each finished product assure constant high quality and product uniformity. We meet the requirements set forth by the ISO 9001 quality management system and the ISO 14001 environmental management system.

### How Do Filter Papers Work?

Filter papers are depth filters. Their efficiency is influenced by various parameters: the mechanical particulate retention, adsorption, pH, surface properties, thickness and strength of the filter paper as well as the shape, density and quantity of particles to be retained. The precipitates deposited on the filter form a "cake layer" which – depending on its density – increasingly affects the progress of an ongoing filtration and decisively affects the retention capability. Therefore, it is essential to select the perfect filter paper to ensure the best filtration results. This choice depends on the filtration method as well as on the amount and properties of the medium to be filtered, the size of the particulate solids to be removed and the required degree of clarification.

### How Do Membrane Filters Work?

Membrane filters retain particles larger than their pore sizes. Smaller particles pass through the membrane or are captured in the membrane. Such filters are used for the filtration of smaller particles and for critical applications such as sterility testing. The choice of the right membrane type depends on the specifications of the solution to be filtered. The most important parameters for this are adsorption, chemical compatibility and the particle size to be retained.










## Ash-free Filter Papers

### For Quantitative and Gravimetric Analyses

These filter papers are used for quantitative and gravimetric analyses as well as for pressure or vacuum filtration. They are made out of 100% cotton linters with an  $\alpha$ -cellulose content of > 98% and are acid-washed to make the papers ashless and achieve high purity.

#### Specifications

| Grade   | Weight (g/m <sup>2</sup> ) | Thickness (mm) | Particle Retention (μm) | Filtration (s) | Precipitates            | Properties                                      |
|---|----------------------------|----------------|-------------------------|----------------|-------------------------|---|
|  388   | 84                         | 0.21           | 12 to 15                | 10             | Coarse crystalline      | Wide-pore, loose structure, fast filtering      |
|  389   | 84                         | 0.19           | 8 to 12                 | 20             | Medium-fine crystalline | Medium- to wide-pore, medium fast filtering     |
|  389 F | 84                         | 0.19           | 8 to 12                 | 20             | Medium-fine crystalline | Medium- to wide-pore, medium fast filtering     |
|  392   | 84                         | 0.17           | 5 to 8                  | 50             | Fine crystalline        | Medium dense, medium fast filtering             |
|  390   | 84                         | 0.16           | 3 to 5                  | 100            | Fine crystalline        | Narrow-pore, dense, slow filtering              |
|  391   | 84                         | 0.15           | 2 to 3                  | 180            | Very fine crystalline   | Fine-pore, dense, very slow filtering           |
|  393   | 100                        | 0.17           | 1 to 2                  | 300            | Very fine crystalline   | Very fine-pore, very dense, very slow filtering |

#### Ordering Information



##### Filter Discs, 100 Pieces

| Dia. in mm | Grade 388    | Grade 389    | Grade 389 F  | Grade 390    | Grade 391    | Grade 392    | Grade 393    |
|------------|--------------|--------------|--------------|--------------|--------------|--------------|--------------|
| 55         | FT-3-101-055 | FT-3-102-055 | FT-3-112-055 | FT-3-103-055 | FT-3-104-055 | FT-3-105-055 | FT-3-127-055 |
| 70         | FT-3-101-070 | FT-3-102-070 |              | FT-3-103-070 | FT-3-104-070 | FT-3-105-070 | FT-3-127-070 |
| 90         | FT-3-101-090 | FT-3-102-090 | FT-3-112-090 | FT-3-103-090 | FT-3-104-090 | FT-3-105-090 | FT-3-127-090 |
| 110        | FT-3-101-110 | FT-3-102-110 | FT-3-112-110 | FT-3-103-110 | FT-3-104-110 | FT-3-105-110 | FT-3-127-110 |
| 125        | FT-3-101-125 | FT-3-102-125 | FT-3-112-125 | FT-3-103-125 | FT-3-104-125 | FT-3-105-125 | FT-3-127-125 |
| 150        | FT-3-101-150 | FT-3-102-150 | FT-3-112-150 | FT-3-103-150 | FT-3-104-150 | FT-3-105-150 | FT-3-127-150 |
| 185        | FT-3-101-185 | FT-3-102-185 | FT-3-112-185 | FT-3-103-185 | FT-3-104-185 | FT-3-105-185 | FT-3-127-185 |
| 240        | FT-3-101-240 | FT-3-102-240 |              | FT-3-103-240 | FT-3-104-240 | FT-3-105-240 | FT-3-127-240 |



##### Folded Filters, 100 Pieces

| Dia. in mm | Grade 388    | Grade 389    | Grade 389 F  | Grade 390    | Grade 391    | Grade 392    |
|------------|--------------|--------------|--------------|--------------|--------------|--------------|
| 110        | FT-4-101-110 | FT-4-102-110 |              | FT-4-103-110 | FT-4-104-110 | FT-4-105-110 |
| 125        | FT-4-101-125 | FT-4-102-125 |              | FT-4-103-125 | FT-4-104-125 | FT-4-105-125 |
| 150        | FT-4-101-150 | FT-4-102-150 |              | FT-4-103-150 | FT-4-104-150 | FT-4-105-150 |
| 185        | FT-4-101-185 | FT-4-102-185 | FT-4-112-185 | FT-4-103-185 | FT-4-104-185 | FT-4-105-185 |
| 240        | FT-4-101-240 | FT-4-102-240 |              |              | FT-4-104-240 |              |



##### Sheets in 580 × 580 mm, 100 Pieces

| Grade 388       | Grade 389       | Grade 390       | Grade 391       | Grade 392       | Grade 393       |
|-----------------|-----------------|-----------------|-----------------|-----------------|-----------------|
| FT-2-101-580580 | FT-2-102-580580 | FT-2-103-580580 | FT-2-104-580580 | FT-2-105-580580 | FT-2-127-580580 |

Other dimensions are available on request

## Wet-strengthened Filter Papers

### For Qualitative Analyses

These qualitative filter papers are essentially used for analytical purposes and routine analyses, whenever no gravimetric analyses are required. They are wet-strengthened and can be used for pressure and vacuum filtration. They are made of refined pulp and linters with an > 95%  $\alpha$ -cellulose content and are very pure with an ash content  $\leq 0.1\%$ .

### Specifications

| Grade | Weight (g/m <sup>2</sup> ) | Thickness (mm) | Particle Retention (μm) | Filtration (s) | Precipitates            | Properties                                      |
|-------|----------------------------|----------------|-------------------------|----------------|-------------------------|---|
| 1288  | 84                         | 0.21           | 12 to 15                | 10             | Coarse crystalline      | Wide-pore, loose structure, fast filtering      |
| 1289  | 84                         | 0.21           | 8 to 12                 | 20             | Medium-fine crystalline | Medium- to wide-pore, medium fast filtering     |
| 1292  | 84                         | 0.17           | 5 to 8                  | 50             | Fine crystalline        | Medium dense, medium fast filtering             |
| 1290  | 84                         | 0.15           | 3 to 5                  | 100            | Fine crystalline        | Narrow-pore, dense, slow filtering              |
| 1291  | 84                         | 0.15           | 2 to 3                  | 180            | Very fine crystalline   | Fine-pore, dense, very slow filtering           |
| 293   | 80                         | 0.15           | 1 to 2                  | 300            | Very fine crystalline   | Very fine-pore, very dense, very slow filtering |

### Ordering Information



#### Filter Discs, 100 Pieces

| Dia. in mm | Grade 1288   | Grade 1289   | Grade 1290   | Grade 1291   | Grade 1292   | Grade 293    |
|------------|--------------|--------------|--------------|--------------|--------------|--------------|
| 55         | FT-3-206-055 | FT-3-207-055 | FT-3-208-055 | FT-3-209-055 | FT-3-210-055 | FT-3-211-055 |
| 70         | FT-3-206-070 | FT-3-207-070 | FT-3-208-070 | FT-3-209-070 | FT-3-210-070 | FT-3-211-070 |
| 90         | FT-3-206-090 | FT-3-207-090 | FT-3-208-090 | FT-3-209-090 | FT-3-210-090 | FT-3-211-090 |
| 110        | FT-3-206-110 | FT-3-207-110 | FT-3-208-110 | FT-3-209-110 | FT-3-210-110 | FT-3-211-110 |
| 125        | FT-3-206-125 | FT-3-207-125 | FT-3-208-125 | FT-3-209-125 | FT-3-210-125 | FT-3-211-125 |
| 150        | FT-3-206-150 | FT-3-207-150 | FT-3-208-150 | FT-3-209-150 | FT-3-210-150 | FT-3-211-150 |
| 185        | FT-3-206-185 | FT-3-207-185 | FT-3-208-185 | FT-3-209-185 | FT-3-210-185 | FT-3-211-185 |
| 240        | FT-3-206-240 | FT-3-207-240 | FT-3-208-240 | FT-3-209-240 | FT-3-210-240 |              |



#### Folded Filters, 100 Pieces

| Dia. in mm | Grade 1288   | Grade 1289   | Grade 1290   | Grade 1291   | Grade 1292   | Grade 293    |
|------------|--------------|--------------|--------------|--------------|--------------|--------------|
| 110        | FT-4-206-110 | FT-4-207-110 | FT-4-208-110 | FT-4-209-110 | FT-4-210-110 |              |
| 125        | FT-4-206-125 | FT-4-207-125 | FT-4-208-125 | FT-4-209-125 | FT-4-210-125 | FT-4-211-125 |
| 150        | FT-4-206-150 | FT-4-207-150 | FT-4-208-150 | FT-4-209-150 | FT-4-210-150 | FT-4-211-150 |
| 185        | FT-4-206-185 | FT-4-207-185 | FT-4-208-185 | FT-4-209-185 | FT-4-210-185 | FT-4-211-185 |
| 240        | FT-4-206-240 | FT-4-207-240 | FT-4-208-240 | FT-4-209-240 | FT-4-210-240 | FT-4-211-240 |



#### Sheets in 580 × 580 mm, 100 Pieces

| Grade 1288      | Grade 1289      | Grade 1290      | Grade 1291      | Grade 1292      | Grade 293       |
|-----------------|-----------------|-----------------|-----------------|-----------------|-----------------|
| FT-2-206-580580 | FT-2-207-580580 | FT-2-208-580580 | FT-2-209-580580 | FT-2-210-580580 | FT-2-211-580580 |

Other dimensions are available on request

## High-Purity Filter Papers

### For Qualitative Analyses

These paper grades are used for analytical purposes that require a low ash content. Grades 292 and 292a are especially suitable for soil analyses because they are low in nitrogen. For phosphate or sodium determination, we recommend grades 131 and 132.

### Specifications

| Grade | Weight (g/m <sup>2</sup> ) | Thickness (mm) | Particle Retention (µm) | Filtration (s) | Material   |
|-------|----------------------------|----------------|-------------------------|----------------|--|
| 292   | 87                         | 0.18           | 5 to 8                  | 45             | Cotton linters, low-nitrogen and nitrates, ash content ≤ 0.06% according to DIN 54370                        |
| 292a  | 97                         | 0.19           | 4 to 7                  | 60             | Cotton linters, low-nitrogen and nitrates, ash content ≤ 0.06% according to DIN 54370                        |
| 132   | 80                         | 0.17           | 5 to 7                  | 55             | Cotton linters and refined pulp, low-phosphate and low-potassium, ash content < 0.02% according to DIN 54370 |
| 131   | 80                         | 0.16           | 3 to 5                  | 100            | Cotton linters and refined pulp, low-phosphate and low-potassium, ash content < 0.02% according to DIN 54370 |

### Ordering Information



#### Filter Discs, 100 Pieces

| Dia. in mm | Grade 131    | Grade 132    | Grade 292    | Grade 292a   |
|------------|--------------|--------------|--------------|--------------|
| 55         |              | FT-3-329-055 | FT-3-205-055 | FT-3-215-055 |
| 70         |              | FT-3-329-070 | FT-3-205-070 | FT-3-215-070 |
| 90         |              | FT-3-329-090 | FT-3-205-090 | FT-3-215-090 |
| 110        |              | FT-3-329-110 | FT-3-205-110 | FT-3-215-110 |
| 125        | FT-3-351-125 | FT-3-329-125 | FT-3-205-125 | FT-3-215-125 |
| 150        |              | FT-3-329-150 | FT-3-205-150 | FT-3-215-150 |
| 185        |              | FT-3-329-185 | FT-3-205-185 | FT-3-215-185 |
| 240        |              | FT-3-329-240 | FT-3-205-240 | FT-3-215-240 |



#### Folded Filters, 100 Pieces

| Dia. in mm | Grade 131    | Grade 132    | Grade 292    | Grade 292a   |
|------------|--------------|--------------|--------------|--------------|
| 110        | FT-4-351-110 | FT-4-329-110 | FT-4-205-110 | FT-4-215-110 |
| 125        | FT-4-351-125 | FT-4-329-125 | FT-4-205-125 | FT-4-215-125 |
| 150        | FT-4-351-150 | FT-4-329-150 | FT-4-205-150 | FT-4-215-150 |
| 185        | FT-4-351-185 | FT-4-329-185 | FT-4-205-185 | FT-4-215-185 |
| 240        |              | FT-4-329-240 | FT-4-205-240 | FT-4-215-240 |



#### Sheets in 580 × 580 mm, 100 Pieces

| Grade 292       | Grade 292a      |
|-----------------|-----------------|
| FT-2-205-580580 | FT-2-215-580580 |

Other dimensions are available on request



## Filter Papers

### For Qualitative-Technical Analyses

These filter papers are used for routine analyses like clarification, determination of substances, but also as discs with a center hole for technical applications. Grades with a wet burst resistance > 30 kPa are referred to as wet-strengthened and are therefore suitable for pressure or vacuum filtration. They are made of refined pulp and linters with an > 95%  $\alpha$ -cellulose content, are very pure with an ash content between <0.1 to 0.15%. Below you will find an overview of the most commonly used grades.

### Specifications

| Grade | Surface | Weight (g/m <sup>2</sup> ) | Thickness (mm) | Particle Retention (μm) | Filtration (s) | Wet Burst Resistance (kPa) | Properties   |
|-------|---------|----------------------------|----------------|-------------------------|----------------|----------------------------|--|
| 3 hw  | Smooth  | 65                         | 0.14           | 8 to 12                 | 20             | 40                         | Medium fast filtering, filter paper for routine work in the lab  |
| 4 b   | Smooth  | 75                         | 0.15           | 8 to 12                 | 22             | > 15                       | Medium fast filtering, filtration of coarse precipitates, wick paper for seed testing                          |
| 603/N | Crêped  | 75                         | 0.25           | > 15                    | 8              | ≥ 50                       | Fast filtering, filtration of sugar solutions  |
| 6     | Smooth  | 80                         | 0.17           | 10 to 13                | 15             | 30                         | Fast filtering, degassing beer before analysis, clarification of spirits                                       |
| 100/N | Smooth  | 85                         | 0.18           | 6 to 8                  | 30             | 80                         | Medium fast filtering, ash content <0.1%, low potassium and sodium content, determination of the sugar content |
| 5 H/N | Crêped  | 85                         | 0.28           | > 40                    | 3              | ≥ 40                       | Very fast filtering, wide-pore, filtration of essential oils   |
| 3 S/h | Smooth  | 200                        | 0.36           | 5 to 7                  | 55             | 15                         | Medium fast to slow filtering, narrow-pore, re-wet test for diapers  |

### Ordering Information

#### Filter Discs

| Dia. in mm | Grade 3 hw (100 Pieces) | Grade 4 b (100 Pieces) | Grade 603/N (100 Pieces) | Grade 6 (100 Pieces) | Grade 100/N (100 Pieces) | Grade 5 H/N (100 Pieces) | Grade 3 S/h (50 Pieces) |
|------------|-------------------------|------------------------|--------------------------|----------------------|--------------------------|--------------------------|-------------------------|
| 55         | FT-3-303-055            | FT-3-309-055           |                          | FT-3-312-055         | FT-3-328-055             |                          | FT-3-307-055            |
| 70         | FT-3-303-070            | FT-3-309-070           |                          | FT-3-312-070         | FT-3-328-070             |                          |                         |
| 90         | FT-3-303-090            | FT-3-309-090           | FT-3-335-090             | FT-3-312-090         | FT-3-328-090             | FT-3-423-090             | FT-3-307-090            |
| 110        | FT-3-303-110            | FT-3-309-110           | FT-3-335-110             | FT-3-312-110         | FT-3-328-110             |                          | FT-3-307-110            |
| 125        | FT-3-303-125            | FT-3-309-125           | FT-3-335-125             | FT-3-312-125         | FT-3-328-125             | FT-3-423-125             | FT-3-307-125            |
| 150        | FT-3-303-150            | FT-3-309-150           | FT-3-335-150             | FT-3-312-150         | FT-3-328-150             | FT-3-423-150             | FT-3-307-150            |
| 185        | FT-3-303-185            | FT-3-309-185           | FT-3-335-185             | FT-3-312-185         | FT-3-328-185             | FT-3-423-185             | FT-3-307-185            |
| 240        | FT-3-303-240            | FT-3-309-240           | FT-3-335-240             | FT-3-312-240         | FT-3-328-240             | FT-3-423-240             | FT-3-307-240            |

#### Folded Filters, 100 Pieces

| Dia. in mm | Grade 3 hw   | Grade 4 b    | Grade 603/N  | Grade 6      | Grade 100/N  | Grade 5 H/N  |
|------------|--------------|--------------|--------------|--------------|--------------|--------------|
| 125        | FT-4-303-125 | FT-4-309-125 | FT-4-335-125 | FT-4-312-125 |              | FT-4-423-125 |
| 150        | FT-4-303-150 | FT-4-309-150 | FT-4-335-150 | FT-4-312-150 | FT-4-328-150 | FT-4-423-150 |
| 185        | FT-4-303-185 | FT-4-309-185 | FT-4-335-185 | FT-4-312-185 |              | FT-4-423-185 |
| 240        | FT-4-303-240 | FT-4-309-240 | FT-4-335-240 | FT-4-312-240 | FT-4-328-240 | FT-4-423-240 |
| 270        | FT-4-303-270 | FT-4-309-270 | FT-4-335-270 | FT-4-312-270 | FT-4-328-270 | FT-4-423-270 |
| 320        | FT-4-303-320 | FT-4-309-320 | FT-4-335-320 | FT-4-312-320 | FT-4-328-320 | FT-4-423-320 |

#### Sheets in 580 × 580 mm, 100 Pieces

| Grade 3 hw      | Grade 4 b       | Grade 603/N     | Grade 6         | Grade 100/N     | Grade 5 H/N     |
|-----------------|-----------------|-----------------|-----------------|-----------------|-----------------|
| FT-2-303-580580 | FT-2-309-580580 | FT-2-335-580580 | FT-2-312-580580 | FT-2-328-580580 | FT-2-423-580580 |

Other dimensions are available on request

## Glass Microfiber Filters

### Without Binder

Binder-free glass microfiber filters are recommended for analytical and gravimetric analyses and also as prefilters. These filters combine fast flow rates with high load capacity and the retention of very fine particles; they are biologically inert, are resistant to most chemicals and withstand temperatures up to 500°C (grade 550-HA up to 550°C).

### Specifications

| Grade     | Weight (g/m <sup>2</sup> ) | Thickness (mm) | Penetration 0.3 µm (%)* | Particle retention in liquids (µm) | Filtration speed (ml/min) | Fulfills the requirements in EN 872:2005 (weigh loss) |
|-----------|----------------------------|----------------|-------------------------|------------------------------------|---------------------------|---|
| MGA       | 55                         | 0.25           | < 0.002                 | 1.6                                | 510                       | Yes   |
| MGB       | 140                        | 0.70           | < 0.002                 | 1.0                                | 210                       |   |
| MGC       | 52                         | 0.26           | < 0.002                 | 1.2                                | 335                       | Yes   |
| MGD       | 120                        | 0.53           | < 0.1                   | 2.7                                | 920                       |   |
| MGF       | 75                         | 0.38           | ≤ 0.001                 | 0.7                                | 110                       |   |
| MGG       | 64                         | 0.28           | ≤ 0.001                 | 1.5                                | 600                       |   |
| 13440     | 88                         | 0.44           |                         | 0.7                                | 120                       | Yes   |
| MG 160    | 75                         | 0.35           | 0.002                   | 1.2                                | 400                       |   |
| MG 550-HA | 65                         | 0.27           |                         | 1.5                                | 400                       |   |
| MG 169    | 68                         | 0.33           |                         | 1.0                                | 130                       |   |

\* Measurement according to EN 143 (0.3 µm, 15 cm/s, paraffin oil)

### Ordering Information



#### Filter Discs

| Dia. in mm | MGA (100 Pieces) | MG 160 (50 Pieces) | MGB (50 Pieces) | MGC (100 Pieces) | MGD (50 Pieces) |
|------------|------------------|--------------------|-----------------|------------------|-----------------|
| 21         |                  |                    | FT-3-1102-021   |                  |                 |
| 25         | FT-3-1101-025    |                    | FT-3-1102-025   | FT-3-1103-025    | FT-3-1104-025   |
| 37         | FT-3-1101-037    | FT-3-01110-037     |                 |                  |                 |
| 47         | FT-3-1101-047    | FT-3-01110-047     | FT-3-1102-047   | FT-3-1103-047    | FT-3-1104-047   |
| 50         | FT-3-1101-050    | FT-3-01110-050     | FT-3-1102-050   | FT-3-1103-050    | FT-3-1104-050   |
| 55         | FT-3-1101-055    |                    | FT-3-1102-055   | FT-3-1103-055    |                 |
| 70         | FT-3-1101-070    | FT-3-01110-070     | FT-3-1102-070   | FT-3-1103-070    | FT-3-1104-070   |
| 80         | FT-3-1101-080    |                    |                 |                  |                 |
| 90         | FT-3-1101-090    | FT-3-01110-090     | FT-3-1102-090   | FT-3-1103-090    | FT-3-1104-090   |
| 100        | FT-3-1101-100    | FT-3-01110-100     | FT-3-1102-100   | FT-3-1103-100    | FT-3-1104-100   |
| 110        | FT-3-1101-110    | FT-3-01110-110     | FT-3-1102-110   | FT-3-1103-110    | FT-3-1104-110   |
| 125        | FT-3-1101-125    |                    | FT-3-1102-125   | FT-3-1103-125    | FT-3-1104-125   |
| 150        | FT-3-1101-150    |                    | FT-3-1102-150   | FT-3-1103-150    | FT-3-1104-150   |

| Dia. in mm | MGF<br>(100 Pieces) | MGG<br>(100 Pieces) | MG 550-HA<br>(100 Pieces) | 13440*          |
|------------|---------------------|---------------------|---------------------------|-----------------|
| 24         |                     |                     | FT-3-01147-024            |                 |
| 25         | FT-3-1105-025       | FT-3-1106-025       |                           |                 |
| 42         |                     |                     |                           | 13440--42-----Q |
| 44         |                     |                     |                           | 13440--44-----Q |
| 47         | FT-3-1105-047       | FT-3-1106-047       | FT-3-01147-047            | 13440--47-----Q |
| 50         | FT-3-1105-050       | FT-3-1106-050       | FT-3-01147-050            | 13440--50-----Q |
| 55         | FT-3-1105-055       | FT-3-1106-055       | FT-3-01147-055            |                 |
| 70         | FT-3-1105-070       | FT-3-1106-070       | FT-3-01147-070            |                 |
| 90         | FT-3-1105-090       | FT-3-1106-090       | FT-3-01147-090            |                 |
| 100        |                     |                     |                           | 13440-100-----K |
| 110        | FT-3-1105-110       | FT-3-1106-110       | FT-3-01147-110            |                 |
| 125        | FT-3-1105-125       | FT-3-1106-125       | FT-3-01147-125            |                 |
| 130        |                     |                     |                           | 13440-130-----K |
| 150        | FT-3-1105-150       | FT-3-1106-150       |                           |                 |

\* Q = 500 pieces | K = 50 pieces

Other dimensions as well as sheets are available on request



## Glass Microfiber Filters

### With Binder

These filters are mostly used either for monitoring air and gas or as a prefilter. They are manufactured with synthetic binding agents to ensure that the filter has a defined strength. They are mechanically and chemically stable, have a temperature resistance up to 180°C and – depending on the binding agent used – are either hydrophobic or hydrophilic.

### Specifications

| Grade       | Weight (g/m <sup>2</sup> ) | Thickness (mm) | Penetration 0.3 µm (%)* | Pressure drop 5.3 cm/s (Pa) | Binding agent |
|-------------|----------------------------|----------------|-------------------------|-----------------------------|---------------|
| MG 227/1/60 | 60                         | 0.32           | < 0.5                   | 260                         | Hydrophobic   |
| 13430       | 220                        | 1.25           | 0.02                    | 360                         | Hydrophilic   |
| 13400       | 73                         | 0.39           | 0.015                   | 363                         | Hydrophilic   |
| MG 227      | 75                         | 0.40           | < 0.01                  | 350                         | Hydrophobic   |
| MG 161      | 75                         | 0.40           | ≤ 0.002                 | ≤ 580                       | Hydrophilic   |
| MG 400      | 80                         | 0.38           | < 0.001                 | 400                         | Hydrophilic   |
| MG 1387/1   | 90                         | 0.37           | ≤ 0.003                 | 400                         | Hydrophilic   |

\* Tested and classified according to the Standard EN 143

### Ordering Information



#### Filter Discs

| Dia. in mm | MG 227/1/60 (100 Pieces) | 13430**         | 13400**         | MG 227 (100 Pieces) | MG 1387/1 (50 Pieces) |
|------------|--------------------------|-----------------|-----------------|---------------------|-----------------------|
| 13         |                          |                 | 13400--13-----S |                     |                       |
| 20         |                          |                 | 13400--20-----S |                     |                       |
| 25         |                          |                 | 13400--25-----Q |                     |                       |
| 42         |                          |                 | 13400--42-----Q |                     |                       |
| 44         |                          |                 | 13400--44-----Q |                     |                       |
| 45         |                          |                 | 13400--45-----Q |                     | FT-3-01125-045        |
| 47         |                          | 13430--47-----S | 13400--47-----Q | FT-3-01120-047      | FT-3-01125-047        |
| 50         |                          |                 | 13400--50-----Q |                     | FT-3-01125-050        |
| 55         |                          |                 |                 | FT-3-01120-055      | FT-3-01125-055        |
| 80         |                          |                 | 13400--80-----N |                     |                       |
| 100        |                          | 13430-100-----K | 13400-100-----K |                     |                       |
| 110        |                          |                 |                 | FT-3-01120-110      | FT-3-01125-110        |
| 120        |                          |                 | 13400-120-----K |                     |                       |
| 124        |                          | 13430-124-----K | 13400-124-----K |                     |                       |
| 125        |                          |                 |                 |                     | FT-3-01125-125        |
| 127        |                          | 13430-127-----K | 13400-127-----K |                     |                       |
| 130        |                          | 13430-130-----K | 13400-130-----K |                     | FT-3-01125-130        |
| 150        | FT-3-01124-150           |                 |                 |                     |                       |

\*\* K = 50 pieces, N = 100 pieces, Q = 500 pieces, S = 200 pieces  
Other dimensions as well as sheets are available on request

## ■ Quartz Microfiber Filters

The quartz microfiber material of the Sartorius pre-heated filters, grade Q3400, is made of high-purity quartz microfibers without any addition of glass microfibers or binding agents. In addition, the Q3400 filter grade is tempered to remove all chemically combined water and to give the filters excellent weight and dimensional stability. Sartorius filters are especially suitable for emissions monitoring at temperatures of up to 900°C and wherever filters of the highest purity are needed.

### □ Specifications

| Grade | Material   | Weight<br>(g/m <sup>2</sup> ) | Thickness<br>(mm) | Penetration,<br>0.3 µm 15 cm/s* | Temperature<br>Resistance |
|-------|--|-------------------------------|-------------------|---------------------------------|---------------------------|
| Q3400 | 100% Quartz microfiber<br>silicium dioxide (SiO <sub>2</sub> ) | 85                            | 0.43              | <0.002                          | up to 900°C               |

\* Tested and classified according to the Standard EN 143

### □ Ordering Information

#### ○ Filter Discs


| Dia. in mm | Q3400           |
|------------|-----------------|
| 20         | Q3400--20-----G |
| 25         | Q3400--25-----G |
| 37         | Q3400--37-----G |
| 45         | Q3400--45-----G |
| 47         | Q3400--47-----G |
| 50         | Q3400--50-----G |
| 82         | Q3400--82-----N |
| 90         | Q3400--90-----N |
| 142        | Q3400-142-----K |
| 150        | Q3400-150-----K |

\* G = 25 pieces, K = 50 pieces, N = 100 pieces


Other dimensions as well as sheets are available on request



## Membrane Filtration – Quick Selection Guide



## Cellulose Nitrate



Cellulose nitrate membrane filters are indicated for many general laboratory applications where a membrane with a high non-specific adsorption is suitable. They are hydrophilic, have high flow rates thanks to their symmetric structure and are compatible with aqueous solutions

(pH 4 to 8), hydrocarbons and several other organic solvents. The cellulose nitrate membranes are available in different pore sizes from 0.2  $\mu\text{m}$  to 8  $\mu\text{m}$ .

### Specifications

| Type  | Pore Size ( $\mu\text{m}$ ) | Thickness ( $\mu\text{m}$ ) | Bubble Point (bar) | Water Flow Rate (ml/min/cm <sup>2</sup> /bar) | Burst Pressure (bar) |
|-------|-----------------------------|-----------------------------|--------------------|---|----------------------|
| 11327 | 0.2                         | 130                         | 4.2                | 25  | $\geq 0.35$          |
| 11306 | 0.45                        | 130                         | 2.4                | 70  | $\geq 0.3$           |
| 11305 | 0.65                        | 130                         | 2                  | 130   | $\geq 0.25$          |
| 11304 | 0.8                         | 130                         | 1.4                | 200   | $\geq 0.2$           |
| 11303 | 1.2                         | 130                         | 1                  | 200   | $\geq 0.2$           |
| 11302 | 3                           | 130                         | 0.5                | 430   | $\geq 0.2$           |
| 11342 | 5                           | 130                         | 0.5                | 570   | $\geq 0.15$          |
| 11301 | 8                           | 130                         | 0.3                | 750   | $\geq 0.1$           |

### Ordering Information


#### Filter Discs

| Dia. in mm | 11301 (8 $\mu\text{m}$ )*    | 11302 (3 $\mu\text{m}$ )*    | 11303 (1.2 $\mu\text{m}$ )* | 11304 (0.8 $\mu\text{m}$ )* |
|------------|------------------------------|------------------------------|-----------------------------|-----------------------------|
| 13         | 11301--13-----N              | 11302--13-----N              | 11303--13-----N             | 11304--13-----N             |
| 20         |                              |                              |                             | 11304--20-----N             |
| 25         | 11301--25-----N              | 11302--25-----N              | 11303--25-----N             | 11304--25-----N             |
| 30         |                              |                              |                             | 11304--30-----N             |
| 37         | 11301--37-----N              |                              |                             | 11304--37-----N             |
| 47         | 11301--47-----N              | 11302--47-----N              | 11303--47-----N             | 11304--47-----N             |
| 50         | 11301--50-----N              | 11302--50-----N              | 11303--50-----N             | 11304--50-----N             |
| 68         |                              | 11302--68-----G              |                             |                             |
| 70         | 11301--70-----G              |                              |                             |                             |
| 90         |                              | 11302--90-----G              | 11303--90-----G             | 11304--90-----G             |
| 100        | 11301-100-----N              | 11302-100-----G              | 11303-100-----G             | 11304-100-----G             |
| Dia. in mm | 11305 (0.65 $\mu\text{m}$ )* | 11306 (0.45 $\mu\text{m}$ )* | 11327 (0.2 $\mu\text{m}$ )* | 11342 (5 $\mu\text{m}$ )*   |
| 13         | 11305--13-----N              | 11306--13-----N              | 11327--13-----N             | 11342--13-----N             |
| 20         |                              | 11306--20-----N              |                             |                             |
| 25         | 11305--25-----N              | 11306--25-----N              | 11327--25-----N             | 11342--25-----N             |
| 30         |                              | 11306--30-----N              |                             |                             |
| 37         |                              | 11306--37-----N              |                             |                             |
| 47         | 11305--47-----N              | 11306--47-----N              | 11327--47-----N             | 11342--47-----N             |
| 50         | 11305--50-----N              | 11306--50-----N              |                             | 11342--50-----N             |
| 55         |                              | 11306--55-----N              |                             |                             |
| 85         |                              | 11306--85-----N              |                             |                             |
| 90         |                              | 11306--90-----N              |                             | 11342--90-----G             |
| 100        | 11305-100-----N              | 11306-100-----N              |                             | 11342-100-----G             |
| 110        |                              | 11306-110-----N              |                             |                             |

\* G = 25 pieces, N = 100 pieces

Other dimensions and packaging units are available on request

## Cellulose Acetate



Cellulose acetate membranes combine high flow rates and thermal stability with very low adsorption characteristics, and are therefore excellently suited for use in pressure filtration devices. They are hydrophilic, have high flow rates thanks to their symmetric structure and are compatible with aqueous solutions

(pH 4–8), oils, alcohols and other organic solvents. The 0.2 µm membrane is the filter of choice for sterile filtration of aqueous solutions, such as nutrient media, buffers and sera. The cellulose acetate membranes are available in different pore sizes from 0.2 to 5 µm.

### Specifications

| Type  | Pore Size (µm) | Thickness (µm) | Bubble Point (bar) | Water Flow Rate (ml/min/cm <sup>2</sup> /bar) | Burst Pressure (bar) |
|-------|----------------|----------------|--------------------|---|----------------------|
| 11107 | 0.2            | 120            | 2.9                | 24  | 0.8                  |
| 11106 | 0.45           | 120            | 1.9                | 69  | 0.7                  |
| 11105 | 0.65           | 120            | 1.5                | 115   | 0.7                  |
| 11104 | 0.8            | 120            | 1                  | 200   | 0.5                  |
| 12303 | 1.2            | 140            | 0.8                | 320   | 0.4                  |
| 12342 | 5              | 140            | 0.4                | 570   | 0.25                 |

### Ordering Information



#### Filter Discs


| Dia. in mm | 11104 (0.8 µm)* | 11105 (0.65 µm)* | 11106 (0.45 µm)* | 11107 (0.2 µm)*  | 12303 (1.2 µm)*  | 12342 (5 µm)*   |
|------------|-----------------|------------------|------------------|------------------|------------------|-----------------|
| 13         | 11104--13-----N |                  | 11106--13-----N  | 11107--13-----N  |                  |                 |
| 25         | 11104--25-----N | 11105--25-----N  | 11106--25-----N  | 11107--25-----N  | 12303--25-----N  | 12342--25-----N |
| 30         |                 |                  | 11106--30-----N  | 11107--30-----N  |                  |                 |
| 37         |                 |                  | 11106--37-----N  |                  |                  |                 |
| 45         |                 |                  |                  |                  |                  |                 |
| 47         | 11104--47-----N | 11105--47-----N  | 11106--47-----N  | 11107--47-----N  | 12303--47-----N  | 12342--47-----N |
| 50         | 11104--50-----N | 11105--50-----N  | 11106--50-----N  | 11107--50-----N  | 12303--50-----N  |                 |
| 70         |                 |                  |                  |                  |                  |                 |
| 85         |                 |                  | 11106--85-----N  |                  |                  |                 |
| 90         | 11104--90-----N | 11105--90-----G  | 11106--90-----G  | 11107--90-----G  |                  |                 |
| 100        |                 |                  | 11106--100-----N | 11107--100-----N | 12303--100-----G |                 |
| 110        |                 |                  | 11106--110-----N |                  |                  |                 |

\* G = 25 pieces, N = 100 pieces

Other dimensions and packaging units are available on request



## ■ Regenerated Cellulose



The very low adsorption membranes are hydrophilic, solvent-resistant (pH 3–12) and therefore suited for the particle removal from solvents. The membrane is asymmetric and reinforced with nonwoven cellulose. They are available in two pore sizes: 0.45  $\mu\text{m}$  and 0.2  $\mu\text{m}$ .

### □ Specifications

| Type  | Pore Size ( $\mu\text{m}$ ) | Thickness ( $\mu\text{m}$ ) | Bubble Point (bar) | Water Flow Rate (ml/min/cm <sup>2</sup> /bar) |
|-------|-----------------------------|-----------------------------|--------------------|---|
| 18407 | 0.2                         | 170                         | 4.4                | 15  |
| 18406 | 0.45                        | 170                         | 2.9                | 30  |

### □ Ordering Information




#### Filter Discs

| Dia. in mm | 18406 (0.45 $\mu\text{m}$ )* | 18407 (0.2 $\mu\text{m}$ )* |
|------------|------------------------------|-----------------------------|
| 13         | 18406--13-----N              | 18407--13-----N             |
| 25         | 18406--25-----N              | 18407--25-----N             |
| 47         | 18406--47-----N              | 18407--47-----N             |
| 50         | 18406--50-----N              | 18407--50-----N             |
| 100        | 18406-100-----G              | 18407-100-----G             |

\* G = 25 pieces, N = 100 pieces

Other dimensions and packaging units are available on request

## ■ Polyethersulfone



Polyethersulfone (PES) membrane filters are hydrophilic, have high flow rates thanks to their symmetric structure, have a low non-specific protein adsorption and are chemically resistant over a pH range of

1 – 14. They are therefore recommended for the filtration of aqueous solutions as well for protein filtration. Furthermore, the low level of extractables makes them suitable for environmental analysis.

### □ Specifications

| Type    | Pore Size (μm) | Thickness (μm) | Bubble Point (bar) | Water Flow Rate (ml/min/cm <sup>2</sup> /bar) | Burst Pressure (bar) |
|---------|----------------|----------------|--------------------|---|----------------------|
| 15458   | 0.1            | 150            | 3.8                | 10  | ≥ 0.6                |
| 15407MI | 0.2            | 150            | 3.5                | 25  | ≥ 0.5                |
| 15406   | 0.45           | 150            | 2.6                | 46  | ≥ 0.5                |

### □ Ordering Information




#### Filter Discs

| Dia. in mm | 15406 (0.45 μm)* | 15407MI (0.2 μm)* | 15458 (0.1 μm)* |
|------------|------------------|-------------------|-----------------|
| 25         | 15406--25-----N  | 15407--25----MIN  | 15458--25-----N |
| 47         | 15406--47-----N  | 15407--47----MIN  | 15458--47-----N |
| 50         | 15406--50-----N  | 15407--50----MIN  | 15458--50-----N |
| 90         |                  | 15407--90----MIK  |                 |

\*K = 50 pieces, N = 100 pieces  
Other dimensions are available on request

## ■ Polyamide



Polyamide membrane filters are hydrophilic and chemically resistant to alkaline solutions and organic solvents. They are therefore recommended for the particle removal from aqueous solutions and

solvents for analytical determination such as HPLC, for the sterile filtration of these liquids as well as for applications where a membrane with a relatively high non-specific adsorption is suitable.

### □ Specifications

| Type  | Pore Size (μm) | Thickness (μm) | Bubble Point (bar) | Water Flow Rate (ml/min/cm <sup>2</sup> /bar) | Burst Pressure (bar) |
|-------|----------------|----------------|--------------------|---|----------------------|
| 25007 | 0.2            | 115            | 3.2                | 15  | ≥ 0.25               |
| 25006 | 0.45           | 115            | 2.3                | 35  | ≥ 0.23               |

### □ Ordering Information




#### Filter Discs

| Dia. in mm | 25006 (0.45 μm)* | 25007 (0.2 μm)* |
|------------|------------------|-----------------|
| 13         | 25006--13-----N  | 25007--13-----N |
| 25         | 25006--25-----N  | 25007--25-----N |
| 47         | 25006--47-----N  | 25007--47-----N |
| 50         | 25006--50-----N  | 25007--50-----N |
| 55         |                  | 25007--55-----N |
| 90         | 25006--90-----G  | 25007--90-----G |

\* G = 25 pieces, N = 100 pieces

Other dimensions and packaging units are available on request

## Hydrophobic PTFE



The main application of these membrane filters is the filtration of air, gases or chemicals. They are made of PTFE (polytetra-fluorethylene) only and are therefore permanently hydrophobic. Unlike other (hydrophilic) filter types, they are not wetted by air humidity, allowing unhindered passage of air at low differential pressures as well. PTFE

membrane filters have an excellent chemical compatibility (pH 1 to 14), so that they are also used for the filtration of solvents and acids, to which other filter types are not resistant. Due to their hydrophobic characteristics, they must be pre-wetted with ethanol or methanol before the filtration of aqueous media.

### Specifications

| Type  | Pore Size (µm) | Thickness (µm) | Bubble Point (bar) | Isopropanol Flow Rate (ml/min/cm <sup>2</sup> /bar) |
|-------|----------------|----------------|--------------------|---|
| 11807 | 0.2            | 65             | 1.4                | 11  |
| 11806 | 0.45           | 80             | 0.9                | 20  |
| 11803 | 1.2            | 100            | 0.45               | 80  |
| 11842 | 5              | 100            | 0.10               | 250   |

### Ordering Information




#### Filter Discs

| Dia. in mm | 11803 (1.2 µm)* | 11806 (0.45 µm)* | 11807 (0.2 µm)* | 11842 (5 µm)*   |
|------------|-----------------|------------------|-----------------|-----------------|
| 13         | 11803--13-----N | 11806--13-----N  | 11807--13-----N |                 |
| 25         | 11803--25-----N | 11806--25-----N  | 11807--25-----N | 11842--25-----N |
| 37         |                 | 11806--37-----N  |                 |                 |
| 42         |                 |                  |                 | 11842--42-----N |
| 47         | 11803--47-----N | 11806--47-----N  | 11807--47-----N | 11842--47-----N |
| 50         | 11803--50-----N | 11806--50-----N  | 11807--50-----N | 11842--50-----N |
| 90         | 11803--90-----G | 11806--90-----G  | 11807--90-----G |                 |
| 100        | 11803-100-----G | 11806-100-----G  | 11807-100-----G | 11842-100-----G |

\* G= 25 pieces, N= 100 pieces

Other dimensions and packaging units are available on request

## ■ Polycarbonate Track-Etched



Those white and hydrophilic polycarbonate track-etched membranes are manufactured from high grade polycarbonate film using track-etch technology. Their capillary pore structure is uniform and precise, with a narrow pore size distribution to retain particles on their surface. Track-etched membranes are an excellent choice for accurate fractionation of particulates because of their precise pore size.

Track-etch technology offers the user distinct performance advantages when excellent surface capture and high sample visibility are required. Their main applications are particulate analysis, epifluorescence microscopy, fluid clarification, cytology, cell biology, bioassays, water microbiology and environmental analysis.

### □ Specifications

| Type  | Pore Size (μm) | Thickness (μm) | Bubble Point (bar) | Water Flow Rate (ml/min/cm <sup>2</sup> /0.7 bar) | Burst Pressure (bar) |
|-------|----------------|----------------|--------------------|---|----------------------|
| 23007 | 0.2            | 25             | 4.8                | ≥ 10  | ≥ 0.7                |
| 23006 | 0.4            | 25             | 2.5                | ≥ 30  | ≥ 0.7                |

### □ Ordering Information



#### Filter Discs, 100 Pieces

| Dia. in mm | 23006 (0.4 μm)* | 23007 (0.2 μm)* |
|------------|-----------------|-----------------|
| 25         | 23006--25-----N | 23007--25-----N |
| 47         | 23006--47-----N | 23007--47-----N |
| 50         |                 | 23007--50-----N |

Other dimensions are available on request

# Blotting | Chromatography Papers



These papers are made of cotton linters only with  $\alpha$ -cellulose content of > 98%. These highly pure papers are not only ideal for blotting and chromatography, but also for a wide range of absorption applications like those common in the life sciences and diagnostics. Below you will find an overview of the most commonly used grades.

## Specifications

| Grade  | Weight (g/m <sup>2</sup> ) | Thickness (mm) | Capillary Rise (mm/30 min) | Capillary Rise (mm/10 min) | Properties  |
|--------|----------------------------|----------------|----------------------------|----------------------------|---|
| FN 4   | 125                        | 0.24           | 95                         |                            | Chromatography paper, ash content < 0.04%                                   |
| FN 7   | 150                        | 0.32           | 145                        |                            | Chromatography paper, ash content < 0.04%                                   |
| FN 30  | 320                        | 0.90           | 240                        |                            | Chromatography paper, ash content < 0.04%, paper for antibiotic test strips |
| FN 100 | 195                        | 0.35           | 115                        | 70                         | The most commonly used chromatography and blotting paper                    |
| BF 3   | 330                        | 0.76           | 30                         | 130                        | Blotting paper to increase and maintain the transport of liquids            |

## Ordering Information



**Sheets in 580 × 600 mm**

| Grade FN 4<br>(100 Sheets) | Grade FN 7<br>(50 Sheets) | Grade FN 30<br>(25 Sheets) | Grade FN 100<br>(50 Sheets) | Grade BF 3<br>(50 Sheets) |
|----------------------------|---------------------------|----------------------------|-----------------------------|---------------------------|
| FT-2-504-580600N           | FT-2-507-580600K          | FT-2-526-580600G           | FT-2-527-580600K            | FT-2-520-580600K          |

Other dimensions and shapes are available on request

## ■ Nitrocellulose Membrane for Blotting



Sartorius nitrocellulose membranes are available in two pore sizes, 0.22  $\mu\text{m}$  and 0.45  $\mu\text{m}$ . Both versions combine the advantages of high protein binding capacity with low background and high membrane stability, which ensures easy handling. Due to its large surface area, the 0.22  $\mu\text{m}$

membrane version is recommended for small proteins. Sartorius blotting membranes are ideal for western blotting, DNA blotting as well as dot or slot blots. They have been optimized for all protein blotting systems, such as electrotransfer, semi-dry or simple capillary blotting.

### □ Specifications

|                              | 0.22 $\mu\text{m}$              | 0.45 $\mu\text{m}$              |
|------------------------------|---------------------------------|---------------------------------|
| Material                     | Cellulose nitrate               | Cellulose nitrate               |
| Thickness                    | 120 $\mu\text{m}$               | 130 $\mu\text{m}$               |
| Water flow rate              | 27 ml/(min.cm <sup>2</sup> bar) | 70 ml/(min.cm <sup>2</sup> bar) |
| Bubble point                 | 4.4 bar                         | 2.4 bar                         |
| Wettability in water         | $\leq 1$ s                      | $\leq 1$ s                      |
| Extractable content in water | $\leq 1\%$                      | $\leq 1\%$                      |
| Burst pressure               | 0.8 bar                         | 0.2 bar                         |
| Binding capacity for IgG     | 200 $\mu\text{g}/\text{cm}^2$   | 200 $\mu\text{g}/\text{cm}^2$   |

### □ Ordering Information

|                       | Roll Size          | Order No.      |
|-----------------------|--------------------|----------------|
| NC 0.22 $\mu\text{m}$ | 30 cm $\times$ 3 m | 11327-----41BL |
| NC 0.45 $\mu\text{m}$ | 30 cm $\times$ 3 m | 11306-----41BL |

All indicated data to be understood as typical average values

## ■ Re-usable, 13 mm Syringe Filter Holders

For the Ultracleaning of Small Volumes Up to About 10 ml



### PTFE Holder for Solvents and Chemicals

Made completely of PTFE, this holder is unaffected by chemicals and contains no trace elements which could be released into the liquid being filtered. It is therefore extremely well suited for particle removal from samples and reagents for analytical methods, such as NMR samples. Other

benefits of this application are the low hold-up volume, the easy cleaning and the drying at a temperature of 180°C. The construction of the holder ensures leak proof sealing without a sealing ring, and avoids twisting of the membrane filter when the top is tightened onto the base.

### □ Specifications

|                          |  |
|--------------------------|--|
| Connectors               | Female Luer Lock inlet, luer slip outlet                               |
| Chemical compatibility   | As for PTFE  |
| Filtration area          | 0.5 cm <sup>2</sup>  |
| Materials                | PTFE top and bottom parts  |
| Max. operating pressure  | 5 bar   500 kPa   72.5 psi   |
| Membrane filter diameter | 13 mm  |
| Sterilization            | By autoclaving (max. 134°C)<br>or by dry heat (max. 180°C)             |
| Hold-up volume           | Less than 0.03 ml after overcoming the<br>bubble point (0.3 ml before) |

### □ Ordering Information

| Description                      | Order No. |
|----------------------------------|-----------|
| 13 mm PTFE Syringe Filter Holder | 16574     |



### Polycarbonate Holder for Aqueous Solutions

This inexpensive filter holder is made of clear, autoclavable polycarbonate. The silicone gasket enables a leak-free

filtration at pressures of up to 7 bar by simply screwing it together manually. Filter supports in the top and bottom parts allow filtration in either direction.

### □ Specifications

|                          |   |
|--------------------------|---|
| Connectors               | Female Luer Lock inlet, luer slip outlet                              |
| Chemical compatibility   | As for polycarbonate and silicone                                     |
| Filtration area          | 0.5 cm <sup>2</sup>   |
| Materials                | Polycarbonate top and bottom part,<br>silicone gasket                 |
| Max. operating pressure  | 7 bar   700 kPa   101.5 psi   |
| Membrane filter diameter | 13 mm   |
| Sterilization            | By autoclaving at 121°C   |
| Hold-up volume           | Less than 0.2 ml after overcoming the<br>bubble point (0.3 ml before) |

### □ Ordering Information

| Description   | Order No.   |
|---|-------------|
| 13 mm Polycarbonate Syringe Filter Holder, pack of 12 | 16514-----E |
| Silicon gasket, 10 × 14.9 × 0.5 mm, pack of 10        | 6980569     |



## ■ Re-usable 25 mm Syringe Filter Holders

For the Ultracleaning and Sterilizing Filtration of Volumes of Up to About 100 ml



### Stainless Steel Holder for Solvents and Chemicals

Made of stainless steel, this holder is heat-resistant, and the chemical compatibility depends only on the inserted filter type.

The top part can easily be mounted on the bottom part using the enclosed tightening tool. Filter supports in the top and bottom parts allow filtration in either direction.

### □ Specifications

|                          |   |
|--------------------------|---|
| Connectors               | Female Luer Lock inlet, luer slip outlet                              |
| Chemical compatibility   | As for stainless steel  |
| Filtration area          | 3 cm <sup>2</sup>   |
| Materials                | Stainless steel (1.4305) top and bottom parts                         |
| Max. operating pressure  | 7 bar   700 kPa   101.5 psi   |
| Membrane filter diameter | 25 mm   |
| Sterilization            | By autoclaving (max. 134°C)<br>or by dry heat (max. 180°C)            |
| Hold-up volume           | Less than 0.1 ml after overcoming the<br>bubble point (0.3 ml before) |

### □ Ordering Information

| Description                   | Order No. |
|-------------------------------|-----------|
| 25 mm Stainless Steel Holder  | 16214     |
| Tightening tool, Polyman 24/5 | 6980595   |



### Polycarbonate Holder for Aqueous Solutions

This inexpensive filter holder is made of clear, autoclavable polycarbonate. The silicone gasket enables a leak-free filtration

at pressures of up to 7 bar by simply screwing it together manually. Filter supports in the top and bottom parts allow filtration in either direction.

### □ Specifications

|                          |   |
|--------------------------|---|
| Connectors               | Female Luer Lock inlet, luer slip outlet                              |
| Chemical compatibility   | As for polycarbonate and silicone                                     |
| Filtration area          | 3 cm <sup>2</sup>   |
| Materials                | Polycarbonate top and bottom parts,<br>silicone gasket                |
| Max. operating pressure  | 7 bar   700 kPa   101.5 psi   |
| Membrane filter diameter | 25 mm   |
| Sterilization            | By autoclaving at 121°C   |
| Hold-up volume           | Less than 0.3 ml after overcoming the<br>bubble point (0.6 ml before) |

### □ Ordering Information

| Description   | Order No.   |
|---|-------------|
| 25 mm Polycarbonate Syringe Filter Holder, pack of 12 | 16517-----E |
| Silicone gasket, 20.5 × 26.5 × 0.5 mm, pack of 10     | 6980570     |

# 25 mm Glass Vacuum Filter Holder

For Hybridization Tests, Particle Testing and Clarification



This filter holder is available in two versions differing from each other only in the type of the filter support. The filter with glass frit ensures uniform distribution of retained particles and is therefore recommended when the residue on the filter surface is of interest. Because it is easy to clean, the device with the PTFE-coated screen support is preferable when the filtrate is required, or when liquids difficult to remove from glass frits must be examined. The PTFE ring, which holds the

glass frit and the screen support, allows for the autoclaving of the devices with a filter in position and protects the edge of the glass frit from breakage and potential leakage. It has a rim around the upper edge to simplify the positioning of the membrane filter when inserted and a silicone O-ring in the underside for a leak-proof seal on the filtrate side. The funnel-shaped top part simplifies filling in the sample.


## Specifications

|                                   |   |
|-----------------------------------|---|
| Outlet spout                      | 12 mm dia.  |
| Parts and materials               | Borosilicate glass funnel and base<br>PTFE   glass filter support (type 16306) or<br>PTFE   stainless steel filter support, coated with<br>PTFE (type 16315)<br>Silicone O-ring 25 × 3 mm<br>Anodized Aluminium clamp |
| Chemical compatibility            | As for glass, PTFE and silicone.<br>The silicone O-ring can be replaced by a<br>fluoroelastomer O-ring (order no. 00118)  |
| Funnel capacity                   | 30 ml   |
| Filtration area                   | 3 cm <sup>2</sup>   |
| Max. operating pressure           | Only for vacuum   |
| Suitable membrane filter diameter | 25 mm (or 24 mm)  |
| Sterilization                     | By autoclaving (max. 134°C) or<br>by dry heat (max. 180°C)  |

## Ordering Information

| Description   | Order No. |
|---|-----------|
| Glass vacuum filtration holder<br>for 25 mm (or 24 mm) membrane filter,<br>with glass frit filter support         | 16306     |
| Glass vacuum filtration holder<br>for 25 mm (or 24 mm) membrane filter,<br>with PTFE-coated screen filter support | 16315     |

Replacement parts are shown in the diagram.



## 50 mm Glass Vacuum Filter Holder

For Particle Testing or Clarification and Sterile Filtration



This filter holder is available in two versions differing from each other only in the type of the filter support. The device with glass frit ensures uniform distribution of retained particles and is therefore recommended, when the residue on the filter surface is of interest. Because it is easy to clean, the device with the PTFE-coated screen support is preferable when the filtrate is required, or when liquids difficult to remove from

glass frits must be examined. The PTFE ring, which holds the glass frit and the screen support, allows the autoclaving of the devices with a filter in position and protects the edge of the glass frit from breakage and potential leakage. It has a rim around the upper edge to simplify the positioning of the membrane filter when inserted and a silicone O-ring in the underside for a leak-proof seal on the filtrate side.




### Specifications

|                                   |  |
|-----------------------------------|--|
| Outlet spouts                     | 15 mm dia.   |
| Parts and materials               | Borosilicate glass funnel and base<br>Silicone caoutchouc lid<br>PTFE   glass filter support (type 16307) or<br>PTFE   stainless steel filter support, coated with<br>PTFE (type 16316)<br>Silicone O-ring 45 x 3 mm<br>Anodized Aluminium clamp |
| Chemical compatibility            | As for glass, PTFE and silicone.<br>The silicone O-ring can be replaced by a<br>fluoroelastomer O-ring (order no. 00124).  |
| Funnel capacity                   | 250 ml   |
| Filtration area                   | 12.5 cm <sup>2</sup>   |
| Max. operating pressure           | Only for vacuum  |
| Suitable membrane filter diameter | 50 mm (or 47 mm)   |
| Sterilization                     | By autoclaving (max. 134°C) or<br>by dry heat (max. 180°C)   |

### Ordering Information

| Description   | Order No. |
|---|-----------|
| Glass vacuum filtration holder<br>for 50 mm (or 47 mm) membrane filter,<br>with glass frit filter support         | 16307     |
| Glass vacuum filtration holder<br>for 50 mm (or 47 mm) membrane filter,<br>with PTFE-coated screen filter support | 16316     |

Replacement parts are shown in the diagram.



# All-Glass Vacuum Filter Holder

For Analytical Determinations, Particle Removal from Solvents



All areas, where liquid and device can come into direct contact, are made of glass or PTFE. The device, in combination with solvent-resistant, hydrophilic RC-membranes, is therefore ideal for ultracleaning and degassing solvents and solvent mixtures for HPLC, GC and AA. Convenience of handling is ensured by several beneficial features. A 6 mm wide

non-ground rim above the ground glass neck of the suction flask prevents the filtrate from contacting grease on the ground glass surface and so avoids its contamination while being poured out of the flask. The hose nipple connector is made of polypropylene for safe connection of the vacuum hose. The filtrate outlet spout ends well below the entrance to this hose nipple.




## Specifications

|                                   |   |
|-----------------------------------|---|
| Parts and materials               | Borosilicate glass funnel, base and flask, sintered glass frit in a PTFE ring and fluoroelastomer O-ring (45 × 3 mm) underneath, anodized aluminium clamp |
| Chemical compatibility            | As for glass and PTFE   |
| Funnel capacity                   | 250 ml  |
| Capacity of the filtrate flask    | 1 liter   |
| Filtration area                   | 12.5 cm <sup>2</sup>  |
| Max. operating pressure           | Only for vacuum   |
| Suitable membrane filter diameter | 50 mm (or 47 mm), 40 or 42 mm prefilter   |
| Sterilization (without connector) | By autoclaving (max. 134°C) or by dry heat (max. 180°C)   |

## Ordering Information

| Description  | Order No. |
|--|-----------|
| All-glass vacuum filter holder for 50 mm (or 47 mm) membrane filter, with vacuum-resistant flask, capacity 1 liter | 16309     |

Replacement parts are shown in the diagram.



## ■ Polycarbonate In-Line Filter Holder

For the Filtration of Liter Volumes of Aqueous Solutions



This holder is made of stable, autoclavable polycarbonate. This practical holder is suitable for many simple laboratory filtrations. It can be connected to a peristaltic pump or a pressure container. The bell-shaped base protects the filtrate from repeated contamination while flowing in a receiver. The holder is characterized

by an excellent resistance to pressure and density setting by simple hand-tightening. The transparent top part allows the visual control of the correct fit of the O-ring. The hose nipples can be replaced by luer connectors to use it as a large area syringe filter holder.




### □ Specifications

|                                   |  |
|-----------------------------------|--|
| Chemical compatibility            | As for polycarbonate, polypropylene and silicone   |
| Filtration area                   | 12.5 cm <sup>2</sup>   |
| Weight                            | 83 g   |
| Threads for connectors            | M 12 × 1 female thread   |
| Materials                         | Polycarbonate top part, base part and hose nipple, polypropylene filter support, silicone O-ring (40 × 5 mm)   |
| Max. operating pressure           | 7 bar   700 kPa   101.5 psi  |
| Suitable membrane filter diameter | 50 mm (40 or 42 mm prefilter)  |
| Sterilization                     | By autoclaving at 121°C<br>The material withstands repeated cycles, provided aggressive cleaning agents are completely washed off and that the boiler water does not contain anti-corrosive or anti-scaling additives. |

### □ Ordering Information

| Description   | Order No.   |
|---|-------------|
| Polycarbonate in-line filter holder for 50 mm membrane filter, pack of 5. | 16508-----B |

Replacement parts are shown in the diagram.



■ 25 mm Stainless Steel Filter Holder

For In-Line Filtration



The G $\frac{1}{4}$  connection threads with density barrel, guarantee leak-proof sealing of the hose nipple and the holder without sealing rings. Other connectors, available as accessories, fit the holder onto reducing

valves or pumps with G $\frac{1}{4}$  female thread (order no. 01030) or G $\frac{3}{8}$  female thread order no. 01029) or onto pressure tanks with G $\frac{3}{8}$  male thread (order no. 00177).


□ Specifications

|                          |   |
|--------------------------|---|
| Connectors               | Hose nipples DN10   |
| Filtration area          | 3 cm <sup>2</sup>   |
| Weight                   | ca. 170 g   |
| Materials                | Stainless steel, except silicone<br>O-ring (21 × 2 mm) and aluminium closing ring |
| Max. operating pressure  | 5 bar   500 kPa   72.5 psi  |
| Suitable membrane filter | 25 mm (20 mm prefilter for the filtration of liquids only)                        |
| Sterilization            | By autoclaving (max. 134°C) or<br>by dry heat (max. 180°C)                        |

□ Ordering Information

| Description  | Order No. |
|--|-----------|
| Stainless steel pressure filter holder for 25 mm dia. membrane filter. | 16251     |

Replacement parts are shown in the diagram.



## 47 mm Stainless Steel Filter Holder

For In-Line Filtration



The filter holder is suitable for a pressure of up to 20 bar. The inlet side valve is convenient for the intermittent run-off of waste water. Other connectors, available as accessories, fit the holder onto reducing

valves or pumps with G $\frac{3}{8}$  female thread (order no. 17089) or onto pressure tanks with G $\frac{3}{8}$  male thread (order no. 17069) or on taps with G $\frac{3}{4}$  male thread (order no. 17068).

### Specifications

|                          |   |
|--------------------------|---|
| Connectors               | Hose nipples DN10   |
| Connection thread        | M12 × 1   |
| Filtration area          | 13 cm <sup>2</sup>  |
| Weight                   | ca. 490 g   |
| Materials                | Stainless steel, except silicone O-ring (42 × 3 mm), PTFE and fluoroelastomer valve seals |
| Max. operating pressure  | 20 bar   2,000 kPa   290 psi  |
| Suitable membrane filter | 47 mm (40 or 42 mm prefilter)   |
| Sterilization            | By autoclaving (max. 134°C) or by dry heat (max. 180°C)                                   |

### Ordering Information

| Description   | Order No.     |
|---|---------------|
| Stainless steel filter holder for 47 mm membrane filter (with adapter M12 × 1 male thread to hose barb DN10, Mat. 316, ref. 6980801) – Replacement parts are shown in the diagram | 16254         |
| Stainless steel filter holder for 47 mm membrane filter (with adapter M12 × 1 male thread to hose barb DN 4 to 5, Mat. 316, ref. 6981132)   | 16278         |
| Stainless steel back pressure screen  | 6980721-----1 |
| Stainless steel filter support screen   | 6980180-----1 |
| Stainless steel underdrain screen   | 00181         |
| Stainless steel connector M12 × 1 male thread to hose barb DN 4-5   | 6981132       |
| Adapter Quick connect nipple length 60 mm male part to male thread M12 × 1, Mat 316   | 17090-----1   |




Diagram for 16254

# Stainless Steel Pressure Filter Holder

For the Filtration of Up to 5 Liter Volumes




A practical filter holder for many laboratory filtrations. It can be attached to a tripod with the help of a steel rod which can be screwed in. The hose nipple is screwed into the side of the top part, leaving room for a large filling opening. This makes pouring in the sample easier, and the sample can be refilled without removing the tube connection to the pressure source. Leak-proof sealing is achieved by

hand-tightening the closing ring. For the filtration of small volumes (up to about 200 ml of soil samples or viscous liquids, such as oils), the holder is connected directly to a pressure source. For the filtration of up to 5 liter volumes of relatively easily filterable liquids (e.g. buffer solutions, solutions for cell counters and tissue culture solutions), it is used in combination with a pressure tank.



## Specifications

|                                   |   |
|-----------------------------------|---|
| Chemical compatibility            | As for stainless steel, PTFE and silicone.<br>If required, the silicone O-ring in the filter support can be replaced by a fluoroelastomer O-ring 00179 or a PTFE O-ring 17038 (by reducing the max. operating pressure to 4 bar   58 psi); the silicone O-ring in the top part can be replaced by a fluoroelastomer O-ring 17145. |
| Filtration area                   | 13 cm <sup>2</sup>  |
| Weight                            | 960 g   |
| Threads for connectors            | M 12 × 1 female thread  |
| Materials                         | Top part, barrel, base part, corrugated iron, closing ring, closure cap, back pressure screen and stainless steel hose nipples 1.4401 (AISI 316), PTFE-coated stainless steel filter support, silicone O-rings, 41 × 2 mm (top part) and 42 × 3 mm (filter support), PTFE-sealing (cap).  |
| Max. operating pressure           | 10 bar   1,000 kPa   145 psi  |
| Suitable membrane filter diameter | 47 mm (40 or 42 mm prefilter)   |
| Sterilization                     | By autoclaving (max 134°C) or by dry heat (180°C)   |



## Ordering Information

| Description                            | Order No. |
|--|-----------|
| Stainless steel pressure filter holder | 16249     |

### Replacement Parts

| Description                                      | Order No. |
|--|-----------|
| Fluoroelastomer O-ring, 42 × 3 mm                | 00179     |
| PTFE O-ring, 42 × 3 mm                           | 17038     |
| Fluoroelastomer O-ring for upper part, 41 × 2 mm | 17145     |

Other replacement parts are shown in the diagram or on page 107.



## ■ Chemical-resistant PTFE Filter Holder

For the Filtration of Aggressive Liquids



The holder hinders the release of trace elements into the filtrate and is resistant to almost all chemicals. The fluoroelastomer O-ring in the top part allows easy hand tightening, and can be replaced by a PTFE

O-ring, order no. 17039. The 6 mm outlet nipple is an integral part of the base, the 10 mm inlet hose nipple can be replaced by a G $\frac{3}{8}$  connector, order no. 17051.

### □ Specifications

|                                   |   |
|-----------------------------------|---|
| Chemical compatibility            | As for PTFE and fluoroelastomer   |
| Filtration area                   | 12.5 cm <sup>2</sup>  |
| Thread for inlet connector        | M 14 × 1.5 male thread  |
| Materials                         | Top part, barrel, base part: corrugated iron, hose nipples and filter support with 40 × 3.5 mm O-ring: PTFE, locking rings: aluminium 39 × 3.5 mm fluoroelastomer O-ring (top part) |
| Max. operating pressure           | 5 bar   500 kPa   72.5 psi  |
| Suitable membrane filter diameter | 47 mm   |
| Sterilization                     | By autoclaving (max 134°C) or by dry heat (180°C)   |


### □ Ordering Information

| Description   | Order No. |
|---|-----------|
| PTFE pressure filter holder, 47 mm, with 200 ml capacity. | 16579     |

### Replacement Parts

| Description              | Order No. |
|--------------------------|-----------|
| PTFE O-ring, 39 × 3.5 mm | 17039     |

Other replacement parts are shown in the diagram.



## Combisart® Manifolds

1-, 3- and 6-Branch



Made of high-grade stainless steel (B.S. 304S3 | AISI 304); accommodates any type of vacuum funnel. Stainless steel three-way valves (taps) allow the vacuum for each filter holder to be individually

controlled and each holder to be sterily vented. The low height of the manifold ports is particularly advantageous for working on a clean bench.

### Ordering Information

| <b>Combisart® Manifolds, without Base Support and Frit</b> | <b>Order No.</b> |
|--|------------------|
| Combisart® 1-branch manifold                               | 16844            |
| Combisart® 3-branch manifold                               | 16842            |
| Combisart® 6-branch manifold                               | 16843            |

| <b>Combisart® Sets, Stainless Steel Capacity</b> | <b>Order No.</b> |
|--|------------------|
| 1-branch 1 × 100 ml                              | 16844-CS         |
| 1-branch 1 × 500 ml                              | 16845-CS         |
| 3-branch 3 × 100 ml                              | 16824-CS         |
| 3-branch 3 × 500 ml                              | 16828-CS         |
| 6-branch 6 × 100 ml                              | 16832-CS         |
| 6-branch 6 × 500 ml                              | 16831-CS         |


In each set stainless steel funnels with lids are preassembled.

### Accessories and Replacement Parts

| <b>Description</b>  | <b>Pack Size</b> | <b>Order No.</b> |
|---|------------------|------------------|
| Plug, conical, to close the venting hole beside the 3-way valve | 10               | 6980225          |
| Silicone O-ring for manifold female threads                     | 3                | 6980235          |
| Rubber tubing, 1 m  | 1                | 16623            |

## ■ Glass Filter Holders; 30, 250 ml

For Particle Counting




### Glass Filter Holders

Two compact vacuum filter holders for easy particulate analysis. Both the top and bottom part of the filter holders are easily and securely fastened together using the

metal clamp. The centering rim on the filter support ensures correct positioning of the membrane filter. The glass frit filter support guarantees uniform distribution of retained particles on the filter surface.

### □ Ordering Information

| Description         |                      | Order No. |
|---------------------|----------------------|-----------|
| Glass filter holder | 30 ml                | 16306     |
| Filter diameter     | 25 mm (or 24 mm)     |           |
|                     | Prefilter, 20 mm     |           |
| Filtration area     | 3 cm <sup>2</sup>    |           |
| Capacity            | 30 ml                | 16307     |
| Outlet              | 12 mm outer diameter |           |
| Glass filter holder | 250 ml               |           |
| Filter diameter     | 47 mm (or 50 mm)     |           |
|                     | Prefilter, 40 mm     | 16307     |
| Filtration area     | 12.5 cm <sup>2</sup> |           |
| Capacity            | 250 ml               |           |
| Outlet              | 15 mm outer diameter |           |




### Adapter, 16836 | Adapter, 16837

For use of a glass filter holder, 16306 or 16307, on a Combisart® stainless steel manifold.


### □ Ordering Information

| Description   | Order No. |
|---|-----------|
| Adapter with 11 mm opening in stopper; for using filter holder 16306 on a Combisart® manifold | 16836     |
| Replacement stopper for 16836   | 00280     |
| Adapter with 14 mm opening in stopper; for using filter holder 16307 on a Combisart® manifold | 16837     |
| Replacement stopper for 16837   | 00281     |



# Polycarbonate Filter Holders


For Particle Counting



**Polycarbonate Filter Holder, 250 ml**  
This reusable, practical filter holder made of autoclavable plastic is ideal for analytical testing outside the laboratory. For use with 47 mm membrane filters.

Outlet: TR 20×2 mm male thread

## Ordering Information



| Description  | Order No. |
|--|-----------|
| Polycarbonate filter holder without receiver flask | 16511     |
| Polycarbonate filter holder with receiver flask    | 16510     |
| Hand vacuum pump with gauge and 60 cm PVC tubing   | 16672     |

## ■ Ready-to-Use Biosart® 250 Funnels

For Particle Counting



### Biosart® 250 Funnel

The Biosart® 250 Funnel has been specially designed for analytical quality assurance. The sterile 250 ml plastic funnel guarantees fast filtration and high sample throughputs

during routine testing. Its large inner diameter allows high flow rates, and the tapered inner walls permit thorough flushing of the funnel, after filtration.

### □ Ordering Information

| Description                                     | Order No.    |
|---|--------------|
| Biosart® 250 Funnel, 50 units, sterile-packaged | 16407-25-ALK |



### Single Support, 16840

For adapting a Biosart® 250 Funnel for use on a Combisart® stainless steel manifold.

### □ Ordering Information

| Description  | Order No. |
|--|-----------|
| Stainless steel filter support for stainless steel manifold. | 16840     |

### Replacement

| Description                                     | Order No. |
|---|-----------|
| Stainless steel frit for 50 mm membrane filters | 6980102   |
| Stainless steel frit for 47 mm membrane filters | 6980103   |
| Silicone flat gasket underneath the frit        | 6980124   |
| PTFE flat gasket underneath the frit            | 6980104   |
| Silicone O-ring for 16840 male thread           | 6980274   |

## ■ Did You know ...

... that you can order more lab catalogs?

Get an overview of our high-quality laboratory instruments, high-grade consumables and excellent services.

Our product portfolio focuses on laboratory instruments, such as lab balances, pipettes and laboratory water purification systems. Moreover, we offer the widest range of consumables like laboratory filters and pipette tips.



**Products for Microbiological Control**  
Publication No.: SL-1530-e  
Order No.: 85034-538-81



**Laboratory Weighing Products**  
Publication No.: WL-0007-e  
Order No.: 98649-018-13



**Sartorius Liquid Handling Products**  
Publication No.: SUL0002-e  
Order No.: 85032-542-25



**Laboratory Filtration Products**  
Publication No.: SLU0006-e  
Order No.: 85032-543-83



**Laboratory Water Purification Products**  
Publication No.: SL-0002-e  
Order No.: 85037-548-54



Download our lab catalogs from our website:  
[www.sartorius.com/lab-catalogs](http://www.sartorius.com/lab-catalogs)

# Sales and Service Contacts

For further contacts, visit [www.sartorius.com](http://www.sartorius.com)

## Europe

### Germany

Sartorius Lab Instruments  
GmbH & Co. KG  
Weender Landstrasse 94–108  
37075 Goettingen

Phone +49.551.308.0  
Fax +49.551.308.3289

### France & Suisse Romande

Sartorius France  
2, rue Antoine Laurent de Lavoisier  
ZA de la Gaudrée  
91410 Dourdan

Phone +33.1.70.62.50.00  
Fax +33.1.64.59.76.39

### Austria

Sartorius Austria GmbH  
Modecenterstrasse 22  
1030 Vienna

Phone +43.1.7965760.0  
Fax +43.1.7965760.24

### Belgium

Sartorius Belgium N.V.  
Leuvensesteenweg, 248/B  
1800 Vilvoorde

Phone +32.2.756.06.71  
Fax +32.2.253.45.95

### Finland & Baltics

Sartorius Biohit Liquid Handling Oy  
Laippatie 1  
00880 Helsinki

Phone +358.9.755.951  
Fax +358.9.755.95.200

### Hungary

Sartorius Hungária Kft.  
Kagyló u. 5.  
2092 Budakeszi

Phone +3623.457.227  
Fax +3623.457.147

### Ireland

Sartorius Ireland Ltd.  
Unit 41, The Business Centre  
Stadium Business Park  
Ballycoolin Road  
Dublin 11

Phone +353.1.8089050  
Fax +353.1.8089388

### Italy

Sartorius Italy S.r.l.  
Viale A. Casati, 4  
20835 Muggiò (MB)

Phone +39.039.4659.1  
Fax +39.039.4659.88

### Netherlands

Sartorius Netherlands B.V.

Phone +31.30.60.53.001  
Fax +31.30.60.52.917  
[info.netherlands@sartorius.com](mailto:info.netherlands@sartorius.com)

### Poland

Sartorius Poland sp.z o.o.  
ul. Wrzesinska 70  
62-025 Kostrzyn

Phone +48.61.6473830  
Fax +48.61.6473839

### Russian Federation

LLC "Sartorius ICR" and LLC "Biohit"  
Uralskaya str. 4, Lit. B  
199155, Saint-Petersburg

Phone +7.812.327.5.327  
Fax +7.812.327.5.323

### Spain & Portugal

Sartorius Spain, S.A.  
Avda. de la Industria, 32  
Edificio PAYMA  
28108 Alcobendas (Madrid)

Phone Spain +34.902.123.367  
Phone Portugal +351.800.855.800  
Fax Spain +34.91.358.96.23  
Fax Portugal +351.800.855.799

### Switzerland

Sartorius Mechatronics Switzerland AG  
Ringstrasse 24a  
8317 Tagelswangen (ZH)

Phone +41.44.746.50.00  
Fax +41.44.746.50.50

### U.K.

Sartorius UK Ltd.  
Longmead Business Centre  
Blenheim Road, Epsom  
Surrey KT19 9QQ

Phone +44.1372.737159  
Fax +44.1372.726171

### Ukraine

LLC "Biohit"  
Post Box 440 "B"  
01001 Kiev, Ukraine  
Phone +380.44.411.4918  
Fax +380.50.623.3162

## Americas

### USA

Sartorius Corporation  
5 Orville Drive, Suite 200  
Bohemia, NY 11716

Phone +1.631.254.4249  
Toll-free +1.800.635.2906  
Fax +1.631.254.4253

### Argentina

Sartorius Argentina S.A.  
Int. A. Avalos 4251  
B1605ECS Munro  
Buenos Aires

Phone +54.11.4721.0505  
Fax +54.11.4762.2333

### Brazil

Sartorius do Brasil Ltda  
Avenida Senador Vergueiro 2962  
São Bernardo do Campo  
CEP 09600-000 - SP- Brasil

Phone +55.11.4362.8900  
Fax +55.11.4362.8901

### Canada

Sartorius Canada Inc.  
2179 Dunwin Drive #4  
Mississauga, ON L5L 1X2

Phone +1.905.569.7977  
Toll-Free +1.800.668.4234  
Fax +1.905.569.7021

### Mexico

Sartorius de México S.A. de C.V.  
Circuito Circunvalación Poniente  
No. 149  
Ciudad Satélite  
53100, Estado de México  
México

Phone +52.5555.62.1102  
Fax +52.5555.62.2942

## Asia | Pacific

### Australia

Sartorius Australia Pty. Ltd.  
Unit 5, 7-11 Rodeo Drive  
Dandenong South Vic 3175

Phone +61.3.8762.1800  
Fax +61.3.8762.1828

### China

Sartorius Scientific  
Instruments (Beijing) Co., Ltd.  
33 Yu An Road, Airport Industrial Park  
Zone B,  
Shunyi District, Beijing 101300,  
P.R.China

Phone +86.10.8042.6300  
Fax +86.10.8042.6486

### Hong Kong

Sartorius Hong Kong Ltd.  
Unit 1012, Lu Plaza  
2 Wing Yip Street  
Kwun Tong  
Kowloon, Hong Kong

Phone +852.2774.2678  
Fax +852.2766.3526

### India

Sartorius Weighing India Pvt. Ltd.  
#69/2-69/3, NH 48, Jakkasandra,  
Nelamangala Tq  
562 123 Bangalore, India

Phone +91.80.4350.5250  
Fax +91.80.4350.5253

### Japan

Sartorius Japan K.K.  
4th Fl., Daiwa Shinagawa North Bldg.  
8-11, Kita-Shinagawa 1-chome  
Shinagawa-ku, Tokyo, 140-0001 Japan

Phone +81.3.3740.5408  
Fax +81.3.3740.5406

### Malaysia

Sartorius Malaysia Sdn. Bhd  
Lot L3-E-3B, Enterprise 4  
Technology Park Malaysia  
Bukit Jalil  
57000 Kuala Lumpur, Malaysia

Phone +60.3.8996.0622  
Fax +60.3.8996.0755

### Singapore

Sartorius Singapore Pte. Ltd  
1 Science Park Road,  
The Capricorn, #05-08A,  
Singapore Science Park II  
Singapore 117528

Phone +65.6872.3966  
Fax +65.6778.2494

### South Korea

Sartorius Korea Ltd.  
8th Floor, Solid Space B/D,  
PanGyoYeok-Ro 220, BunDang-Gu  
SeongNam-Si, GyeongGi-Do, 463-400

Phone +82.31.622.5700  
Fax +82.31.622.5799

### Thailand

Sartorius (Thailand) Co. Ltd.  
129 Rama 9 Road,  
Huaykwang  
Bangkok 10310

Phone +66.2643.8361-6  
Fax +66.2643.8367



◀ [www.sartorius.com](http://www.sartorius.com)